# Application News

No. SSI-LCMS-082K

## Liquid Chromatograph Mass Spectrometer, LCMS™-8060

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### Separation and Quantitation of Oxysterol Compounds Using LC-MS/MS

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Compound name	Abbreviation
24 (S)-hydroxycholesterol	24HC
(D7)22-hydroxycholesterol	22HC(d7)
25-hydroxycholesterol	25HC
27-hydroxycholesterol	27HC
(D7)7α-hydroxycholesterol	7αHC(d7)
7α-hydroxycholesterol	7αHC
7β-hydroxycholesterol	7βHC
7-Ketocholesterol	7KC
(D7) 7-Ketocholesterol	7KC(d7)
7α-hydroxycholestenone	7αHCn
(D3) Vitamin D3	VitD3(d3)
Zyosterol	Zymo
Desmosterol	Desmo
7α,27 dihydroxycholestenone	7α,27diHC,3one
Cholesterol	CH
7dehydrocholesterol	7DHC

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SHIMADZU Scientific Korea Corp.  
www.shimadzu.co.kr

For Research Use Only. Not for use in diagnostic procedures. Not available in the USA, Canada, and China. This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.