Application News

No. 04-AD-0262-ENK

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Instrument	: TOC-L _{CSH} TOC analyzer and SSM-5000A solid sample combustion unit
Cell length	: Short cell
SSM Carrier Gas	: Oxygen, 500 mL/min
Oxidation Method	: Cobalt Oxide/Platinum catalytic combustion 980 °C
Method	: NPOC
Calibration Curve	: 2 % TC standard solution prepared from Glucose(C ₆ H ₁₂ O ₆)

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Instrument	: TOC-L _{CSH} TOC analyzer
Accessory	: High Suspension Kit(High Conc TC)
Oxidation Method	: Standard Platinum catalyst, combustion, 680 °C
Method	: NPOC
Calibration Curve	: 2-20 mg/L NPOC standard solution prepared from Potassium Hydrogen Phthalate(C ₈ H ₅ O ₄)

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[1] Avramidis, P., Nikolaou, K. and Bekiari, V. (2015). TOC and Total Nitrogen in Sediments and Soils. Agriculture and Agricultural Science Procedia 4: 425 – 430.
 [2]. BS EN 15936 (2012). Sludge, Treated Biowaste, Soil & Waste – Determination of TOC by Dry Combustion.
 [3]. Goto, A. (2020). Measurement of total carbon in soil and compost by TOC solid sample measurement system. Shimadzu Application News no 078.
 [4]. SCA-130-105 (2019). TOC determination in solid samples using the suspension method. Application Handbook – Sum Parameter (SEG-A-099).
 [5]. TOC-L for Monitoring a Microscopic Algae Biomass (2014). Shimadzu C391-E096.
 [6]. Ikezawa, Y. (2017). Carbon analysis of blue-green algae suspensions. Shimadzu Application News no 061.