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■ Improved Selectivity of Skin Sensitization Test ADRA with RF-20Axs and i-PDeAll

ADRA, as a skin sensitization test, that evaluates the allergic reaction of chemicals to the skin using N-(2-(1-naphthyl)acetyl)-L-cysteine(NAC) and α -N-(2-(1-naphthyl)acetyl)-L-lysine(NAL) which contains naphthalene rings and cysteine or lysine. An example analysis according to the OECD Test Guideline is given in AN_01-00258. Depending on the test substance, coelution peaks may appear near the NAC and the NAL peaks. In this case, LC separation conditions should be considered.

The ADRA test must be performed within 72 hours from sample preparation.

If the co-elution peak appears, the test must be repeated from the beginning. In this report, we introduce two methods to avoid the interferences of coelution peaks by increasing the selectivity of detection without changing the analysis. One example is the use of a fluorescence detector(RF), which is more selective than a UV detector, to avoid the interferences of coelution. Another example is an analysis using a photodiode array detector(PDA) with the deconvolution function i-PDeA II to obtain separated peaks.*)

■ Simultaneous Analysis of Dipotassium Glycyrrhizinate and Tranexamic Acid

Glycyrrhizic acid is a medicinal compound found in Chinese crude medicines such as licorice and is generally known to have pharmacological effects such as anti-allergy, anti-inflammation, and detoxification. Dipotassium glycyrrhizate (GK2), a salt from glycyrrhizic acid is widely used for products such as cosmetics, shampoos, and toothpastes, in addition to over-the-counter drugs. Tranexamic acid (TA), an artificially synthesized amino acid, is generally known for its anti-inflammatory, hemostatic, and other pharmacological effects. Additionally, it is also used in cosmetics as an active whitening agent.

In general, it is difficult to retain TA on C18 columns due to its high polarity, and it often show tailing peak shape due to its basicity. On the other hand, GK2 is a compound that is relatively well retained on C18 columns. Therefore, this article introduces of TA and GK2 in cosmetics employing sodium perchlorate as a nt of mobile phase for simultaneous analysis of these two compounds.

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■ Quantitative Analysis of 57 Fragrance Allergens in Cosmetics Using Twin Line MS System

Under the European Cosmetics Directive (EC 1223/2009). 24 of the fragrances used in cosmetics are regulated as fragrance allergens. If a leave-on product contains more than 0.001 % or a rinse-off product contains more than 0.01 %, it must be labeled on the product. Recently, however, following an opinion by the Scientific Committee on Consumer Safety (SCCS), the number of regulated compounds has been scheduled to be expanded to more than 80. Therefore, manufacturers of cosmetics need to know exactly what fragrance compounds are contained in finished products or raw materials.

This Application News describes an example of analysis of 57 fragrance allergens in a commercial hair oil using the twin line MS system, based on the analysis method reported by IFRA.

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■ Evaluation of Structural Strength of Semisolids Conforming to USP 915 and USP 1912 by Penetration Method

Methods for characterization of the structural strength or consistency of semisolid materials which are scheduled for inclusion in the United States Pharmacopeia USP 915 and USP 1912, are useful for understanding the performance and quality stability of various pharmaceutical products. These USP chapters have shown that comparative evaluations of characteristics between semisolid materials are possible through evaluation of structural strength. As representative characteristics of pharmaceutical products related to structural strength, the ease of application and quality of spreading of ointments and lotions, resistance to dripping of liquid preparations for external use, and ease of application of adhesives for transdermal patches can be mentioned.

In the past, majority of evaluations of the structural strength of pharmaceutical products were conducted by sensory evaluation. However, the need for quantitative evaluation of the structural strength of pharmaceuticals in order to ensure a stable supply of products with appropriate viscosity has become apparent in recent years, but until now, there had been no officially recognized method for quantitative evaluations of structural strength in the pharmaceutical industry. By providing measurement methods related to the characterization of pharmaceutical products in USP 915 and USP 1912, it is hoped that techniques which enable quantitative evaluation of the viscoelasticity of pharmaceutical products will gain official recognition, and a stable supply of products and more accurate comparative evaluations of products under development will be possible.

Four evaluation methods for pharmaceutical products, namely, Strain ramp measurement, Shear rate ramp measurement, Oscillation amplitude sweep measurement, and Penetrometry measurement, are scheduled for inclusion in USP 915. Of these four methods, this article introduces the constant-speed measurement method as one evaluation method by penetrometry measurement. Because a predetermined conical test jig geometry is used with this method, it is not necessary to select a jig, and since the test is to measure the cone penetration by the penetration depth of the conical test jig, the test analyst can evaluate the data by a simply, easy-to-understand method.

In connection with penetration depth, details are provided in USP 1912, and the various characteristic values, including penetration depth, are explained from the theoretical viewpoint referring to the standard test methods ASTM D217 and ASTM D937, which were established for the semisolid substances such as industrial grease and petrolatum, respectively.

This article introduces examples of application of the penetration method in tests of a medical hand cream, which is sold commercially as a medicine, and a toothpaste known for displaying speed-dependency of viscosity as the pharmaceutical products.

■ Particle Size Evaluation in Hair Conditioner –Undiluted Measurement of Highly Concentrated Emulsions–

Emulsions are liquid mixtures of two immiscible solvents, such as oil and water, consisting of a dispersion of microparticles of one liquid in the other. Familiar examples of emulsions include hair conditioners, suntan lotion, and mayonnaise. The effectiveness of cosmetics, the texture of food products, and other characteristics can vary depending on the size of the particles in the emulsion. Therefore, product quality and other features can be quantitatively evaluated by measuring the particle size distribution in emulsions.

However, due to the thermodynamic instability of many emulsions, slight differences in how they are prepared can change their state or stability. This article describes an example of using a Shimadzu SALD-2300 laser diffraction particle size analyzer to measure undiluted and diluted emulsions.

■ Webinar@ The Newest Alternative to Animal Testing for Skin Sensitization

In recent years, there has been a growing movement to promote alternative methods to animal testing and the EU and a few other countries have even moved to ban the import and sale of cosmetics that have been tested on animals. The Amino acid Derivative Reactivity Assay (ADRA) Method developed by Fujifilm, is an innovative alternative test method for skin sensitization without the use of animals and it was adopted by the OECD in 2022

as a reliable test method.

In this webinar, we will present the technical background of the ADRA method, the use of reagent kits to perform ADRA testing, and our solutions for ADRA analysis that deliver both enhanced selectivity and high-speed testing. Register for this webinar and learn more about how the ADRA method will help you better test the safety of cosmetics and pharmaceuticals.

Key Learning Objectives:

- Outline of the OECD TG442C (ADRA: Amino acid Derivative Reactivity Assay)
- Overview of ADRA kits and related products
- How to conduct ADRA testing using HPLC

Who Should Attend:

- R&D and Manufacturer of Cosmetics
- Companies considering Alternative Methods to Animal Testing
- Industries interested in Alternative Methods to Animal Testing

► For more information and Applications.

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OTC Drug Manufacturer대상 Cosmetic Campaign Price 실시

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**SHIMADZU** | Shimadzu Scientific Korea Corp.