

## Application News

# No. C146A

### ■ nSMOL<sup>TM</sup> Antibody BA Kit Features

nSMOL is Shimadzu's completely new and breakthrough technology that enables selective proteolysis of the Fab region of monoclonal antibodies. This technique facilitates method development independent of a variety of antibody drugs and achieves a paradigm shift in the bioanalysis of antibody drugs.

nSMOL<sup>™</sup> Antibody BA Kit

- Bevacizumab analysis -

Furthermore, this is the only method with respect to antibody drugs that has fulfilled the criteria of "Guideline on Bioanalytical Method Validation in Pharmaceutical Development" for low MW drug compounds issued by the Japanese Ministry of Health, Labour and Welfare. Shimadzu also offers optimization methods and protocols, and nSMOL can be applied to clinical research at various institutions.

### Antibody Drug Classification and Selection of Quantitation Peptides

Monoclonal antibodies are produced from the hybridoma with mouse spleen lymphocyte and myeloma cells. While mice are predominantly used as hosts, in recent years a variety of hosts are now available to produce monoclonal antibodies.

Furthermore, production of the variable Fv by phage display technology and high-throughput screening of affinity sequence has become alternative standard procedure.

Antibody drugs are classified into four classes according to specific structure.

The complementarity-determining region (CDR) of antibody specificity against human IgGs becomes smaller according to the "mouse"  $\rightarrow$  "chimeric"  $\rightarrow$  "humanized"  $\rightarrow$  "fully human" antibody. More precise selection of quantitation peptides becomes particularly important in the nSMOL, which is used to perform structure specificity-indicated analysis.

The nSMOL enables selective proteolysis in variable regions. This allows selection of quantitation peptides that reflect the structural characteristics of antibodies. Antibodies have three CDRs on each heavy and light chain, CDR2 is known as the region that makes first contact with an antigen. The signature peptide by nSMOL are mainly from CDR2 containing peptides.



Fig. 1 Antibody Drug Classification

### Analysis Conditions for Bevacizumab Using the nSMOL

<Sample Processing Protocol>

LCMS Bioanalysis of Antibody Drugs Using

Fab-Selective Proteolysis nSMOL-Part 2

With the nSMOL technique, the same sample processing protocol can be applied to all antibody drugs. For details, refer to Shimadzu Application News (Trastuzumab analysis).

<LCMS Analysis Conditions>

[I C] NexeraX2 System	n
Column	: Shim-pack GISS C18 (50 mm $\times$ 2.1 mm)
Column oven	: 50 °C
Solvent A	: 0.1 % formic acid/water
Solvent B	: 0.1 % formic acid/acetonitrile
Gradient	: 1 %B (1.5 min)/1-35 %B (3.5 min)/
	95 %B (1 min)/1 %B (1 min)
Flow rate	: 0.4 mL/min
Injection	: 10 μL
[MS] LCMS-8050, 806	0
lonization	: ESI Positive
DL	: 250 °C
Heat Block	: 400 °C
Interface	: 300 °C
Nebulizer gas	: 3 L/min
Drying gas	: 10 L/min
Heating gas	· 101/min

#### Bevacizumab Quantitation Peptides

Averaged accuracy

Peptide	MRM transition	Purpose	
P <sub>14</sub> R	512.1>292.3 (b3+) 512.1>389.3 (b4+) 512.1>660.4 (b6+)	For quantitation (IS) For structure confirmation For structure confirmation	
FTFSLDTSK	523.3>797.4 (y7+) 523.3>898.5 (y8+) 523.3>650.3 (y6+)	For quantitation For structure confirmation For structure confirmation	
STAYYLQMN SLR	642.3>748.4 (y6+) 642.3>861.5 (y7+) 642.3>620.3 (y5+)	For quantitation For structure confirmation For structure confirmation	
VLIYFTSSLH SGVPSR	588.3>775.9 (y14++) 588.3>602.3 (y6+) 588.3>939.5 (y9+)	For quantitation For structure confirmation For structure confirmation	
* Ouantitation range in human plasma		: 0.15 to 300 µg/ml	

: 101.3 %



Fig. 4 MRM Chromatograms of VLIYFTSSLHSGVPSR (in Human Plasma)

### Full Validation Results for Bevacizumab

<Precision and accuracy>

Set Concentration [µg/ml]	Data Average (N = 15)	Accuracy (%)	CV (%)
0.439	0.464	106	11.7
240	235	98.1	4.45
< Franza that tast			

<Freeze-thaw test>

Set Concentration [µg/ml]	Data Average (N = 5)	Accuracy (%)	Temperature (°C)
0.439	0.395	89.9	-20
240	253	106	-20

### <Long-term stability test>

Set Concentration [µg/ml]	Data Average (N = 5)	Accuracy (%)	Temperature (°C)
0.439	0.477	109	-20
240	223	93.1	-20

### <Processed sample stability for 48 h>

-				
	Set Concentration [µg/ml]	Data Average (N = 5)	Accuracy (%)	Temperature (°C)
Ī	0.439	0.445	101	5
	240	245	102	5





### Observations, Conclusions, and References

Although Bevacizumab quantitation peptide using nSMOL was obtained from the same region of Trastuzumab, the optimal peptide sequences for bioanalysis will depend on the interference with endogenous IgGs.

With respect to multiplexed quantitation of three sequences, nSMOL bioanalysis for Bevacizumab fulfilled the full validation criteria.

The lower limit of quantitation is  $0.15 \,\mu$ g/ml and the same assay method can be used from preclinical to clinical trials.

<References>

Shimadzu Corporation

www.shimadzu.com/an/

lwamoto N et al. Analyst, 2014, DOI:10.1039/c3an02104a

Iwamoto N et al., Drug Metab Pharmacokinet., 2016, DOI:10.1016/j.dmpk.2015.11.004 <Chief Scientists>

Noriko Iwamoto, Ph.D. and Takashi Shimada, Ph.D., Technology Research Laboratory

Notes: The product described in this document has not been approved or certified as a medical device under the Pharmaceutical and Medical Device Act of Japan.

It cannot be used for the purpose of medical examination, treatment or related procedures.



### For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "©". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

First Edition: Mar. 2017