

Application News

No. **G312**

Gas Chromatography

Testing the Gasoline Dilution Rate of Gasoline in Engine Oil in Accordance with ASTM D3525 and Japan Petroleum Institute Standard JPI-5S-24

If gasoline or diesel mixes into the engine oil, it decreases the oil viscosity and prevents achieving the proper performance as a lubricating oil. Measuring the fuel dilution rate serves as a key indicator during oil replacement, because it can determine the degradation status of engine oil. The test methods used to measure the fuel dilution rate are specified in standards such as U.S. ASTM D3524, D3525, and D7593. The standard for gasoline is ASTM D3525. It corresponds to the Japan Petroleum Institute standard JPI-5S-24 in Japan. This article describes an example of measuring the dilution rate of gasoline in engine oil in accordance with ASTM and JPI standards.

A. Miyamoto, R. Kubota, T. Wada

Preparing Samples

 $20~\mu L$ of the internal standard substance n-C $_{14}{}^{*1}$ was added to a 1 g sample of engine oil that contains gasoline.

Mixture samples for measuring the column resolution were prepared by adding 1 % each, by volume, of n-C $_{14}$ and n-C $_{16}$ *3 to n-C $_{8}$ *2.

Each sample was analyzed without dilution with a solvent.

■ Analytical Conditions

Analytical conditions specified based on respective standards are indicated in Table 1.

Table 1 Analytical Conditions

Model: Nexis™ GC-2030 AF/AOC-20i

Column: SH-RtxTM-1 (5 m × 0.53 mm I.D., df = 1.00 μ m)

Column Temp.: 50 °C (0 min) - 45 °C/min - 290 °C (10 min)

Total: 5.33 min

Injection Temp.: 255 °C
Carrier Gas: N₂, 15 mL/min
Purge Flow: 3 mL/min
Injection Method: Split 1:5

Carrier Gas Controller: Constant linear velocity mode


Detector: FID
Detector Temp.: 300 °C
Injection Volume: 0.1 µL

*1: Tokyo Chemical Industry Co., Ltd., 99.0 % or higher

*2: FUJIFILM Wako Pure Chemical Corporation, 98.0 % or higher

*3: Tokyo Chemical Industry Co., Ltd., 98.0 % or higher

■ Chromatogram of Engine Oil Containing Gasoline



■ Precautions for Injecting Samples

The given analytical conditions do not require any pretreatment steps, such as dilution with n-C $_{16}$, carbon disulfide (CS $_2$), or other solvent. A small quantity of highly viscous samples is injected, but the results confirm that adequate accuracy can be achieved by modifying the consumables or injection actions involved. Sample injection conditions are summarized in Table 2.

Table 2 Sample Injection Conditions

| Syringe: | 0.5 µL volume syringe (P/N 000445 in Fig. 2)*4 |
|---------------------------|--|
| Rinse Solvent: | CS ₂ |
| Plunger Aspiration Speed: | Low |
| Sample Rinsing: | None |
| Number of Pump Strokes: | Zero times |
| Insert: | P/N 227-35007-01 (Fig. 3) Wool is positioned 18 mm from top |

^{*4:} If new, repeatedly aspirate and discharge solution to ensure the plunger slides smoothly before use.



Fig. 2 0.5 µL Syringe (P/N 000445)



Fig. 3 Deactivated Insert for Split Injection (P/N 227-35007-01)

■ Confirming the Column Resolution




Fig. 4 Chromatogram of Column Resolution Measurement Sample

Results from analyzing column resolution measurement samples are shown in Fig. 4. The column resolution between n-C₁₄ and n-C₁₆ was 7.3. Results confirmed that the resolution between n-C₁₄ and n-C₁₆ in samples for measuring column resolution is at least 3 and not more than 8 (USP), as required by ASTM standards.

■ Repeatability of Gasoline Dilution Rates

The repeatability of gasoline dilution rates is shown in Table 3. Excellent repeatability %RSD (n = 10) results were obtained. The results also confirmed that all within-laboratory accuracy values satisfy tolerances required by standards.

Table 3 Repeatability %RSD (n = 10) of Gasoline Dilution Rates (%)

| | Sample 1 | Sample 2 | Sample 3 |
|---------|----------|----------|----------|
| 1 | 0.66 | 4.91 | 9.91 |
| 2 | 0.66 | 4.87 | 9.89 |
| 3 | 0.66 | 4.86 | 9.97 |
| 4 | 0.65 | 4.89 | 9.90 |
| 5 | 0.65 | 4.87 | 9.80 |
| 6 | 0.65 | 4.84 | 9.76 |
| 7 | 0.65 | 4.82 | 9.76 |
| 8 | 0.65 | 4.88 | 9.73 |
| 9 | 0.65 | 4.86 | 9.71 |
| 10 | 0.65 | 4.85 | 9.64 |
| Average | 0.65 | 4.87 | 9.81 |
| %RSD | 0.74 | 0.52 | 1.08 |

Summary

This analysis achieved accuracy levels required by ASTM D3525 and JPI-5S-24 using the indicated analytical conditions with a nitrogen carrier gas, without involving dilution with a solvent or other pretreatment steps. These standards specify requirements for using the backflush method to eliminate residual components. For more information about high-throughput analysis for gasoline dilution rate testing using backflushing, refer to Application News No. G313 regarding ASTM D7593. Other Application News bulletins related to fuel dilution rates are indicated in the list of references below.

List of References

| Standard | Item Analyzed | Application News No. |
|-----------|---------------------|----------------------|
| D3524 | Diesel | G310 |
| JPI-5S-23 | Diesel | G311 |
| D3525 | Gasoline | G312 |
| JPI-5S-24 | Gasoline | G312 |
| | Gasoline | G313 |
| D7593 | Diesel Biodiesel | G314 |

Reference Documents ASTM D3525-04 JPI-5S-24-2017

Nexis is a trademark of Shimadzu Corporation in Japan and/or other countries.

Rtx is either a trademark or a registered trademark of Restek Corporation in the United States and/or other countries.



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See http://www.shimadzu.com/about/trademarks/index.html for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shirmadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shirmadzu on or before the date of publication, and subject to change without notice.

First Edition: Nov. 2019