

# Application News

High Performance Liquid Chromatography

No.L416

## Analysis of Hexafluorophosphate Ion in Electrolytic Solution for Lithium-Ion Rechargeable Batteries by Ion Chromatography

Lithium hexafluorophosphate is a common electrolyte used in rechargeable batteries. Here, we introduce an example of the analysis of hexafluorophosphate ion in

the electrolytic solution of a lithium-ion rechargeable battery with the Shimadzu Prominence HIC-SP suppressor-type ion chromatograph.

### ■ Analysis of Standard Solution

Fig. 1 shows the chromatogram obtained from analysis of a hexafluorophosphate ion standard solution, and Table 1 shows the analytical conditions.

Due to the extremely strong retention of the hexafluorophosphate ion on the Shim-pack IC-SA2 column which is the standard column for the Prominence HIC-SP, the Shim-pack IC-SA2(G) guard column was also used. Moreover, instead of using the difficult to handle lithium hexafluorophosphate as the standard solution, ion analysis was conducted using potassium hexafluorophosphate instead.

Table 1 Analytical Conditions

Column	: Shim-pack IC-SA2(G) (10 mm L. × 4.6 mm I.D.)
Mobile Phase	: 12 mmol/L NaHCO <sub>3</sub> , 0.6 mmol/L Na <sub>2</sub> CO <sub>3</sub>
Flow Rate	: 1.0 mL/min
Injection Vol.	: 10 μL
Column Temp.	: 30 °C
Detection	: CDD-10A <sub>SP</sub> (Suppressor)




Fig. 1 Chromatogram of Hexafluorophosphate Ion

### ■ Linearity

Fig. 2 shows the calibration curve generated using hexafluorophosphate ion standard solutions with concentrations from 2.5 to 250 mg/L. Excellent linearity was obtained with an  $R^2$  value greater than 0.999.




Fig. 2 Linearity

### ■ Repeatability

Table 2 shows the retention time and peak area relative standard deviation ( $n = 6$ ) obtained from analysis of a 100 mg/L standard solution. Excellent repeatability of both retention time and peak area were obtained.

Table 2 Repeatability of Retention Time and Peak Area

	Retention Time	Peak Area
1st	4.78	102135
2nd	4.78	102422
3rd	4.79	100941
4th	4.77	103548
5th	4.77	101649
6th	4.77	102475
%RSD	0.16	0.86

### ■ Analysis of Electrolytic Solution for Lithium-Ion Rechargeable Battery

Fig. 3 shows the chromatogram obtained from the analysis of an electrolytic solution in a commercially available lithium-ion rechargeable battery. The electrolytic solution was analyzed after diluting it 1000 to 1 with pure water, and then filtering it through a membrane filter.

Table 3 shows the quantitation result. The concentration of this particular lithium ion rechargeable battery electrolytic solution was reported to be 1 mol/L, and the actual value obtained using the above-mentioned analytical conditions was 1.00 mol/L (144.6 mg/L). Lithium hexafluorophosphate is hydrolyzable, but good quantitation results were obtained under these conditions.

Table 3 Quantitative Value

	Displayed Rating Value	Analytical Value
LiPF <sub>6</sub> Concentration	1 mol/L	1.00 mol/L




Fig. 3 Chromatogram of Electrolytic Solution for a Lithium-Ion Rechargeable Battery

[Referenc]

K.Tasaki, K.Kanda, S.Nakamura, M.Ue : Journal of The Electrochemical Society, 150 (12) A1628-A1636 (2003)

 **SHIMADZU**

SHIMADZU CORPORATION. International Marketing Division

3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan Phone: 81(3)3219-5641 Fax: 81(3)3219-5710