

Application  
Data Sheet

No.20

System Gas Chromatograph

High Sensitive CO, CO<sub>2</sub>, CH<sub>4</sub> Analysis  
Nexis GC-2030CCC4  
GC-2014CCC4

This system is designed to measure a trace amount of carbon monoxide (CO), methane (CH<sub>4</sub>) and carbon dioxide (CO<sub>2</sub>) in a gas sample, such as He, H<sub>2</sub>, N<sub>2</sub> and Ar. The sample is injected automatically through a 10-port valve. First, a Porapak-N pre-column is used to cut the above C<sub>2</sub> compounds. The Porapak functions to separate CO/CH<sub>4</sub> and CO<sub>2</sub>. CO and CH<sub>4</sub> are separated by an MS-13X column, while CO<sub>2</sub> moves through the Porapak-Q. CO/CH<sub>4</sub> and CO<sub>2</sub> are then combined before a methanizer. CO and CO<sub>2</sub> are reduced to CH<sub>4</sub> by means of a nickel catalyst and detected by a flame ionization detector (FID). If the matrix contains O<sub>2</sub>, this concentration should be less than 0.1% to protect the catalyst from damage. The system includes LabSolutions GC workstation software.

**Analyzer Information**

**System Configuration:**

Two valves / four packed columns /  
Methanizer with FID detector

**Sample Information:**

CO, CO<sub>2</sub>, CH<sub>4</sub>

**Concentration Range:**

No.	Name of Compound	Concentration Range	
		Low Conc.	High Conc.
1	CO	1.0ppm	100ppm
2	CO <sub>2</sub>	1.0ppm	100ppm
3	CH <sub>4</sub>	1.0ppm	100ppm

Detection limits may vary depending on the sample.  
Please contact us for more consultation.

**System Features**

- Single channel with packed columns
- Hydrocarbons and water are backflushed by the pre-column while trace CO, CO<sub>2</sub>, CH<sub>4</sub> reach FID.
- Good separation between CH<sub>4</sub> and CO with MS-13X packed column
- 13 minutes analysis time

**Typical Chromatograms**




Fig. 1 Chromatogram of FID

First Edition: November, 2017