

Application  
Data Sheet

No.84

System Gas Chromatograph

H2O in Gas Sample Analysis System  
Nexis GC-2030H2O  
GC-2014H2O

This instrument is designed for the determination of H2O by gas chromatography (GC) and detection by TCD. A Porapak-Q pre-column is used, Porapak-Q and Sunpak S are used to separate the water from the matrix. LabSolutions chromatography software handles all aspects of GC control, automation, and data handling.

**Analyzer Information**

**System Configuration:**

Two valves / two packed columns with one TCD detector

**Sample Information:**

H<sub>2</sub>O

**Concentration Range:**

No.	Name of Compound	Concentration Range		Detector
		Low Conc.	High Conc.	
1	H2O	0.05%	50.00%	TCD-1

Detection limits may vary depending on the sample. Please contact us for more consultation.

**System Features**

- 7 minutes analysis for H2O analysis can be carried out
- One TCD channel
- Good repeatability

**Typical Chromatograms**




Fig. Chromatogram of TCD

First Edition: November, 2017