




# Nexera Application Data Sheet No.12

## Ultra-High-Resolution Analysis of Polycyclic Aromatic Hydrocarbons

As the analysis of polycyclic aromatic hydrocarbons in the environment demands the highly efficient resolution of multiple components, long sub-2 µm columns or multiple columns are often used. Such cases demand a UHPLC system with a high system pressure tolerance. The 130 MPa pressure tolerance of the Shimadzu Nexera can safely and flexibly accommodate such demands. This Application Data Sheet introduces the ultra-high-resolution simultaneous analysis of polycyclic aromatic hydrocarbons using Nexera with two Phenomenex Kinetex columns.

### Simultaneous Analysis of 16 Polycyclic Aromatic Hydrocarbon Components

Resolution was investigated using a standard mixture of 16 polycyclic aromatic hydrocarbon components (1 to 20 mg/L, acetonitrile solution). The performance was compared for one and two Phenomenex Kinetex columns with 2.6 µm particle size. (Kinetex columns are core-shell columns with a 0.35 µm porous membrane bonded to a 1.9 µm solid core.) Connecting two columns achieved 1.5 minimum resolution for all components. The fluorescence detector achieved high sensitivity detection of all components, except Component 2.



Peaks :	
1. Naphthalene, 2. Acenaphthylene, 3. Fluorene, 4. Acenaphthene, 5. Phenanthrene, 6. Anthracene, 7. Fluoranthene, 8. Pyrene, 9. Chrysene, 10. Benzo(a)anthracene, 11. Benzo(b)fluoranthene, 12. Benzo(k)fluoranthene, 13. Benzo(a)pyrene, 14. Dibenz(a,h)anthracene, 15. Indeno(1,2,3-cd)pyrene, 16. Benzo(g,h,i)perylene	