




Nexera Application Data Sheet No.13

Ultrafast Analysis of Aldehydes and Ketones

Aldehydes and ketones in the environment can be analyzed using HPLC by derivatizing them with 2,4-dinitrophenylhydrazine (2,4-DNPH). For multicomponent separation of aldehyde/ketone DNPH derivatives, reverse-phase gradient elution is generally used. Nexera enables stable ultrafast gradient elution through accurate solution delivery and the use of a high-efficiency gradient mixer. This document introduces an example of ultrafast analysis performed on 13 aldehyde/ketone DNPH derivatives by gradient elution using Nexera and a Phenomenex Kinetex C18 column (particle size 2.6 μ m, a core-shell column where a 0.35 μ m porous membrane is combined with a 1.9 μ m solid core).

Batch analysis of 13 aldehyde/ketone DNPH derivatives by binary high-pressure gradient elution

A batch analysis was performed on standard mixtures (1 μ g/mL and 5 μ L injection each) by binary high-pressure gradient elution with water/THF mixture and acetonitrile. In this analysis, the repeatability (n = 6, %RSD) of the retention time and peak area of formaldehyde-DNPH derivatives were 0.087 % and 0.16 % respectively, showing Nexera's high-accuracy gradient performance.



Column : Phenomenex Kinetex 2.6 µm C18 100 Å

 $\label{eq:mobile phase} \begin{array}{ll} & (50 \text{ mmL. x } 3.0 \text{ mmI.D., } 2.6 \text{ } \mu\text{m}) \\ \text{Mobile phase} & : A: \text{Water/THF=8/2} \end{array}$

B: Acetonitrile


Gradient : B 20 % (0 min)→50 % (4 min) Flow rate : 1.5 mL/min (Mixer 180 μL)

Column temp. : $50\,^{\circ}\text{C}$ Injection volume : $5\,\mu\text{L}$ Detection : UV 360 nm
Flow cell : Semi-micro cell

Batch analysis of 13 aldehyde/ketone DNPH derivatives by ternary high-pressure gradient elution

Nexera's modular design enables various systems to be configured, such as a high-pressure gradient elution system with up to three solutions or a low-pressure gradient elution system with up to four solutions. Below is an example where standard mixtures (1 μ g/mL and 2 μ L injection each) were efficiently analyzed in a shorter time using ternary high-pressure gradient elution.

Gradient



Column : Phenomenex Kinetex 2.6 μm C18 100 Å

(50 mmL. x 2.0 mml.D., 2.6 μm)
Mobile phase : A: Water

B: Acetonitrile

C: THF

: B 20 % (0 min)→50 % (2.0-2.4 min) C 15 % (0 min)→10 % (2.0-2.4 min)

Flow rate : 0.9 mL/min (Mixer 180 µL)

Column temp. : $50\,^{\circ}\text{C}$ Injection volume : $2\,\mu\text{L}$ Detection : UV 360 nm
Flow cell : Semi-micro cell

Peaks:

1. Formaldehyde, 2. Acetaldehyde, 3. Acetone, 4. Acrolein, 5. Propionaldehyde, 6. Crotonaldehyde, 7. 2-Butanone, 8. Methacrolein, 9. *n*-Butyraldehyde, 10. Benzaldehyde, 11. Valeraldehyde, 12. *m*-Tolualdehyde, 13. Hexaldehyde (as their 2,4-DNPH derivatives)

