

# Application News

## No. C230

LC-MS

### Analysis of Favipiravir in Human Plasma Using Fully Automated Sample Preparation LC/MS/MS System

#### Introduction

Favipiravir (brand name: Avigan®), which was developed by FUJIFILM Toyama Chemical Co., Ltd, is one of the RNA polymerase inhibitors used for treating influenza. In Application News C229, we introduced a robust, highly-sensitive analysis using LC/MS/MS with manual pretreatment. However, manual pretreatment of plasma samples entails a certain level of workload. This report introduces a method of analyzing favipiravir using a fully automated sample preparation LC/MS/MS system that can reduce variation between procedures, sample mix-ups, and risk of exposure to the samples (Fig. 1).

E. Imoto, D. Kawakami



Fig. 1 Fully Automated Sample Preparation LC/MS/MS System (CLAM™+LC/MS/MS)

#### Fully Automated Sample Preparation of Favipiravir in Plasma

For analysis of low-molecular weight compounds in plasma using a LCMS™, it is common to use supernatant collected following deproteinization by adding an organic solvent. With the fully automated sample preparation LC/MS/MS system, these preparatory steps are done automatically just by placing a blood collection tube in the system after plasma separation (Fig. 2). Pretreatment of the next sample can also be performed in parallel with LC/MS/MS analysis, which can greatly reduce the time required to analyze each sample.

This analysis was performed in a per-sample cycle time of 6.5 minutes from plasma pretreatment to the analysis of favipiravir using LC/MS/MS.




Fig. 2 Workflow of Fully Automated Sample Pretreatment

#### Analytical Conditions and Pretreatment of Samples

Favipiravir (PN: C8720\*1), as the target compound, and [<sup>13</sup>C,<sup>15</sup>N]-favipiravir (PN: C8853\*1), as its stable isotope, were purchased from Alsachim, one of the companies of the Shimadzu Group. [<sup>13</sup>C,<sup>15</sup>N]-favipiravir was used as the internal standard (ISTD). Favipiravir was spiked with commercially available human plasma treated with EDTA 2K to prepare calibration curves and QC samples. Analysis was performed using the LC and MS analytical conditions shown in Table 1 and the MRM transition in Table 2. Shim-pack Scepter™ C18-120 (50 mm×2.1 mm I.D., 1.9 µm, P/N: 227-31012-03) was used as the analytical column. Fig. 3 shows the MS chromatograms and structural formulas of the compounds.

A calibration curve was prepared using calibration points at plasma concentrations of 1, 2, 5, 10, 20, 50 and 100 µg/mL for favipiravir (n = 5 for each calibration point). [<sup>13</sup>C,<sup>15</sup>N]-favipiravir (20 µg/mL) solution was prepared using acetonitrile and used as ISTD. The pretreatment steps for the plasma sample spiked with favipiravir are shown in Fig. 2. Samples were automatically prepared through the following series of steps: mixing 20 µL of 75% isopropyl alcohol (IPA), 50 µL of plasma, 10 µL of ISTD and 200 µL of acetonitrile, shaking the mixture, and then filtration of the mixture using a PTFE membrane filter. Finally, the prepared sample was used for analysis.


\*1 Shimadzu GLC and Alsachim Product numbers

Table 1 LC and MS Analytical Conditions



| <LC analytical conditions>               |  | <MS analytical conditions> |            |
|--|--|----------------------------|------------|
| UHPLC                                    | Nexera™ X2   | LC/MS/MS system            | LCMS-8060  |
| Analysis column                          | Shim-pack Scepter C18-120 (50 mm×2.1 mm I.D., 1.9 µm)  | Interface                  | Heated ESI |
| Mobile phase                             | A: 0.05 % Formic acid – water<br>B: 0.05 % Formic acid – acetonitrile                          | MS analysis mode           | MRM (+)    |
| Gradient program (%B)                    | 5 % (0 – 0.30 min) →<br>30 % (0.35 min) →<br>90 % (1.50 – 2.50 min) →<br>5 % (2.60 – 4.00 min) | Heat block temperature     | 400 °C     |
| Flow rate                                | 0.4 mL/min   | DL temperature             | 250 °C     |
| Column oven temperature                  | 40 °C  | Interface temperature      | 300 °C     |
| Injection volume                         | 1.0 µL   | Nebulizing gas flow rate   | 3 L/min    |
| Rinse solution (for external rinse only) | MeOH   | Drying gas flow rate       | 10 L/min   |
|  |  | Heating gas flow rate      | 10 L/min   |

Table 2 MRM Transitions of Favipiravir and [<sup>13</sup>C,<sup>15</sup>N]-Favipiravir

| Compound  | Ion            | Precursor ion (m/z) | Product ion (m/z) |
|---|----------------|---------------------|-------------------|
| Favipiravir                                     | Quantifier ion | 157.70              | 85.10             |
|   | Qualifier ion  | 157.70              | 113.20            |
| [ <sup>13</sup> C, <sup>15</sup> N]-Favipiravir | Quantifier ion | 159.70              | 85.10             |
|   | Qualifier ion  | 159.70              | 113.20            |



Favipiravir  
Formula : C<sub>5</sub>H<sub>4</sub>FN<sub>3</sub>O<sub>2</sub>



[<sup>13</sup>C,<sup>15</sup>N]-Favipiravir  
Formula : C<sub>4</sub><sup>13</sup>CH<sub>4</sub>FN<sub>2</sub><sup>15</sup>NO<sub>2</sub>




Fig. 3 MS Chromatograms and Structural Formulas of Favipiravir (Left) and [<sup>13</sup>C,<sup>15</sup>N]-Favipiravir (Right)

### Preparation of Calibration Curve

The calibration curve prepared using the fully automated sample preparation LC/MS/MS are shown in Table 3. Good linearity with R<sup>2</sup> of 0.9987 was obtained in the set calibration range. The precision of favipiravir (%RSD) was 1.0 % – 5.6 % over the entire concentration range, including the quantitative lower limit. The accuracy of favipiravir ranged between 95 % – 105 %, within acceptance limits of 100 ± 15 %.

Table 3 Linearity, Precision and Accuracy of Favipiravir in Plasma Using Fully Automated LC/MS/MS Obtained from Analysis

| Favipiravir |                      |                       |                |            | Calibration Curve |
|-------------|----------------------|-----------------------|----------------|------------|-------------------|
| ID          | Spiked Conc. (µg/mL) | Average Conc. (µg/mL) | Precision %RSD | Accuracy % |                   |
| Blank       | ---                  | ---                   | ---            | ---        |                   |
| Level 1     | 1                    | 1.04                  | 4.3            | 104        |                   |
| Level 2     | 2                    | 2.00                  | 1.9            | 100        |                   |
| Level 3     | 5                    | 4.77                  | 2.4            | 95         |                   |
| Level 4     | 10                   | 9.87                  | 1.9            | 99         |                   |
| Level 5     | 20                   | 19.7                  | 3.5            | 99         |                   |
| Level 6     | 50                   | 52.7                  | 1.0            | 105        |                   |
| Level 7     | 100                  | 97.9                  | 5.6            | 98         |                   |

### Validation of Analytical System Using QC Samples

Favipiravir was prepared at the following plasma concentrations as QC samples: 3, 50, 90 µg/mL to evaluate its repeatability (Table 4) and between-days reproducibility comparing results obtained over three days (Table 5). Based on the repeatability test result, the precision of favipiravir (%RSD) was 1.6 % – 3.0 %. The accuracy ranged between 94 % – 97 % with acceptance limit of 100 ± 15 %. Based on the test results for between-days reproducibility, the precision of favipiravir (%RSD) was 0.2 % – 7.6 %. The accuracy ranged between 88 % – 99 % within acceptance limit of 100 ± 15 % during QC sample analyses on each of the three days.

Table 4 Repeatability of Favipiravir in Plasma

| Compound    | QC Sample | Spiked Conc. (µg/mL) | Intra-Assay (n=6)     |                |            |
|-------------|-----------|----------------------|-----------------------|----------------|------------|
|             |           |                      | Average Conc. (µg/mL) | Precision %RSD | Accuracy % |
| Favipiravir | Low       | 3                    | 2.90                  | 2.2            | 97         |
|             | Medium    | 50                   | 48.3                  | 3.0            | 97         |
|             | High      | 90                   | 84.8                  | 1.6            | 94         |

Table 5 Between-Days Reproducibility of Favipiravir in Plasma

| Compound    | QC Sample | Spiked Conc. (µg/mL) | Day 1st (n=3)         |                |            | Day 2nd (n=3)         |                |            | Day 3rd (n=3)         |                |            |
|-------------|-----------|----------------------|-----------------------|----------------|------------|-----------------------|----------------|------------|-----------------------|----------------|------------|
|             |           |                      | Average Conc. (µg/mL) | Precision %RSD | Accuracy % | Average Conc. (µg/mL) | Precision %RSD | Accuracy % | Average Conc. (µg/mL) | Precision %RSD | Accuracy % |
| Favipiravir | Low       | 3                    | 2.95                  | 0.2            | 98         | 2.94                  | 7.6            | 98         | 2.64                  | 3.7            | 88         |
|             | Medium    | 50                   | 49.3                  | 1.0            | 99         | 47.1                  | 2.8            | 94         | 46.0                  | 3.6            | 92         |
|             | High      | 90                   | 85.8                  | 1.3            | 95         | 81.3                  | 1.0            | 90         | 79.4                  | 1.7            | 88         |

### Conclusion

Using favipiravir spiked with plasma, a fully automated sample preparation LC/MS/MS analytical system has been developed. The prepared calibration curve showed good linearity. The repeatability and between-days reproducibility of favipiravir were evaluated using QC samples. Good accuracy and reproducibility were obtained.

The product described in this document has not been approved or certified as a medical device under the Pharmaceutical and Medical Device Act of Japan. It cannot be used for the purpose of medical examination and treatment or related procedures.

CLAM, LCMS, Shim-pack Scepter and Nexera are trademarks of Shimadzu Corporation in Japan and/or other countries.

AVIGAN is a registered trademark of Global Response Aid Inc. in the United States.

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with the trademark symbol "TM" or "®".