

# Application News

No. L577

lon Chromatography

## Analysis of Cyanide Ion and Cyanogen Chloride in Mineral Water by Ion Chromatograph-Post Column Method

Mineral water has become an indispensable beverage in our daily lives, as annual per-capita consumption of mineral water in 2019 was 31.7 liters in Japan, 119 liters in the United States, 125.2 liters in Germany, and 147.4 liters in France <sup>(1)</sup>.

Mineral water is classified as a soft drink, and safe products are distribution to the market based on compositional standards in Japan. In a revision of the standards for soft drinks and similar products, Ministry of Health, Labour and Welfare (MHLW) issued notification *Shokuan* 1222 No. 1 "Partial amendment to the Specifications and Standards for Foods, Food Additives, *Etc.* and the Ministerial Ordinance on Milk and Milk Products Concerning Compositional Standards, *Etc.*" dated December 22, 2014 <sup>(2)</sup>, which set the content of cyanide ion and cyanogen chloride in mineral water at no more than a total of 0.01 mg/L, whether sterilized/disinfected or not. The CODEX international standard <sup>(3)</sup> for foods stipulates that the content of cyanide must not exceed 0.07 mg/L.

This article introduces an example of analysis of the cyanide ion and cyanogen chloride in three types of mineral water using a Shimadzu Nexera<sup>™</sup> cyanic analysis system. The analysis complied with the Soft Drinks Test Method in MHLW notification *Shokuan* 1222 No. 4 (hereinafter, notification of enforcement)<sup>(4)</sup>.

A. Morita

#### Analysis Method

In this analysis, detection is conducted at the wavelength of 638 nm using post-column derivatization by the 4-pyridinecarboxylatepyrazolone method after separation of the cyanide ion and cyanogen chloride by the ion exclusion mode, in which a sodium tartrate buffer solution was used as the mobile phase. A twostage reaction is used in this post-column method. The first reaction is a chlorination reaction with a chloramine T solution, and the second reaction is a coloring reaction using a 1-phenyl-3-methyl-5-pyrazolone/4-pyridinecarboxylate solution.

Fig. 1 shows the flow channel diagram of the Nexera cyanic analysis system compliant with the notification of enforcement. Table 1 shows the analysis conditions. To prevent cyanogen chloride from evaporating, the autosampler vials were cooled to  $4 \,^{\circ}$ C in this analysis.

Table 1 Analysis Conditio
---------------------------

<separation> Column Guard column Mobile phase Flow rate Column temp. Injection vol. Vial</separation>	<ul> <li>Shim-pack<sup>™</sup> Amino-Na (100 mm×6.0 mm l.D., 5 μm)<sup>*1</sup></li> <li>Shim-pack CN(G) (10 mm×6.0 mm l.D., 5 μm)<sup>*2</sup></li> <li>10 mmol/L Sodium tartrate buffer</li> <li>0.6 mL/min</li> <li>40 °C</li> <li>100 μL</li> <li>Shimadzu Vials, LC, Polypropylene <sup>*3</sup></li> </ul>			
<post-column reaction=""></post-column>				
First reaction				
Reagent	: Phosphate buffer containing Tg/L Chioramine T			
Flow rate	: 0.5 mL/min			
Second reaction	: 40 C			
Reagent	: 28.7 mmol/L 1-Phenyl-3-Methyl-5-Pyrazolone +96.5 mmol/L Sodium 4-Pyridinecarboxylate			
Flow rate	: 0.5 mL/min			
Reaction temp.	: 100 °C			
Detection	: UV-VIS detector at 638 nm (Lamp: W)			

\*1 : P/N 228-18837-91, \*2 : P/N 228-18837-93, \*3 : P/N 228-31537-91



Fig. 1 Flow Channel Diagram

#### Analysis of Standard Solution

Fig. 2 shows the result when 100  $\mu L$  of a standard solution of the cyanide ion and cyanogen chloride (0.01 mg/L each) was injected. The notification of enforcement requires separate preparation of the cyanide ion standard solution and the cyanogen chloride standard solution.



(0.01 mg/L each)

## Linearity of Calibration Curve

Fig. 3 shows the calibration curves of the cyanide ion and cyanogen chloride standard solutions. The calibration curves were prepared for the concentration range of 0.0025 to 0.025 mg/L, as specified in the notification of enforcement. The calibration curves showed good linearity, as the coefficient of determination ( $r^2$ ) was 0.999 or higher.



(Left: Cyanide Ion, Right: Cyanogen Chloride)

### Repeatability

Using 0.0025 mg/L of the standard solutions of the cyanide ion and cyanogen chloride, 6 continuous analyses were carried out, and the repeatability (peak area) of repeated injections was verified. The relative standard deviations of the cyanide ion and cyanogen chloride were 0.26% and 0.55%, respectively, showing satisfactory repeatability. Thus, it was found that system performance is stable.

### Analysis of Mineral Water Samples

Figs. 4 to 6 show the results when 100  $\mu$ L of three types of mineral water having different hardnesses was injected. The notification of enforcement does not allow addition of a phosphate buffer, as is done in analyses of tap water. The figures also show the results when 0.001 mg/L of the cyanide ion, i.e., 1/10 of the concentration of the standard value, was added to the mineral water. Fig. 7 shows the chromatograms of the 0.001 mg/L standard solution of the cyanide ion and ultrapure water. Table 2 shows the results of a spike-and-recovery test (Average of n = 3 each).

Table 2 Hardness and Recov	very Rate of Mineral Water Samp	oles

Name	Hardness [mg/L]	Recovery rate [%]
Mineral water A	10	102
Mineral water B	304	97
Mineral water C	1468	99

## Conclusion

This article introduced an example of analysis of the cyanide ion and cyanogen chloride in three types of mineral water using a Shimadzu Nexera cyanic analysis system. The recovery rate was confirmed by adding the cyanide ion at a concentration of 1/10 of the standard value to the mineral water samples. Satisfactory recovery rates were obtained even with mineral waters having different hardnesses.

<References>

Shimadzu Corporation www.shimadzu.com/an/

- (1) Statistical data of the Mineral Water Association of Japan, Transition of Per-Capita Consumption of Mineral Water (April 2, 2020) https://minekyo.net/relays/download/5/123/3/444/?file=/files/libs/444/202004021630262016.pdf (In Japanese only)
- Notification of the Ministry of Health, Labour and Welfare Shokuan 1222 No. 1, Partial amendment to the Specifications and Standards for Foods, Food Additives, Etc. and the Ministerial Ordinance on Milk and Milk Products Concerning Compositional Standards, Etc. dated December 22, 2014.
   (3) CODEX Standard for Natural Mineral Waters: CXS 108-1981 (Adp. 1981, Rev. 1997, 2008, Amd, 2001, 2011, 2019)
- (3) CODEX Standard for Natural Mineral Waters: CXS 108-1981 (Adp. 1981, Rev. 1997, 2008, Amd. 2001, 2011, 2019)
  (4) Notification of the Ministry of Health, Labour and Welfare Shokuan 1222 No. 4, Test Methods Related to Partial Amendment to Standards for Soft Drinks, Etc. dated December 22, 2014.

Nexera and Shim-pack are trademarks of Shimadzu Corporation in Japan and/or other countries. Shimadzu Nexera cyanic analysis system is only sold in Japan.



#### For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <a href="http://www.shimadzu.com/about/trademarks/index.html">http://www.shimadzu.com/about/trademarks/index.html</a> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.



CN

1.00 mAL

0.75



(Top: 0.001 mg/L Standard Solution of Cyanide Ion, Bottom: Ultrapure Water)

First Edition: Nov 2020