

Application News

**Gas Chromatograph Mass Spectrometry** 

## No. **M287**

# Analysis of Dioxins in Foods and Feeds Using GC-MS/MS with BEIS

Persistent Organic Pollutants (POPs) in food and animal feed products are analyzed by a variety of methods. In particular, quantitative analysis of dioxins to low concentration levels is necessary because dioxins have high toxicity, even among POPs. Until now, magnetic sector-type (double focusing) GC-MS (GC-HRMS) had been used in analysis of dioxins, as this method provides high quantitation accuracy. However, because the quantitation performance of triple quadrupole GC-MS (hereinafter, GC-MS/MS) has improved remarkably in recent years, the EU issued Commission Regulations (EU) 589/2014 and 644/2017 establishing analysis methods using GC-MS/MS as official test methods with the same positioning as methods using GC-HRMS. For this reason, analysis of dioxins by GC-MS/MS has attracted interest.

In ionization by the El method, Boost Efficiency Ion Source (BEIS) is an ion source that maximizes ion efficiency by optimizing the focal point of the electron beam. BEIS realizes sensitivity approximately four times\* higher than that of conventional ion sources.

In this article, we analyzed dioxins in about 250 samples of approximately 40 types of food and animal feed products using a GC-MS/MS with BEIS. (However, the dioxins were limited to polychlorinated dibenzo-p-dioxin (PCDD) and polychlorinated dibenzofuran (PCDF).) Quantitation performance was evaluated by comparing the analysis results obtained by GC-HRMS and GC-MS/MS. We also evaluated the number of analyses possible while maintaining sensitivity at low concentrations in order to verify the durability of the GC-MS/MS instrument.

M. Takakura

### Samples and Analysis Conditions

All food and feed samples were prepared using automatic pretreatment devices (extraction: SpeedExtractor (BÜCHI Labortechnik AG), sample clean-up: GO-xHT (Miura Co., Ltd.)).

Nonane was used as the final solvent of the samples, and the final solvent amount was 10  $\mu L.$  Standard samples were prepared by mixing DF-ST and DF-LCS (Wellington Laboratories Inc.).

The GC-MS/MS analysis conditions registered in EU Regulations Compliant GC-MS/MS Method Package for Dioxins in Foods were used as the GC-MS/MS analysis conditions. Table 1 shows the detailed conditions.

#### Analysis Results of Standard Sample

As the concentration range of the calibration curve, standard samples were prepared for concentrations from 0.025 pg/ $\mu$ L to 1 pg/ $\mu$ L (double concentration for Octa-PCDD/PCDF).

In the EU Regulations, all compounds must satisfy the two criteria shown below (partially excerpted from EU 589/2014 and 644/2017) at the LOQ (limit of quantitation). This study confirmed that both criteria can be satisfied at the lowest point of the calibration curve. The results are shown in Table 2 on the following page.

#### <u>1. S/N ratio</u>

The concentration of an analyte in the extract of a sample which produces an instrumental response at two different ions to be monitored with a S/N (signal/noise) ratio of 3:1 for the less intensive raw data signal.

or, if for technical reasons the signal-to-noise calculation does not provide reliable results,

#### 2. Lowest concentration point on a calibration curve

The lowest concentration point on a calibration curve that gives an acceptable ( $\leq$ 30%) and consistent (measured at least at the start and at the end of an analytical series of samples) deviation to the average relative response factor calculated for all points on the calibration curve in each series of samples.

		5 Analysis Conditions			
Instrument composition		Detailed analysis conditions (GC)			
mple preparation (extraction): SpeedExtractor (BÜCHI Labortechnik AG)		Inlet liner:	Topaz Single Gooseneck Inlet Liner, w/ Wool (Restek Corp., P/N: 23336)		
Sample preparation (clean-up): Auto sampler:	GO-xHT (Miura Co., Ltd.) AOC-20i/s	Column:	SH-Rxi <sup>™</sup> -5Sil MS (60 m, 0.25 mm l.D., 0.25 μm), (SHIMADZU, P/N: 227-36036-02)		
GC-MS/MS:	GCMS-TQ™ 8050	Injection Mode:	Splitless		
Software:	GCMSsolution <sup>™</sup> Ver. 4.50SP1	Sampling Time:	1.00 min		
	LabSolutions Insight <sup>™</sup> Ver. 3.6	Injection Temp.:	280 °C		
	EU Regulation Compliant GC-MS/MS Method	Column Oven Temp.:	150 °C (1 min)→(20 °C/min)→220 °C→(2 °C/min) →260 °C (3 min)→(5 °C/min)→320 °C (3.5 min)		
	Package for Dioxins in Foods	High Pressure Injection	n: 450 kPa (1.5 min)		
	5	Flow Control Mode:	Linear Velocity (45.6 cm/sec.)		
	6.001()	<ul> <li>Purge Flow:</li> </ul>	20 mL/min		
Detailed analysis conditions (AOC-20i/s)		<ul> <li>Carrier Gas:</li> </ul>	Helium		
# of Rinses with Solvent (Pre-run					
# of Rinses with Solvent (Post-run): 3					
# of Rinses with Sample:	0	Analysis conditions (M	,		
Plunger Speed (Suction):	Low	Ion Source Temp.:	230 °C		
Viscosity Comp. Time:	0.2 sec.	Interface Temp.:	300 °C		
Plunger Speed (Injection):	High	Detector Voltage:	1.8 kV (Absolute)		
Syringe Insertion Speed:	High	Loop time:	0.8 sec. (for native compounds)		
Pumping Times:	5		0.2 sec. (for labeled compounds)		
Inj. Port Dwell Time:	0.3 sec.	Transition:	Conditions of Method Package		
Terminal Air Gap:	No	lon source:	BEIS		
Plunger Washing Speed:	High	Emission current:	150 μΑ		
Washing Volume:	6 µL		•		
Injection Volume:	2 µL				

Table 1 GC-MS/MS Analysis Conditions

\* The increase in sensitivity is dependent on the compound.

	Table 2 Evaluation Results of EOQ in Analysis of Standard Samples									
I.D.	Compound name	Average RRF	RRF (Level 1)	RRFDev (%) (Level 1)	S/N ratio (Level 1)	LOQ				
1	2,3,7,8-Tetrachlorodibenzo-p-dioxin	1.07	1.15	8.10	552	0.025 pg/μL				
2	1,2,3,7,8-Pentachlorodibenzo-p-dioxin	1.09	0.97	10.56	411	0.025 pg/μL				
3	1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	1.14	1.39	22.26	269	0.025 pg/μL				
4	1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin	0.95	0.92	2.72	254	0.025 pg/μL				
5	1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin	1.03	1.25	21.44	260	0.025 pg/μL				
6	1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin	0.92	0.82	11.46	421	0.025 pg/μL				
7	Octachlorodibenzo-p-dioxin	1.19	1.04	12.21	915	0.050 pg/µL				
8	2,3,7,8-Tetrachlorodibenzofuran	1.10	1.05	4.66	793	0.025 pg/μL				
9	1,2,3,7,8-Pentachlorodibenzofuran	1.04	1.00	3.23	483	0.025 pg/μL				
10	2,3,4,7,8-Pentachlorodibenzofuran	0.97	0.89	7.59	474	0.025 pg/μL				
11	1,2,3,4,7,8-Hexachlorodibenzofuran	1.03	0.82	20.72	447	0.025 pg/μL				
12	1,2,3,6,7,8-Hexachlorodibenzofuran	1.09	1.36	24.62	446	0.025 pg/μL				
13	2,3,4,6,7,8-Hexachlorodibenzofuran	1.09	1.39	27.83	286	0.025 pg/μL				
14	1,2,3,7,8,9-Hexachlorodibenzofuran	1.06	1.23	16.10	315	0.025 pg/μL				
15	1,2,3,4,6,7,8-Heptachlorodibenzofuran	1.17	1.05	10.37	705	0.025 pg/μL				
16	1,2,3,4,7,8,9-Heptachlorodibenzofuran	1.02	0.97	4.97	650	0.025 pg/µL				
17	Octachlorodibenzofuran	1.00	0.84	15.80	689	0.050 pg/μL				




Fig. 1 Chromatograms of Representative Compounds in Analysis of Standard Samples (Concentration: 0.025 pg/µL)

#### Evaluation of Sensitivity in Analysis of Actual Samples

Sensitivity at low concentration levels in analysis of actual samples was verified. Fig. 2 on the next page shows the chromatograms of representative compounds. Satisfactory sensitivity near the limit of quantitation (LOQ) was also obtained in analysis of the actual samples.

### Evaluation of Quantitation Accuracy in Analysis of Actual Samples


More than 250 samples of approximately 40 kinds of food and animal feed products were analyzed using GC-MS/MS. The quantitation accuracy of GC-MS/MS was evaluated by analyzing the same GC-MS/MS samples by GC-HRMS and comparing the results. The results were compared by converting the quantitation values of each sample to TEQ (Toxicity Equivalent Quantity). The results are shown in Fig. 3 on the next page.


In the figure on the left in Fig. 3, the horizontal and vertical axes are shown by linear scales. In the figure on the right, logarithmic scales are used to enable detailed confirmation of the results of samples with small quantitation values. Both figures show the quantitation values by GC-HRMS on the horizontal axis, and those by GC-MS/MS on the vertical axis. In both figures, when a correlation exists between the GC-HRMS and GC-MS/MS values, the values are close to a line with a slope of 1 (blue broken line).


In the figure on the left, the quantitation values by GC-MS/MS and GC-HRMS were similar in all samples in which quantitation values were detected at the level of 1 ng/kgTEQ or more.

The figure on the right shows lines (orange broken lines) where the ratio of the quantitation values by GC-MS/MS and GC-HRMS were 50% or 200% as a standard for difference between the quantitation value obtained by the two methods. In samples in which the quantitation value by GC-HRMS was 0.1 ng/kgTEQ or less, some scattered results in which the quantitation value ratio fell outside the range of 50% to 200% were seen.

The sample with the smallest maximum permissible level (ML) was pork fat, with ML of 1.0 ng/kgTEQ. For this reason, a large difference in the quantitation values (i.e., a quantitation value ratio outside the 50% to 200% range) is possible at concentration levels at least 10 times lower than ML. However, no significant difference could be seen in the quantitation performance of GC-MS/MS and GC-HRMS at the concentration level required in analyses.







Yellow Pigment
 other

Fig. 3 Comparison of Quantitation Values of GC-MS/MS and GC-HRMS (Left: Linear Scale, Right: Logarithmic Scale)

### Evaluation of Durability in Analysis of Actual Samples

As an evaluation of durability in analysis of dioxins in food products, actual samples and standard samples (concentration: 0.05 pg/ $\mu$ L) were analyzed alternately, and the number of analyses possible while maintaining sensitivity was evaluated based on the transition of sensitivity for low concentration standard samples. A total of more than 500 analyses of standard samples and actual samples were carried out. Fig. 4 shows the results.

In Fig. 4, the horizontal axis shows the number of analyses and the vertical axis shows the peak area at each analysis number. No large decrease in sensitivity occurred after more than 500 analyses. Next, Table 3 shows the average peak area and repeatability from the 1<sup>st</sup> to the 530<sup>th</sup> analysis. Repeatability was less than 20%RSD for all compounds, indicating that sensitivity could be maintained through the entire test.

## Conclusion

In this experiment, dioxins in more than 250 samples of approximately 40 types of food and animal feed products were analyzed by GC-MS/MS using BEIS, and the quantitation performance of GC-MS/MS was evaluated by comparing the analysis results by GC-MS/MS and GC-HRMS. As a result of this evaluation, no difference was found in the quantitation performance of GC-MS/MS and GC-HRMS at the concentration level necessary in analyses.

Durability in analysis of dioxins in actual samples was also evaluated, and no decrease in sensitivity at low concentration levels occurred after more than 500 analyses.

Based on these results, we found that BEIS has high sensitivity for dioxins and also enables analysis with high durability.

#### Acknowledgement

We wish to take this opportunity to express our sincere thanks to the Laboratoire d'Etude des Résidus et Contaminants dans les Aliments for cooperation in providing samples, guidance, and other assistance in the preparation of this paper.




Fig. 4 Transition of Peak Area in Repeated Analysis

Table 3 Average Peak Area and Repeatability for Standard Samples (Concentration: 0.05 pg/µL) in Durability Test

I.D.	Compound name	Average peak area	SD	%RSD(n = 17)
1	2,3,7,8-Tetrachlorodibenzo-p-dioxin	73596	8321	11.31
2	1,2,3,7,8-Pentachlorodibenzo-p-dioxin	60713	8803	14.50
3	1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	55956	7025	12.55
4	1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin	58167	8034	13.81
5	1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin	56035	9095	16.23
6	1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin	51663	7452	14.43
7	Octachlorodibenzo-p-dioxin	84614	10814	12.78
8	2,3,7,8-Tetrachlorodibenzofuran	105930	11847	11.18
9	1,2,3,7,8-Pentachlorodibenzofuran	80339	9592	11.94
10	2,3,4,7,8-Pentachlorodibenzofuran	88317	10631	12.04
11	1,2,3,4,7,8-Hexachlorodibenzofuran	67814	11761	17.34
12	1,2,3,6,7,8-Hexachlorodibenzofuran	74759	9636	12.89
13	2,3,4,6,7,8-Hexachlorodibenzofuran	75794	9605	12.67
14	1,2,3,7,8,9-Hexachlorodibenzofuran	67878	6056	8.92
15	1,2,3,4,6,7,8-Heptachlorodibenzofuran	67665	10199	15.07
16	1,2,3,4,7,8,9-Heptachlorodibenzofuran	62914	9356	14.87
17	Octachlorodibenzofuran	103483	13911	13.44

GCMS-TQ, GCMSsolution and LabSolutions Insight are trademarks of Shimadzu Corporation in Japan and/or other countries.

Topaz and Rxi are either trademarks or registered trademarks of Restek Corporation in the United States and/or other countries.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®"

First Edition: Jan. 2020



Shimadzu Corporation www.shimadzu.com/an/

#### For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See <a href="http://www.shimadzu.com/about/trademarks/index.html">http://www.shimadzu.com/about/trademarks/index.html</a> for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.