

Application News

Gas Chromatography

No. **G317**

Analysis of Aroma Components in Beer Using HS-20 and Nexis™ GC-2030

Beer contains an extremely large number of aromatic compounds, or aroma components. In this article, nine aroma components, which have a large influence on the flavor of beer, including alcohols and esters, were analyzed in five types of commercial beer. A multivariate analysis (principal component analysis, hierarchical clustering analysis) of those analysis results was also conducted, and the beers were compared.

Y. Takemori

Analysis Method

A standard solution for the calibration curve was prepared by diluting a mixed stock solution of aroma components with 4% ethanol. 10 g of this solution was introduced into a headspace vial, an internal standard was added, and the vial was sealed, after which a headspace analysis was conducted and the calibration curve was prepared. Real samples were prepared in the same manner as described above using 10 g of non-deaerated beer and the internal standard and then sealed in a headspace vial. A headspace analysis was conducted, and quantitation was done using the prepared calibration curve.



Appearance of Nexis[™] GC-2030 + HS-20

Analysis Conditions

Table 1 shows the device composition and analysis conditions.

Table 1 Device Composition and Analysis Conditions

Headspace sampler (HS-20)		Gas chromatograph (Nexis GC-2030 / FID-2030)		
Mode_	: Loop	Injection Mode	: Split	
Oven Temperature	: 40 ℃	Split Ratio	: 1:5	
Sample Line Temperature	: 70 ℃	Carrier Gas	: He	
Transfer Line Temperature	: 80 °C	Carrier Gas Control	: Column flowrate (5 mL/min)	
Vial Pressure	: 150 kPa	Column	: SH-Rtx™-Wax	
Vial Heat-retention Time	: 45 min		(30 m×0.53 mm l.D., 1 μm)	
Vial Pressurization Time	: 1 min	Column Temp	: 40 °C (5 min) - 10 °C /min – 190 °C (5 min)	
Vial Pressurization Equilibrating Time	: 0.1 min	Detector Temp	: 200 °C	
Loading Time	: 0.5 min	Detector Gas	: H ₂ 32 mL/min	
Loading Pressurization Time	: 0.1 min		Make-up (He) 24 mL/min	
Injection Time	: 0.5 min		Air 200 mL/min	
Needle Flush Time	: 5 min			

Analysis Results of 5 Commercial Beers

Fig. 1 shows a representative chromatogram, and Table 2 shows the quantitation results (average of n=3) of the aroma components of the five commercial beer samples.




Fig. 1 Chromatogram of Beer

Table 2 Quantitation Results for Aroma Components (Unit: mg/L)

Compound	Beer A	Beer B	Beer C	Beer D	Beer E
Ethyl acetate	21.065	32.552	12.584	30.386	13.825
n-propyl alcohol	11.342	11.217	11.985	23.662	9.360
Isobutyl alcohol	10.490	10.941	14.126	0.190	9.198
Isoamyl acetate	2.188	2.725	0.936	0.459	1.120
Isoamyl alcohol	62.552	52.284	59.169	59.692	54.575
Ethyl caproate	0.067	0.091	0.052	0.088	0.094
Ethyl caprylate	0.114	0.159	0.067	0.124	0.137
Phenylethyl acetate	0.350	0.808	0.125	Not detected	0.231
Phenylethyl alcohol	16.891	21.655	17.025	13.420	16.583

Multivariate Analysis Results and Study

Using eight components (excluding phenylethyl acetate) that were detected in all the beers in this analysis (n=3), a principal component analysis and hierarchical clustering analysis were carried out with the multivariate analysis software SIMCA® 15 (Infocom Corporation). Fig. 2 shows a score plot of the results. The five beer samples were clearly separated on the score plot, and the differences among the beers could be visualized. Fig. 3 shows a dendrogram. The degree of similarity could be visualized from the dendrogram.

The score plot and dendrogram show that Beer D differs greatly from the other four beers, suggesting that the concentration ratio of the eight components is different in Beer D from the other beers. Moreover, only Beer D was a top fermentation ale beer, and the other beers were bottom fermentation lager beers.




Fig. 2 Score Plot

Fig. 3 Dendrogram

Conclusion

Nine aroma components in beer could be analyzed with high sensitivity, and differences in the component concentrations among the beers could be recognized. Simple visualization of the comparison of the beers was possible by a multivariate analysis of those analysis results.

<Reference>

Brewers Association of Japan, BCOJ Beer Analysis Methods (2013)

Nexis is a trademark of Shimadzu Corporation in Japan and/or other countries.

Rtx is a trademark or registered trademark of Restek Corporation in the United States and/or other countries.

SIMCA is a registered trademark of Sartorius Stedim Data Analytics AB.

First Edition: Dec. 2019



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Shimadzu disclaims any proprietary interest in trademarks and trade names used in this publication other than its own. See http://www.shimadzu.com/about/trademarks/index.html for details.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.