Application News

No. SSK_GCMS_2101

Gas Chromatograph Mass Spectrometer

GCMS-O를 이용한 식품 중 향 물질의 정성분석 (Qualitative Analysis of Flavor compounds using Gas Chromatography Mass Spectrometer- Olfactometry in Food)

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- 1) Shimadzu Technical Report, New Approach for Off-Flavor-Complaint Product Analysis Using GC/MS Off-Flavor Analyzer (C146-E293)

SHIMADZU Scientific Korea Corp.
www.shimadzu.co.kr

For Research Use Only. Not for use in diagnostic procedures. Not available in the USA, Canada, and China. This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Copyright © 2018 SHIMADZU group. All rights reserved.