Application News

No. SSK_ICPMS_2102

Pharmaceutical Analysis/ ICPMS-2030

USP에 따른 MCT의 원소불순물 분석 : ICPMS-2030 (Analysis of Elemental Impurities in Medium-chain Triglycerides According to USP)

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RF power	1.20 kW
Sampling depth	5.0 mm
Plasma gas flow	8.0 L/min
Auxiliary gas flow	1.10 L/min
Carrier gas flow	0.70 L/min
Cell gas (He) flow	6.0 mL/min
Torch type	Mini torch
Sampling & Skimmer cone	Nickel
Quantification method	Internal Standard Correction Method

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□ Reference

1. Medium-Chain Triglycerides, USP43-NF38, 2020.11.01
2. Triglycerides, Medium-Chain, 01/2020:0868, European Pharmacopoeia (Ph. Eur.) 10th Edition, 2019. 07
3. USP43-NF38 <232> ELEMENTAL IMPURITIES—LIMITS, <233> ELEMENTAL IMPURITIES—PROCEDURES, 2020.11.01
4. ICH(The International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use) Q3D(R1) Guideline for Elemental Impurities, Step 4, 2019.03.22

Shimadzu Corporation
www.shimadzu.com/an/

SHIMADZU Scientific Instruments Inc.
www.shimadzu.co.kr

For Research Use Only, Not for use in diagnostic procedures. |

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.