Application News

No. AD-0235K

Pharmaceutical / HS-20 Nexis GC-2030

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	n-부탄올 회수율 (%)			
	Run 1	Run 2	Run 3	평균
Sample 1	99.5	101.2	100.4	100.4
Sample 2	99.2	99.0	100.6	99.6
Sample 3	99.3	99.2	99.3	99.3

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- 1) <https://www.cdc.gov/coronavirus/2019-ncov/hcp/hand-hygiene.html>
- 2) 식품의약품안전처 고시 의약외품 표준제조기준 제 2018-14호
- 3) Shimadzu Application News No. G333, Alcohol Determination of Sanitizer Gel in accordance with USP<611>.
- 4) Lee. J., et al, hand Sanitizers: A Review on Formulation Aspects, Adverse Affects and Regulation, Int J Environ Res Public Health, 2020 May 11; 17(9):3326.
- 5) US EPAMETHOD 5021A, VOLATILE ORGANIC COMPOUNDS IN VARIOUS SAMPLE MATRICES USING EQUILIBRIUM HEADSPACE ANALYSIS, 2014.
- 6) Ettre, L.S., Kolb, B., (2006), STATIC HEADSPACE-GAS CHROMATOGRAPHY. Canada: John Wiley & Sons.

SHIMADZU Scientific Korea Corp.
www.shimadzu.co.kr

For Research Use Only. Not for use in diagnostic procedures. Not available in the USA, Canada, and China. This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "@". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "@". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Copyright © 2018 SHIMADZU group. All rights reserved.