Application News

No. 06-SAIP-GC-034-ENK

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Similarity Search Results

Report View Compound Info Process Help

#	Similarity	Register	Ret. Index	Compound Name	Mol Wt	Formula	Library Name
1	92	<input checked="" type="checkbox"/>	1260	Phosphoric acid-TMS : phosphoric acid	314	C9H27O4P3	OA_TMS_BPX5_23min_V3.lib
2	45	<input type="checkbox"/>	1256	Glycerol-3TMS : propane-1,2,3-triol	308	C12H32O3S3	OA_TMS_BPX5_23min_V3.lib
3	40	<input type="checkbox"/>	1269	Ethylmalonic acid-2TMS : 2-ethylpropanedioic acid	276	C11H24O4S2	OA_TMS_BPX5_23min_V3.lib
4	34	<input type="checkbox"/>	1260	Octanoic acid-TMS : octanoic acid	216	C11H24O2S	OA_TMS_BPX5_23min_V3.lib
5	32	<input type="checkbox"/>	1271	2-Ketooisocaproic acid-oxime-2TMS : 4-methyl-2	289	C12H27NO3S2	OA_TMS_BPX5_23min_V3.lib
6	26	<input type="checkbox"/>	1255	Benzoic acid-TMS : benzoic acid	194	C10H14O2S	OA_TMS_BPX5_23min_V3.lib

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ID #	Compound Name	ID #	Compound Name	ID #	Compound Name	ID #	Compound Name
1	Lactic acid	17	Phosphoric acid	33	Creatinine	49	Homovanillic acid
2	Glycolic acid	18	Ethylmalonic acid	34	2-Hydroxyphenylacetic acid	50	Citric acid
3	Pyruvic acid	19	Succinic acid	35	3-Hydroxyglutaric acid	51	Hippuric acid
4	2-Keto-isovaleric acid	20	Methylsuccinic acid	36	3-Phenyllactic acid	52	Dopamine
5	2-Hydroxybutyric acid	21	Uracil	37	3-Hydroxy-3-methylglutaric acid	53	Homogentisic acid
6	Oxalic acid	22	Fumaric acid	38	Phenylpyruvic acid	54	Methylcitric acid
7	4-Cresol	23	Thymine	39	2-Ketoglutaric acid	55	3-(3-Hydroxyphenyl)-3-hydroxypropionic acid
8	2-Hydroxyisovaleric acid	24	Glutaric acid	40	Tartaric acid	56	Vanilmandelic acid
9	3-Hydroxybutyric acid	25	3-Methylglutaconic acid	41	4-Hydroxybenzoic acid	57	4-Hydroxyphenyllactic acid
10	3-Methyl-2-oxovaleric acid	26	3-Methylglutaric acid	42	N-Acetylaspartic acid	58	Pyridoxine
11	Malonic acid	27	Mandelic acid	43	Arabinose	59	Sebacic acid
12	Acetoacetic acid	28	Citramalic acid	44	Suberic acid	60	Ascorbic acid
13	Methylmalonic acid	29	Malic acid	45	Quinolinic acid	61	2-Hydroxyhippuric acid
14	2-Ketoisocaproic acid	30	Adipic acid	46	Aconitic acid	62	Pantothenic acid
15	2-Hydroxyisocaproic acid	31	2-Phenyllactic acid	47	Orotic acid	63	Kynurenic acid
16	4-Hydroxybutyric acid	32	5-Hydroxymethyl-2-furoic acid	48	Tricarballic acid	64	Margaric acid
						65	3-Indoleacetic acid

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■ References

1. Naraian, R., & Kumari, S. (2017). Microbial production of organic acids. *Microbial Functional Foods and Nutraceuticals*, 93.
2. Suiryanrayna, M. V., & Ramana, J. V. (2015). A review of the effects of dietary organic acids fed to swine. *Journal of animal science and biotechnology*, 6(1), 1-11.
3. Yan, R., Jiang, H., Gu, S., Feng, N., Zhang, N., Lv, L., & Liu, F. (2020). Fecal metabolites were altered, identified as biomarkers and correlated with disease activity in patients with systemic lupus erythematosus in a GC-MS-based metabolomics study. *Frontiers in immunology*, 2138.

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For Research Use Only, Not for use in diagnostic procedures. Not available in the USA, Canada, and China. This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Copyright © 2022 SHIMADZU group. All rights reserved.