Application News

No. GCMS-2202-K

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1. Ishimura, T.; Iwai, I.; Matsui, K.; Mattonai, M.; Watanabe, A.; Robberson, W.; Cook, A.; Allen, H.L.; Pipkin, W.; Teramae, N.; Ohtani, H.; Watanabe, C. Qualitative and quantitative analysis of mixture of microplastics in the presence of calcium carbonate by Pyrolysis-GC/MS. *Journal of Analytical and Applied Pyrolysis*, 2021 (157): 105188.
2. WK67788 - New Test Method for Identification of Polymer Type and Quantity (Mass) Measurement of Microplastic Particles and Fibers in Waters with High-to-Low Suspended Solids Using Pyrolysis-Gas Chromatography/Mass Spectrometry: Py-GC/MS.
3. Joint ISO/TC 147/SC 2 - ISO/TC 61/SC 14 WG: Plastics (including microplastics) in waters and related matrices.
4. Pipkin, W.; Belganeh, R.; Robberson, W.; Allen, H.L.; Cook, A.; Watanabe, A. Identification of Microplastics in Environmental Monitoring Using Pyrolysis-GC-MS Analysis. *LCGC North America*, 2021 April 01; 39(4):179-186.
5. How to Kickstart Your Micro- and Nanoplastic Py-GC/MS Analysis. Shimadzu Application News. GCMS-2201.

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Part Number	Item Name	Item Description
UAMP-K01	UAMP Column Bracket	UA Precolumn 50: Ultra Alloy-50 (2 m x 0.25 mm ID x 1.0 μm) and separation column: Ultra Alloy-5 (30 m x 0.25 mm ID x 0.5 μm)
PY1-3346	Packed GC glass insert	Packed GC glass insert for Shimadzu
220-90906-00	Eco-Cup LF	PY-2020 Eco-Cup LF
220-94824-18	Quartz Pyrolysis Tube	Pyrolysis tube, 120MM PY
220-90917-00	Graphite Vespel ferrule	Graphite Vespel Ferrule for Py-2020 (pk of 5)
220-94824-14	Needle Set	Syringe, deactivated needle (PY) (3/pk)
PY1-4940	MP calibration standard set	MP calibration standard consisting of 12 polymers
MS402280	Vent-free GC/MS Adapter	Vent-free GC/MS Adapter N-50
220-94792-03	Quartz wool	Phthalate-free quartz wool
220-94792-00	Septa	Septa for Pyrolysis
221-49662-91	O-ring	O-ring for sealing glass liners when using a Pyrolyzer

No. GCMS-2202-K

Shimadzu Corporation
www.shimadzu.com/an/
 Shimadzu Scientific Korea
www.shimadzu.co.kr

For Research Use Only. Not for use in diagnostic procedures. Not available in the USA, Canada, and China.
 This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Copyright © 2023 SHIMADZU group. All rights reserved.