Application News

J107K

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- 1) Official Specifications Related to Typical Fertilizer Based on the Fertilizer Control Act (Ministry of Agriculture, Forestry and Fisheries Notification No. 284, February 22, 1986, Revised on August 8, 2012 by the Ministry of Agriculture, Forestry and Fisheries Notification No. 1985, Enforced from September 7, 2012)
- 2) Enforcement Regulations Regarding the Soil Contamination Countermeasures Law (Ministry of the Environment Ordinance No. 29, December 26, 2002)
- 3) Ordinance on Test Method for Water Quality of Sewage (Ministry of Health and Welfare/Ministry of Construction Ordinance No. 1, December 17, 1962, Revised on May 23, 2012 by the Ministry of Land, Infrastructure, Transport and Tourism/Ministry of the Environment Ordinance No. 2)
- 4) Ordinance for Determination of Effluent Standards (Prime Minister's Office Ordinance No. 35, June 21, 1971, Revised on September 4, 2013 by the Ministry of the Environment Ordinance No. 20)
- 5) JIS K0102-2013 (Industrial Wastewater Test Method)

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For Research Use Only. Not for use in diagnostic procedures. Not available in the USA, Canada, and China. This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Copyright © 2023 SHIMADZU group. All rights reserved.