Application News

No. 01-00662-K

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1. Simple Analysis of Impurities in Oligonucleotide Therapeutics Using a Single Quadrupole Mass Spectrometer, Application News No.01-00656-EN
2. An Oligonucleotide Impurity Analysis Workflow Using LabSolutions Insight™ Biologics Software, Application News No.01-00595A-EN

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	Enable	Detect	m/z	Charge	Actual	Mass	Weight	Intensity
1	<input checked="" type="checkbox"/>	Auto	1261.8500	17	16.9974	21468.57376	0.058	116658
2	<input checked="" type="checkbox"/>	Auto	1191.3500	18	18.0024	21462.43104	0.077	174175
3	<input checked="" type="checkbox"/>	Auto	1128.6500	19	19.0016	21463.48832	0.232	583305
4	<input checked="" type="checkbox"/>	Auto	1072.3500	20	19.9983	21467.14560	0.262	728457
5	<input checked="" type="checkbox"/>	Auto	1021.1500	21	21.0000	21465.30288	0.335	1025763
6	<input checked="" type="checkbox"/>	Auto	974.7500	22	21.9986	21466.66016	0.036	119888

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Shimadzu Corporation
www.shimadzu.com/an/

Shimadzu Scientific Korea
www.shimadzu.co.kr

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. See <http://www.shimadzu.com/about/trademarks/index.html> for details.

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.