

Are you sticking with drinking Bottled Water? Assessment of PFAS content in commercial samples

Eberhardt Kuhn, Brahm Prakash, Ruth Marfil-Vega

Shimadzu Scientific Instruments, Inc. 7102 Riverwood Drive, Columbia, MD 21046

Introduction

The presence of Per- and Polyfluorinated Alkyl Substances (PFAS) in drinking water is being thoroughly studied due to the persistence of these compounds in the environment and their potential health effects. However, there is limited knowledge about the occurrence of these chemicals in bottled water, despite the increasing concerns about PFAS in the food supply. This poster shows results from a fast and simple direct injection method similar to draft EPA method 8237, using the Shimadzu LCMS-8050 to analyze seven commercially available samples of bottled water for 24 PFAS. The results demonstrate that the instrument's performance exceeds the Quality Assurance requirements in FDA draft method C-010.01 for other matrices, including milk (which is the most similar to water), as well as the limits established by the EPA for the analysis of PFAS drinking water.

While the origin of the water itself maybe the source of PFAS in bottled water, we also wanted to investigate the importance of the type of materials. Migration of PFAS from Food Contact Materials (FCM) is known to occur in all kinds of food containers. In this study, we procured bottled water in several different bottle materials, as well as two types of water source. These included spring and purified water, and bottles made from 5 different kinds of container materials: plastic (virgin and recycled), glass, metal, and cardboard. Preliminary results indicate that observed PFAS levels seem to depend on both the bottle material and the water source.

Instrument Operating Conditions

Table 1: Chromatography and mass spectrometer conditions (for more details, please see Poster 626-18P)

Parameter		Value					
LCMS	Shimadzu LCMS-8050						
Analytical Column	Restek Raptor C18 2.1 mm ID. x 150 mm L., 2.7 µm)						
·	Part No 9304A62						
Solvent Delay Column	Restek PFAS Delay Column (2.1 mm ID. x 50 mm L)						
	Part No. 27854						
Column Oven Temperature	40 ° C						
Injection Volume	30 μL						
Mobile Phase	A: 20 mmol Ammonium Acetate in 5 % (v/v) Acetonitrile in reagent water						
	B: 10 mmol Ammoni	um Acetate in 95 % (v/v) Acetonitrile in reagent water					
Gradient Flow rate	0.3 mL/ Min						
Gradient	Time (minutes)	<u>% B</u>					
	0	0					
	1	20					
	6	50					
	14	100					
	17	100					
	18	0					
	21	0					
Run time	21 minutes						
Nebulizing gas flow	5 L/min						
Heating gas flow	15 L /Min						
Interface temperature	300 °C						
Desolvation Line temperature	100 °C						
Heat Block temperature	200 °C						
Drying gas flow	5 L /min						
Acquisition cycle time	21 min						
Total MRMs	66						

Experimental

We analyzed 24 target PFAS compounds and 19 surrogates in various types of water. The analysis of PFAS was performed using a Shimadzu Nexera X2 SIL-30AC autosampler and a LCMS-8050 triple quadrupole mass spectrometer. An injection volume of 30 µL was used in this study. A detailed description of the LC/MS/MS parameters is included in Table 1.

Chromatography was adjusted to obtain maximum resolution between peaks in the shortest time possible with minimum coelution of isomers. The total run time of 21 minutes includes a final wash out with concentrated acetonitrile to flush the column, remove background residuals contaminants and restore column performance before starting the next run. The method could easily be modified to include isotopic dilution or internal calibration if needed for quantifying the concentrations.

Calibration Standards

Standards available from Wellington Laboratories were used for these studies (Catalog no. PFAC-24PAR and MPFAC-24ES). These standards were then diluted to working standards using 95:5 acetonitrile:water as the diluent. The working standards were used to create a calibration curve ranging from 5-200 ppt with the injection solvent consisting of 50:50 water:methanol with 0.1% acetic acid in order to match the injection solvent for the extracted samples. Filtration was not performed on the calibration standards.

Sample Preparation

Seven types of bottled water as sample matrices were tested using reagent water as the blank. Each sample was diluted 50:50 with MeOH and 0.1% acetic acid, spiked with isotopically labeled surrogates and vortexed for 2 min. The samples were then filtered through 0.2 µm syringe filters and analyzed by LC/MS/MS.

All compound parameters, including precursor ion, product ion, and collision energies, were optimized. There are at least two multiple reaction monitoring (MRM) transitions for most of the analytes. For more details, please see Pittcon 2020 Poster 626-18P.

Results and Discussion

It is known that PFAS can be present in reagents, glassware, pipettes, tubing, degassers and other parts from the LC-MS/MS instruments. PFAS contamination coming from the LC system is eliminated using a delay column placed between the reagents and the sample valve. This separates PFAS in the sample from the PFAS in the LC system. All supplies used to conduct the study were free from PFAS contamination. To monitor the lack of contamination two blanks were injected at the beginning of each batch: system null injection (air injection) and reagent blank (0.1% acetic acid in high purity water:methanol (50:50)). Figure 1 shows the schematic of the delay column set up, and Figure 2 shows the importance of having a delay column and its impact on data quality: the background PFAS contamination elutes at a later retention time than the targeted compounds present in the samples.

Fig. 1: Schematic of Delay Column System




Fig. 2: Comparison of Chromatograms with and without a Delay Column




Table 2: Recoveries and Reproducibility of 80 ppt Standard (n=3)

PFAS Name	Average ppt	% Recovery	%RSD
PFBA	103.6	129.5	21.2
MPFBA	81.9	102.3	2.3
PFPeA	79.3	99.1	2.5
M5PFPeA	83.5	104.4	1.4
4-2 FTS	85.5	106.9	6.4
M2-4-2 FTS	86.3 107.9		9.4
PFHxA	78.1	97.6	3.7
M5PFHxA	80.4	100.6	2.8
PFBS	79.6	99.4	3.4
M3PFBS	79.9	99.8	1.3
PFHpA	81.4	101.7	3.4
M4PFHpA	82.3	102.9	3.0
PFHxS	79.9	99.8	3.9
PFPeS	76.3	95.4	8.3
6-2 FTS	80.9	101.1	10.8
M2-6-2 FTS	77.4	96.8	23.0
PFOA	76.7	95.9	6.6
M8PFOA	82.7	103.4	4.9
M3PFHxS	77.7	97.1	5.2
PFOS	83.9	104.9	15.6
PFNA	76.4	95.5	2.8
M9PFNA	82.4	103.0	1.7
PFHpS	83.0	103.7	12.3
8-2 FTS	70.1	87.6	14.6
M2-8-2 FTS	74.3	92.8	3.8
N-EtFOSAA	73.1	91.4	2.9
N-MeFOSAA	79.2	99.1	2.7
PFDA	80.0	100.0	4.1
M6PFDA	77.0	96.2	2.9
d3-NMeFOSAA	73.6	92.0	7.0
M8PFOS	79.8	99.7	7.7
d5-NEtFOSAA	81.6	102.0	15.6
PFUnA	79.0	98.7	4.6
M7PFUnA	77.8	97.3	3.8
PFNS	80.6	100.8	19.6
PFDoA	76.6	95.7	4.7
M2PFDoA	74.0	92.6	2.8
PFDS	84.2	105.2	10.0
PFTriA	73.3	91.6	1.9
FOSA		80.2 100.3	
M8FOSA	73.5 91.9		6.5 4.0
PFTreA	76.5	95.6	3.8
M2PFTreA	73.5	91.9	5.0
HFPO-DA	82.4	103.0	4.2
3C-HFPO-DA SUR	79.3	99.1	5.9
ADONA	79.5 79.9	99.9	3.2
9CI-PF3ONS	79.9 77.9	99.9	3.7
11CI-PF3OUdS	77.9 79.4	99.2	9.2
11047130003	13.4	33.2	⋑.∠

Table 3: PFAS Data by Water Bottle Material (n=3)

DEAC Name	blank	alass	a a walle a a wal	ma et al	nlostis 4	nlostic 2	nlostic 2	waanalaad u laatin
PFAS Name	<u>blank</u>	glass	cardboard		plastic 1	plastic 2	plastic 3	recyclced plastic
PFBA	NQ	NQ	NQ	23.8 ± 1.4	NQ	15.3 ± 2.1	104.3 ± 4.2	18.1 ± 2.0
PFPeA	ND	ND	ND	ND	ND	ND	ND	ND
4-2 FTS	ND	NQ	NQ	NQ	NQ	NQ	NQ	NQ
PFHxA	ND	ND	ND	ND	ND	ND	ND	ND
PFBS	NQ	ND	ND	NQ	ND	ND	ND	ND
PFHpA	ND	NQ	NQ	ND	NQ	NQ	ND	ND
PFHxS	ND	ND	ND	ND	ND	ND	ND	ND
PFPeS	ND	ND	ND	ND	ND	ND	ND	ND
6-2 FTS	NQ	NQ	ND	NQ	81.5 ± 3.9	NQ	ND	253.9 ± 26.1
PFOA	ND	NQ	NQ	NQ	NQ	NQ	NQ	NQ
PFOS	ND	ND	ND	ND	ND	ND	ND	ND
PFNA	ND	NQ	NQ	NQ	NQ	NQ	NQ	NQ
PFHpS	ND	ND	ND	ND	ND	ND	ND	ND
8-2 FTS	NQ	ND	ND	ND	ND	ND	ND	ND
N-EtFOSAA	ND	ND	ND	ND	ND	ND	ND	ND
N-MeFOSAA	ND	ND	ND	ND	ND	ND	ND	ND
PFDA	NQ	ND	NQ	NQ	ND	NQ	NQ	NQ
PFUnA	ND	ND	ND	ND	ND	ND	ND	ND
PFNS	ND	ND	ND	ND	ND	ND	ND	ND
PFDoA	NQ	NQ	NQ	ND	NQ	ND	ND	NQ
PFDS	ND	ND	ND	ND	ND	ND	ND	ND
PFTriA	ND	ND	ND	ND	ND	ND	ND	ND
FOSA	ND	ND	ND	NQ	ND	ND	ND	ND
PFTreA	ND	ND	ND	ND	ND	ND	ND	ND
HFPO-DA	ND	ND	ND	NQ	ND	ND	ND	ND
ADONA	ND	ND	ND	ND	ND	ND	ND	ND
9Cl-PF3ONS	ND	ND	ND	ND	ND	ND	ND	ND
11Cl-PF3OUdS	ND	ND	ND	ND	ND	ND	ND	ND

ND = not detected; NQ = not quantitated

Summary and Conclusions

This study evaluated the direct injection analysis of 24 PFASs and 19 mass-labeled surrogates in bottled water using Shimadzu UFMS™ LCMS-8050. The data shows excellent performance of the LCMS-8050 for PFAS analysis in bottled water matrices with minimal sample preparation. Of the seven types of bottled water containers tested, plastic had the highest amount of PFAS present. In particular, the bottle made from recycled plastic showed by far the highest amount of PFAS. Glass and cardboard bottles had no detectable PFAS levels. While this data would suggest that the plastic recycling process introduces additional PFAS into the water bottles, it is too early to draw that conclusion with any certainty. The data also suggests that the source of the water can contribute PFAS to the bottled water. More experiments are necessary to confirm these preliminary results.

References

FDA Makes Available Testing Method for PFAS in Foods

https://www.fda.gov/food/cfsan-constituent-updates/fda-makes-available-testing-method-pfas-foods-and-final-results-recent-surveys

FDA Makes Available Results from Second Round of Testing for PFAS in Foods from the General Food Supply https://www.fda.gov/food/cfsan-constituent-updates/fda-makes-available-results-second-round-testing-pfas-foods-general-food-supply

Environmental Science and Technology. "Polyfluorinated Compounds: Past, Present, and Future" http://www.greensciencepolicy.org/wp-content/uploads/2014/10/Lindstrom-Strynar-and-Libelo-2011.pdf

International Bottled Water Association – Recommended PFAS Levels in bottled water https://www.bottledwater.org/and-polyfluoroalkyl-substances-pfas-and-bottled-water

US EPA. Basic Information about Per- and Polyfluoroalkyl Substances (PFAS) https://www.epa.gov/pfas/basic-information-about-and-polyfluoroalkyl-substances-pfass

Analysis of Per-and Polyfluoroalkyl Substances (PFAS) in Non-Drinking Water Matrices Using the LC-Triple Quad Mass Spectrometer Brahm Prakash, Ruth Marfil-Vega, Christopher Gilles, William Lipps Pittcon 2020 Poster 626-18P Shimadzu Scientific Instruments, Inc Columbia, MD 21046

There are two data points that stand out:

Recoveries of an 80 ppt standard are shown

in Table 2. The data represent an average

Recoveries ranged from 87.6% to 129.5%,

The LOQ was determined at 10 ppt in the

compounds, but only two were found in any

of the samples – PFBA (perfluoro butanoic

acid) and 6-2 FTS (fluorotelomer sulfonate)

The highest levels were found in plastic

bottles. Much of the attention in PFAS

analysis has been on the longer chain

did not detect either one of these in the

effects of PFBA and 6-2 FTS, although

Table 3 shows a summary of the PFAS

residues in each water sample; each

claimed "spring water" as their water

in consumer products.

samples tested. Limited information is

analogs, especially PFOS and PFOA. We

currently available about the potential health

interest on these compounds has increased

in recent years as they are frequently used

sample type was analyzed in triplicate. The

sample "Plastic 2" was labelled as "purified

water" on the bottle, while all other bottles

source. Only two of the samples had no

bottle and the cardboard container. The

other containers had at least one PFAS

above the levels recommended by The

International Bottled Water Association

PFAS, and 10 ppt for total PFAS

exposure limits.

(IBWA). This organization has guidelines

for its members of 5 ppt for any individual

concentration. The EPA has set guidelines of total PFAS concentration below 70 ppt in

drinking water, while the FDA is currently

contamination and will use this data to set

testing many types of foods for PFAS

detectable PFAS concentration – the glass

This demonstrates the accuracy and

reproducibility of the measurements

and %RSD was below 10% for most

Our method screened for 24 PFAS

of 3 individual runs.

compounds.

- 1. by far the highest level of PFAS was measured in the bottle made from recycled plastic (highlighted in yellow) (it was stated on the label that the bottle was made with 50% recycled plastic). Taken at face value, this would imply that the recycling process introduces additional PFAS into the produced plastic bottles.
- 2. The bottle labeled "Plastic 3" contained by far the highest amount of PFBA (highlighted in blue). Since the water source was identified as "Spring Water", it could indicate local PFAS contamination near the source of that water.

However, we only have one data point each so far, so it is too early to come to any meaningful conclusions. The data does warrant further investigation, though, on both the source of the water and the bottle material.