

# Improvement of total analytical workflow by using online SFE-SFC

**Pittcon 2017** 2100-1P

Kenichiro Tanaka<sup>1</sup>; Keiko Matsumoto<sup>2</sup>; Hidetoshi Terada<sup>2</sup>;  
Takato Uchikata<sup>2</sup>; Yasuhiro Funada<sup>2</sup>

<sup>1</sup>Shimadzu Scientific Instruments Inc., Columbia, MD,

<sup>2</sup>Shimadzu Corporation, Kyoto, Japan

## Introduction

Total analysis time is composed of sample preparation, a chromatographic separation of the target compounds and column equilibration. Online SFE-SFC can decrease the total analysis time, because after simply placing a sample in a special extraction vessel, it automatically performs the entire process from extraction of target components to acquisition of data. Furthermore, SFC can be expected to separate isomers in a shorter time than HPLC. Online SFE-SFC is particularly time-saving in the pharmaceutical industry. More than half of low molecular-weight drugs have stereoisomers, and pharmacological activities of each enantiomer are different. Therefore, it is important that the efficacy and safety of compounds are accurately evaluated as enantiomers, especially in pharmaceutical formulations and its related industries. Chiral separation using SFC and HPLC is one of the typical methods for purifying enantiomers from

racemic mixtures. A suitable column and mobile phase for targeted chiral separation have to be evaluated before starting the analysis. To determine the optimized analytical conditions, a large number of candidate conditions have to be examined and this process requires extensive method development. A more prompt and simplistic system for determining the optimized analytical conditions has been needed. A HPLC/SFC switching system meets and exceeds those needs.

We demonstrated high-resolution separations using a HPLC/SFC switching system, which automatically switches analytical modes between HPLC and SFC in a single sequence. Here, we report the process of a high-efficiency method development workflow for chiral compounds by using SFC and HPLC in a single sequence. In addition, an online SFE-SFC method was developed and a plasma sample was analyzed under the optimized conditions.

## Materials and Method

### Sample and column

Two standard chiral compounds (Omeprazole, Warfarin) were analyzed, as shown in Fig. 1. For HPLC / SFC chiral analytical columns, the i CHIRAL-6 series (CHIRALPAK® IA/IB/IC/ID/IE/IF, Daicel Corporation) were used. Column specifications, such as stationary phase and particle size, are shown in Table 1 and Fig. 2.





Fig. 1: Structures of Chiral Compounds


## Improvement of total analytical workflow by using online SFE-SFC

Table 1: Chiral Columns used in this study

Columns name	Stationary Phase	Particle size	Diameter	Length	
CHIRALPAK® IA	Amylose tris (3, 5-dimethylphenylcarbamate)	3 $\mu$ m	SFC	SFC	
CHIRALPAK® IB	Cellulose tris (3, 5-dimethylphenylcarbamate)		3.0 mm	100 mm	
CHIRALPAK® IC	Cellulose tris (3, 5-dichlorophenylcarbamate)		HPLC	4.6mm	50mm
CHIRALPAK® ID	Amylose tris (3-chlorophenylcarbamate)				
CHIRALPAK® IE	Amylose tris (3, 5-dichlorophenylcarbamate)				
CHIRALPAK® IF	Amylose tris (3-chloro-4-methylphenylcarbamate)				



**Amylose Derivatives**



**Cellulose Derivatives**




Fig. 2: Chiral Columns

## System

Fig. 3 shows a flow diagram of the "Nexera UC online-SFE-HPLC/SFC switching system" that was developed in this experiment. This system consists of the online SFE-SFC system "Nexera UC", which combines supercritical fluid extraction and supercritical fluid chromatography, and ultra high-performance liquid chromatography "Nexera X2". The single instrument can be used for both SFC and HPLC by switching pumps for

delivering solvent or CO<sub>2</sub>, and by regulating the backpressure or not. Furthermore, solvent switching valves and column switching valves are assembled into this system, allowing combinations of columns and mobile phases for the scouting analyses to be changed automatically.

It allows for online SFE-SFC analysis without reconnecting flow lines.

# Improvement of total analytical workflow by using online SFE-SFC





Fig. 3: Flow diagram "Nexera UC online-SFE-HPLC/SFC switching system" (for online SFE-SFC analysis)

## Software

Software performs the work of switching between HPLC and SFC. Purging of mobile phases, which is required when switching between SFC and HPLC, is accomplished by simply executing a batch table automatically generated

by the dedicated software: "Nexera UC Method Scouting Solution Ver.2". It also enables multiple mobile phases and columns to be used for method scouting.



Nexera UC Method Scouting Software

Analysis	Visit	Tray Name	Sample Name	Method File	Data File
1	-1	1		di#Base method for SFC analysis.icm	\\Data\Project\1\EquilData\011.icf
2	-1	1		di#WDisplacementMethodforSFC.icm	\\Data\Project\1\EquilData\012.icf
3	-1	0		di#Base method for SFC analysis.icm	\\Data\Project\1\EquilData\013.icf
4	1	1	Blank	di#Base method for SFC analysis.icm	UCRP_MeOH_20_40_014.icf
5	1	1	STD-1	di#Base method for SFC analysis.icm	UCRP_MeOH_20_40_015.icf
6	1	1	STD-2	di#Base method for SFC analysis.icm	UCRP_MeOH_20_40_016.icf
7	4	1	Sample-1	di#Base method for SFC analysis.icm	UCRP_MeOH_20_40_017.icf
8	4	1	Sample-2	di#Base method for SFC analysis.icm	UCRP_MeOH_20_40_018.icf
9	-1	0		di#SwitchingMethodforUHPLC.icm	\\Data\Project\1\EquilData\019.icf
10	-1	1		di#WDisplacementMethodforLC.icm	\\Data\Project\1\EquilData\020.icf
11	-1	1		di#Base method for LC analysis.icm	\\Data\Project\1\EquilData\021.icf
12	-1	0		di#Base method for LC analysis.icm	\\Data\Project\1\EquilData\022.icf
13	1	1	Blank	di#Base method for LC analysis.icm	VPODS_MeOH_Water_20_40_018.icf
14	2	1	STD-1	di#Base method for LC analysis.icm	VPODS_MeOH_Water_20_40_019.icf
15	3	1	STD-2	di#Base method for LC analysis.icm	VPODS_MeOH_Water_20_40_019.icf
16	4	1	Sample-1	di#Base method for LC analysis.icm	VPODS_MeOH_Water_20_40_018.icf

Automatically Generated Batch File

- **Automation of method scouting process**  
12 column switching for HPLC/SFC system  
Automatic examination with multiple gradient patterns
- **Quickly Check Analysis Results**  
The results for all the conditions considered can be automatically output as an Excel file

Switches from UHPLC to SFC  
SFC analysis  
Switches from SFC to UHPLC  
UHPLC analysis

Fig. 4: Method scouting software

## Improvement of total analytical workflow by using online SFE-SFC

### Analytical conditions

HPLC / SFC analytical conditions are shown in Table 2 / Table 3. By using 3 mobile phases and 12 columns (HPLC 6 columns and SFC 6 columns), a total of 36 analytical conditions were automatically examined. Online-SFE-SFC analytical conditions for a rat plasma sample are shown in Table 4.

Table 2: HPLC chiral analytical conditions

No.	Mobile phase (A/B)	Others
1	Hexane/ Ethanol	B Conc.(%) : 20% (Isocratic) Time program : 20 % (0 -6min) - 40 % (6-8 min : wash) - 20 % (8 -12min)
2	Hexane / Isopropyl alcohol	Flow Rate : 2mL/min Column Temp : 40 °C
3	Methyl tertiary butyl ether / Ethanol	Inj. Vol. : 1 uL Detection : PDA@220 nm

Table 3: SFC chiral analytical conditions

No.	Modifier	Others
1	Methanol	Modifier Conc. (%) : 20% (Isocratic) Time program : 20 % (0 -5min) - 40 % (5-7 min : wash) - 20 % (7 -10min)
2	Ethanol	Flow Rate : 3mL/min Column Temp. : 40 °C
3	Acetonitrile Ethanol = 75 / 25 (v/v)	Inj. Vol. : 1 uL BPR Press : 10 MPa Detection : PDA@220 nm

Table 4: Online-SFE-SFC analytical conditions

SFE
Extraction Time : 15 min
Mobile Phase : A: CO <sub>2</sub> B: Ethanol
B Conc. : 0 % (v/v)
Flowrate : 5.0 mL/min
Back Pressure : 10 MPa

SFC
Column : CHIRALPAK® IA (100 mm L. × 3.0 mm I.D., 3 μm)
Mobile Phase : A: CO <sub>2</sub> B: Ethanol
Time Program : 0 (0min) - 20 % (1 to 4 min) → 40 % (4 to 9min)
Flowrate : 3.0 mL/min
Column Temp. : 40 °C
Back Pressure : A) 10 MPa, B) 40 MPa
Detector : PDA@220 nm

## Results

### Column and modifier scouting

All chromatograms of warfarin are shown in Fig. 5. In some conditions, isomers were successfully separated. A comparison of SFC/HPLC chromatograms using 36 analytical conditions indicated that SFC parameters provided better separation for warfarin.

## Improvement of total analytical workflow by using online SFE-SFC





Fig. 5: Chromatograms of Warfarin with all scouting conditions.  
(The number in chromatograms shows rank of Fig. 6)


Data processing software “Multi-data Report (MDR)” was able to pick the best separation chromatogram quickly by comparing the resolutions, number of detected peaks, and other variables. With this software, it is possible to compare the data quantitatively and thus it makes data processing more efficient (Fig. 6).

## Improvement of total analytical workflow by using online SFE-SFC

### LabSolutions LC/GC + CLASS-Agent + Agent Report



### LabSolutions DB/CS + Multi-Data Report (MDR)



11 min

Register To Agent  
Select data  
Select Template  
Open Excel  
Print out

0 min

Agent Report

MDR

<Time required for one report output>

Ranking	Run No.	Analytical Condition	Resolution	Selectivity	Symmetry factor		Retention factor		Area %		Detected Peaks
					Peak 1	Peak 2	Peak 1	Peak 2	Peak 1	Peak 2	
1	2	SFC_IA_EtOH	4.61	1.89	1.65	1.49	1.53	2.89	49.93	50.07	2
2	3	SFC_IA_MeOH	4.08	1.86	1.35	1.28	1.17	2.17	49.90	50.10	2
3	14	SFC_IE_EtOH	3.16	1.43	1.19	1.13	2.51	3.57	50.01	49.99	2
4	6	SFC_IB_MeOH	2.96	1.73	1.38	1.28	0.72	1.25	49.83	50.17	2
5	5	SFC_IB_EtOH	2.85	1.60	1.59	1.43	1.18	1.88	49.79	50.21	2
6	17	SFC_IF_EtOH	2.72	1.38	1.29	1.20	1.96	2.71	49.84	50.16	2

Fig 6: Workflow of outputting report which shows estimation result of the separation for Warfarin

## Online SFE-SFC analysis

Glass filter was prepared by dropping 5  $\mu$ L of rat blood plasma spiked with a warfarin standard. Fig. 7 shows the results from analyzing the warfarin isomer by online SFE-SFC (analytical conditions in Table 4). Warfarin chiral compounds were detected and separated enough after only fifteen minutes of SFE extraction.




Fig. 7: Chromatogram of rat blood plasma spiked with a warfarin standard

## Conclusion

- With the “Nexera UC HPLC/SFC switching system for chiral screening”, analytical conditions suitable for chiral compounds could be quickly determined. Furthermore, the data processing software “LabSolutions DB/CS + Multi-Data Report” achieved higher efficiency.
- Online SFE-SFC system can be used to reduce the work involved in pretreatment processes that were previously performed manually for research in pharmaceutical fields. The evaluation of quantitative analysis will be examined in future work.

First Edition: March, 2017