


Gel Permeation Chromatography System

# Nexera GPC System



# Nexera™ GPC System

Gel Permeation Chromatography System



## Analytical Intelligence



- Automated support functions utilizing digital technology, such as M2M, IoT, and Artificial Intelligence (AI), that enable higher productivity and maximum reliability.
- Allows a system to monitor and diagnose itself, handle any issues during data acquisition without user input, and automatically behave as if it were operated by an expert.
- Supports the acquisition of high quality, reproducible data regardless of an operator's skill level for both routine and demanding applications.

# Integrating Analytical Intelligence with Exceptional Instrument Performance

The Nexera GPC System provides highly reliable and highly extensible performance in a wide variety of applications. Hardware provides rapid instrument stabilization and excellent reproducibility of analytical results, while software includes analytical workflow automation and overlap injection functions and features a user-friendly analysis screen view, all leading to higher productivity.

In addition, Analytical Intelligence functionality based on state-of-the-art digital technology provides highly reliable analytical results by preventing potential problems and ensuring the instrument is always in a stable state, regardless of the skill, knowledge, or experience level of the analyst.




# Ensures Highly Reliable Analysis Regardless of an Operator's Skill Level

## Protects Expensive GPC Columns: FlowPilot



HPLC columns can be damaged by sudden pump starts and stops or extreme gradient changes. The Nexera automatically uses FlowPilot (Smart Flow Control) to increase the flow rate gradually to the set point. There is no need to create startup protocols for each analysis.




## Mobile Phase Levels Measured in Real Time



Reservoir tray weight sensors (optional) can be used to monitor the volume of mobile phase or autosampler rinse solution in up to twelve\* containers. The containers can also be checked remotely from a smart device.

You will no longer need to worry about running out of mobile phase mid-analysis, because the device will notify you before starting the run if the volume remaining is too low.

\* Up to 12 solutions can be monitored using 1-liter containers or up to 4 solutions using larger containers (2-liter or up to 5-liter containers).



Checks whether the quantity to be consumed by the analysis is available.




The mobile phase level can even be monitored from outside the laboratory.

## Auto-Diagnostics and Recovery



In rare cases, air bubbles can form in the mobile phase and cause problems if inhaled into the pump. The Nexera has the ability to monitor baseline changes and pressure fluctuations to check for abnormalities.

When it detects an unusual fluctuation, it can automatically pause the analysis, purge the flow path, and restart analysis once it has confirmed recovery to normal pressure.



# Automated Workflow





ANALYTICAL INTELLIGENCE

Functionality for assisting the entire series of analytical process steps, from instrument startup to shutdown, significantly increases work efficiency. FlowPilot mobile phase flowrate control functionality is linked to column oven temperature, so that columns are protected during automatic column equilibration. By automatically checking for baseline drift, instrument stability is maintained without the operator having to remain next to the instrument.

## Automated and Easy System Suitability Testing (SST)


The LabSolutions™ software is designed with a function to use a control sample or standard sample to check system compatibility prior to measurement of an unknown sample.

This function can be included in batch analysis to allow the actual sample to be measured only when a positive test result is obtained. The software can also output test results in text or CSV format.




## Shorter Analysis Times — Overlap Injection Function

Using the overlapped injection function, included standard in Nexera autosamplers (SIL-40 series), significantly shortens the time required to analyze multiple samples, which increases throughput.



Example of Using Overlapped Injection to Improve Analysis Cycle Time by About 50%



As with normal analysis, molecular weight can be calculated from the analysis window.

# Assembling Elements Demanded in GPC Analysis

## Excellent Reproducibility in Molecular Weight Distribution Analysis



For the analysis of molecular weight distribution, since the molecular weight of compounds eluted between the exclusion limit and permeation limit are considered an exponential function of elution time (calibration curve), even a small variation in elution time causes a substantial change in the result.

The Nexera GPC System utilizes high-speed micro-plunger actuation and automatic pulse compensation to provide pulse-free solvent delivery and excellent elution time reproducibility.

## Excellent Stability Maintained Even During Subtle Changes in Room Temperature

Nexera series detectors feature outstanding temperature-control functionality that provides superior baseline stability.

Since differential refractive index detectors are particularly sensitive to room temperature fluctuations, the RID-20A detector features dual temperature-controlled optics to ensure superior baseline stability by minimizing the effects of slight temperature changes.






## Identifying Unseparated Additives as Single Peaks



If LabSolutions software is used with a PDA detector, then peak convolution functionality (i-PDeA II) based on MCR-ALS (multivariate curve resolution alternating least square) analysis can be used to qualitatively and quantitatively identify peaks that are not completely separated by the column.

Peak convolution can also be used to confirm the purity of target peaks.



i-PDeA II chromatogram identifies the two separate peaks

# The RID-20A Differential Refractive Index Detector Allows Productivity Improvements in GPC Analysis


Inheriting the stability and extensibility that are the strengths of the Prominence series, the new RID-20A differential refractive index detector is designed with a new reference-cell auto-purge feature and validation support function that dramatically improve GPC analysis productivity.



## Reduced Stabilization Time and Improved Baseline Stability: After 30 Minutes Power ON

The RID-20A achieves shorter baseline stabilization time after turning ON the power through improved dual-temperature control of the optical system and superior lamp performance.

The stable baseline ensures reliable molecular weight distribution analysis.



## 30%\* Savings in Solvent Usage and Reduced Environmental Burden

The amount of mobile phase consumed can be saved by returning column eluate to the mobile phase bottle during intervals when no component peaks are eluted. This reduces the cost of analysis per sample and mitigates the environmental impact.


The recycle valve kit can be attached to the RID-20A and the SPD-40/40V.

When performing a 15-minute analysis 50 times (flow rate: 1 mL/min)

Without recycle valve kit ... 750 mL

With recycle valve kit ... 500 mL

→ Mobile phase usage reduced by approx. 30%\* (\* under specified conditions)




The use of the recycling valve cannot affect the baseline and the peak shape.

## Shimadzu's Proprietary Technology Supports Highly Sensitive Analytical to Preparative Analysis Applications

The four-partition photodetector in the RID-20A allows a wide refractive index range (0.01 to 5000 μRIU). The single detector supports all applications from highly sensitive measurements to preparative measurements using the three operation modes shown below.


A (Analytical) Mode	High-sensitivity to general-purpose analysis
P (Preparative) Mode	High-concentration analysis, semi-preparative analysis (up to 20 mL/min)
L (Large-scale prep.) Mode	Flow selection block allows large-volume preparative analysis (up to 150 mL/min)



# LabSolutions GPC Software

## Easy Analysis of Molecular Weight Distribution via Graphical Interface

### Create a Calibration Curve in Just Three Steps



**STEP 1** Drag and drop standard-sample data into the window.

**STEP 2** Load the retention times.

**STEP 3** Enter the molecular weight of the standard sample.


A maximum of 64 data points are available. Virtual points are also easy to set, and calibration curve appropriateness can be checked visually while choosing from a wide variety of approximation equations. Calibration curves can be corrected using the Mark-Houwink equation, and other correction methods based on Q-factors or degree of polymerization are available.

### Graphical GPC Data Analysis Window

- Manipulation of peak integration possible by means of the graphical interface
- Management of data from multiple detectors within a single file


Because the molecular weight distribution curve is updated whenever peak integration is performed, results for mean molecular weight, intrinsic viscosity, polydispersity, and other parameters can be confirmed immediately.

Time and detector sensitivity can also be corrected based on an internal standard peak or control sample.



### The Data Comparison Window Allows Simultaneous Evaluation of Multiple Samples


- The elution curves and derivative and integrated molecular weight distribution curves for up to 10 samples can be overlaid on a graph.
- Statistical results can be displayed for mean molecular weight, intrinsic viscosity, and polydispersity.



## Diverse Range of Report Formats

LabSolutions is equipped with report templates for a variety of analysis results. The software can accommodate a diverse range of reports thanks to a wealth of report items and a highly flexible layout.

In addition, a PDF output function is included as standard, so analysis results reports can be managed by automatic import to a database, helping your laboratory to go paperless and promoting an eco-friendly analysis process.




## Reduce the Work Involved in Creating a Final Report

Do you move your analytical results to a spreadsheet program (e.g., Excel®) to create a final report?

LabSolutions includes a multi-data report feature, which reduces the work involved in report creation. Analytical results are automatically entered into a spreadsheet equivalent to the one used in Excel®, eliminating the need to move the data.


### [Conventional Workflow for Report Creation]



Transfer of results reports to a spreadsheet program

- Time and effort associated with transfer tasks, entry errors, and double checking
- Problems of paper/file management
- Time and effort associated with regular verification work

### [Workflow of Report Creation with the Multi-Data Report]




Transfer of results is unnecessary

- Elimination of transfer errors
- Easy management of reports via the database
- Appropriate file protection provided by an audit trail feature

# Superior Extensibility of the Nexera Series


## Choose a Detector Designed to Your Objectives

A lineup of detectors is available, comprising the RID-20A differential refractive index detector, photodiode array detectors, UV-VIS detectors, and evaporative light-scattering detectors. Select the detector according to application, from organic solvent-based GPC applications for analysis of industrial polymers to water-based GPC applications for analysis of biopolymers. Also, using AccuSpot™ allows you to construct an automated analysis system capable of easy online processing—everything from GPC separation to MALDI-TOF MS measurement.



## Perform Diverse Analyses Within a Single Measurement

Using a differential refractive index detector together with a photodiode array detector enables high-sensitivity detection of peaks that are indistinct when analyzed with the differential refractive index detector alone. This makes it possible not only to measure molecular weight distribution, but also to identify and quantitate additive agents and impurities contained in polymers, all in a single analysis.



## Shim-pack™ GPC Series Columns for GPC Measurement **CoreFocus**

Shim-pack GPC series columns include polystyrene polymers with different degrees of cross-linking, optimized based on target molecules, from polymers to oligomers.

### Shim-pack GPC-800 Series High-Performance GPC Columns

#### For use with tetrahydrofuran (800 Series)

P/N	Product Name	Exclusion Limit Molecular Weight (polystyrene)	Size (length × inner diameter, mm)
228-20803-91	Shim-pack GPC-801	$1.5 \times 10^3$	300 × 8.0
228-20804-91	Shim-pack GPC-802	$5 \times 10^3$	300 × 8.0
228-20805-91	Shim-pack GPC-8025	$2 \times 10^4$	300 × 8.0
228-20806-91	Shim-pack GPC-803	$7 \times 10^4$	300 × 8.0
228-20807-91	Shim-pack GPC-804	$4 \times 10^5$	300 × 8.0
228-20808-91	Shim-pack GPC-805	$4 \times 10^6$	300 × 8.0
228-20810-91	Shim-pack GPC-80M (Mixed gel)	$2 \times 10^7$	300 × 8.0
228-20812-91	Shim-pack GPC-800P	Guard column	—

\* Please refer to the below URL and QR code when you use other organic solvent such as chloroform. This guide page explains a compatible column which is usable for your analysis and how to exchange the shipping solvent (contained in a new column) for the Shim-pack GPC series to a different solvent required for analysis.


[https://www.shimadzu.com/an/products/liquid-chromatography/hplc-consumables/shim-pack-gpc-series/features.html#anchor\\_1](https://www.shimadzu.com/an/products/liquid-chromatography/hplc-consumables/shim-pack-gpc-series/features.html#anchor_1)



### Shim-pack GPC-80M Mixed Gel Column

Shim-pack GPC-80M is a mixed gel column capable of analyzing samples over a wide range of molecular weight distributions.

#### Analysis of Polystyrene Standard



#### ■ Conditions

Column : Shim-pack GPC-80M  
(2 columns in series)  
(P/N: 228-20810-91)  
Mobile phase : Tetrahydrofuran  
Flow rate : 1.0 mL/min  
Column temp.: Ambient  
Detection : UV 254 nm



### Shim-pack GPC-2000 Series Preparative Columns

Shim-pack GPC-2000 Series are dedicated to preparative LC of tetrahydrofuran and chloroform. They ensure separation efficiency comparable to analytical columns as well as large preparative capacity.

## Nexera GPC System Standard Configuration

P/N	Description	Model	Qty
228-65502-58	System Controller	CBM-40* <sup>1</sup>	1
228-65002-58	Solvent Delivery Unit	LC-40D	1
228-65018-58	Degassing Unit	DGU-403	1
228-65101-58	Autosampler	SIL-40C	1
228-71762-46	Vial plate 1.5 mL		2
228-67723-42	Solenoid valve for drain (PTFE)		1
228-65202-58	Column Oven	CTO-40C	1
228-65306-58	Refractive Index Detector	RID-20A	1
228-65508-58	Reservoir Tray		1
228-53184-54	SUS Piping (I.D. 0.3 mm × 600 mm)		2
228-70247-42	40 Series Wiring Kit, B		1
—	Power cable		1
—	Power outlet unit		1
228-57647-43	Hand Tool Kit		1
228-38583-42	Glass Reservoir Bottle, w/ Cap, 5 pcs		1
228-77057-41	Bottle cap for GPC		2
228-65525-58	Mobile Phase Monitor (Controller)		1
228-65526-58	Mobile Phase Monitor (bottle holder)		1
074-83019-01	Mobile Phase Monitor AC Adapter		1
—	Analysis column	Shim-pack GPC Series* <sup>2</sup>	—
223-62700-92	LC Workstation	LabSolutions* <sup>3</sup>	1
223-19189-92	GPC Software	LabSolutions GPC Software* <sup>4</sup>	1

\*<sup>1</sup>: The CBM-40Alite (228-65501-58) and SCL-40 (228-65500-58) can also be chosen.

\*<sup>2</sup>: Choose according to the range of molecular weights to be measured.

\*<sup>3</sup>: A printer and cable are required separately.

\*<sup>4</sup>: LabSolutions DB GPC software (223-13812-92) and LabSolutions CS GPC software (223-18995-92) can also be chosen.

\* The addition of the dedicated HFIP kit may be required depending on the mobile phase composition. Please contact the Shimadzu sales staff for further information.

## Nexera GPC System Standard Specifications

Item	Specifications
System configuration	Modular
Measurement method	Single Flow
Mobile phase delivery method	Parallel Double Plunger
Mobile phase flow rate setting range	0.0001 to 10.0000 mL/min
Mobile phase flow rate accuracy	No more than 1% or 2 µL/min, whichever is greater (0.01 to 2 mL/min)
Mobile phase flow rate precision	No more than 0.06% RSD or 0.02 min SD, whichever is greater
Mobile phase degassing method	Vacuum membrane
Degassing line flow rate	400 µL
Sample injection method	Variable loop weighing
Sample injection volume setting range	0.1 to 100 µL (standard)
Number of samples processed	162 (1.5 mL vials)
Column temperature control method	Forced Air Circulation

Item	Specifications
Column temperature control range	(Room temperature + 10°C) to 85°C
Detector noise	2.5 × 10 <sup>-9</sup> RIU max.
Detector drift	1 × 10 <sup>-7</sup> RIU/h max.
Detector cell volume	9 µL
Detector cell temperature control range	30 to 60°C
Calibration curve approximation formulae	Linear, 3rd-order, 3rd-order + hyperbolic curve, 5th-order, 5th-order + hyperbolic curve, 7th-order, 7th-order + hyperbolic curve, or broken line
Calibration curve correction functions	Internal standard correction, Q-factor, RID sensitivity correction
Peak Integration	Automatic processing according to parameter settings. (Manipulation possible)
Molecular weight calculations	Mn, Mw, Mz, Mz1, Mv, Mw/Mn, Mv/Mn, Mz/Mw, and intrinsic viscosity
Data output	ASCII format

Nexera, Analytical Intelligence logo, LabSolutions, AccuSpot, AXIMA, Shim-pack and CoreFocus are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries. Excel is either a registered trademark or a trademark of Microsoft Corporation in the United States and/or other countries.



Shimadzu Corporation  
www.shimadzu.com/an/

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.