


Fully Automated Sample Preparation Module for LCMS





Fully Automated Sample Preparation Module for LCMS<sup>™</sup> Based on Shimadzu's Extensive Experience with Clinical Examination Systems for Blood Coagulation Analysis

Simply Place Blood Collection Tubes in Position and the Module Automatically Performs Everything from Sample Pretreatment to LCMS Analysis

## Achieves Reliable Data Acquisition

- Fully automates processes from pretreatment to LCMS analysis.
- Reduces variability in internal standard substance concentration (using vials with septa installed).
- Equipped with improved and expanded accuracy control functionality.
- Reduces operator errors and risk of infection.

# Eliminates the Need to Review Settings for Routine Work

- Automates processes from searching toxicological drug libraries to generating reports.\*1
- Provides the optimal analytical method.\*2

# Further Improves Efficiency of Research or Business Processes

- Makes it easy to create or review new pretreatment protocols.
- Samples can be added or processing interrupted in an emergency.\*3
- Text files can be used to load analysis requests from an LIS or output concentration values to a specified folder.

<sup>\*1:</sup> Requires the LC/MS/MS forensic toxicological database or rapid toxicology screening system ver. 2, which are sold separately.

<sup>\*2:</sup> Limited to commercial LCMS analysis kits and other methods available for supply.

<sup>\*3:</sup> The number of interruptions or samples that can be added is limited by the number of pretreatment vials placed.

# Innovation in Automation



Fully Automated Sample Preparation Module for LCMS

CLAM-2030

# CLAM-2030 Features

# Automates All Process Steps, from Pretreating to Measuring Blood, Urine, Serum and Plasma Samples

- Enables seamless analysis and management of samples throughout the entire process from pretreatment to LC/MS/MS measurement.
- Automation means samples can be analyzed automatically at nights or on days off.

2

# Supports Measuring a Wide Range of Compounds

- Flexible HPLC and MS configuration to meet application need.
- Measure hundreds of compounds in a single analysis.
- Easily change between different methods.

3

# **Checking Various Measurement Results Is Easy**

• Facilitates data review and flags inconsistent data and instrument conditions.



- Analysis of blood serum/plasma components requires
- pretreating samples by centrifuging in advance.

Prepare a blood collection tube.





#### Full automation

- Add blood or urine samples.
- Add internal standard substance.
- Add organic solvent to denature proteins.
- Stir.
- Filter using suction.
- Derivatize by heating (optional).
- Transfer samples to the HPLC autosampler.









- 6. Place vials in a rack.
- 7. Request analysis.



- Complete integration of CLAM-2030 with LCMS enables seamless analysis.
- Successfully tested using a wide range of compounds (full list available on page 8).



- Check quantitative results (concentration values) and the mass chromatograms for each component.
- Display calibration curve accuracy % (precision) and data alarms with respect to the measurement results for substances subject to accuracy control.
- As qualitative analysis, output reports on library search results.




# Reduces Operator Errors and Risk of Infection During Manual Operations

## **Reduces Operator Errors**

By simply placing uncapped blood collection tubes (or sample cups) and pretreatment vials in the system, and requesting analysis, the system performs all other process steps automatically, from pretreatment to LCMS analysis. With an LCD touch panel screen and a user interface that requires no instruction manual, the system helps ensure reliable and rapid operations with minimal operator errors due to manual operations.




Carousel for samples and reagents
Racks for pretreatment vials
Dispensing probe for samples
Dispensing probe for reagents



### **Reduces Infection Risk**

The CLAM-2030 minimizes human contact with potentially hazardous samples through automation. Following pretreatment, solid and liquid waste are contained within the body of the CLAM-2030 and collected. Integrated sensors alert users through the software that waste levels are high and require emptying.



Fully Automated Sample Preparation Module for LCMS

7

# Applications

The follow are examples of applications using a CLAM module.

### Target Components

6-acetylmorphine, Aldosterone, Amikacin, Amiodarone, Amisulpride, Amphetamine, Androstenedion,

Anhydroecgonine methyl ester, Arbekacin, Aripiprazole

8 Barbiturate drug and Bromovalerylurea, BDB, Benzoylecgonine, 39 Benzodiazepines and their metabolites, Bepridil, Buprenorphine

Carbamazepine (CBZ), Carbamazepine10-11-epoxide, 2C-B, 2-Cl, CBZ-Diol, Cibenzoline, Ciprofloxacin,

Clopazine, Cocaethylene, Cocaine, Codeine, Corticosterone, Cortisol, Cyclosporine

Daptomycin, Dehydroaripiprazole, 11-Deoxycortisol, Desmethylclopazine, Dextromethorphan, DHEA,

DHEA-S, Diazepam, Dihydrocodeine, DOACs (Apixaban, Rivaroxaban, Edoxaban, Dabigatran)

Ecgonine methylester, EDDP, Ethosuximide, Ethylmorphine, Everolimus

Felbamate, Flecainide

Gabapentin, Gentamicin

Haloperidol, Hydoromorphone, Hydrocodone, 25-hydroxyvitamin D3, 25-hydroxyvitamin D2, Hygromycin C

Kanamycin, Ketamine

Lacosamide, Lamotrigine, Levetiracetam, Levofloxacin, Levomepromazine, Linezolid

6-MAM, MBDB, MDA, MDEA, MDMA, Melperone, Mephedrone, Metanephrin, Methadone, Methamphetamine,

Methaqualone, Methcathinone, Methiopropamine, 3,4-Methylenedioxypyrovalerone, Methylphenidate, Mexiletine Morphine, *m*-CPP, 4-MTA

Naloxone, Naltrexone, N-Desmethylolanzapine, NDMS, Neomycin B, Nimetazepam, Norbuprenorphine, Norephedrine, Norfenfluramine, Norketamine Noroxycodone, Normetanephrin, Norpseudoephedrine, Norguetiapine, n-Desethylamiodarone

10-OH-CBZ, 17-OH-Progesterone, 9-OH-Risperidone, Olanzapine, Oxacarbazepine, Oxycodone

Paromycin, PEMA, Perampanel, Perazine, Phencyclidine, Phenobarbital, Phenytoin, Pholcodine, Pilsicainide,

Pipamperone, Pregabalin, Primidone, Progesterone, Promethazine, Propoxyphene,

Pyridoxal-5'-Phoshate (P5P)

Quetiapine Risperidone

Rapamycin (Sirolimus), Retigabine, Ritalinic acid, Rufinamide

Sertindole Sulpiride, Sotalol, Stiripentol, Streptomycin, Sultiame

Tacrolimus (FK 506), Teicoplanin, Teststerone, Thiamin Pyrophosphate (TPP), Thioridazine, Tiagabine,

Topiramate, 12 Tri-/Tetra-cyclic antidessant

Valpronic acid, Vancomycin

Ziprasidone, Zonisamide, Zotepine Zuclopenthixol, (Z)-Chorprothixene

# High-Speed High-Separation Analysis Connectible to LCMS Systems

CLAM-2030 modules can be connected to Nexera LC-MS/MS ultra high performance liquid chromatograph systems, which offer exceptional speed and resolving power. By connecting a CLAM-2030 module to a Nexera LC-MS/MS system, the CLAM-2030 control software can be used to acquire data by automatically analyzing pretreated samples in the LCMS system. The MS detector can be selected based on the target components, analysis frequency, and sensitivity requirements of the customer.

### **MS** Detectors

Low

#### LCMS-8040

This triple quadrupole mass spectrometer offers increased sensitivity (compared to previous Shimadzu models), based on advancements in ion optics.

### LCMS-8045

Triple quadrupole mass spectrometer features a newly designed heated ESI probe and collision cell.





#### Sensitivity

### LCMS-8050

With a heated ESI probe and the UFsweeper<sup>™</sup> III collision cell, the LCMS-8050 achieves outstanding measurement speeds and quantitative performance.



### LCMS-8060

This triple quadrupole mass spectrometer features Shimadzu's patented "UF Technologies."

By including an optimized ion loading unit and new UF-Qarray<sup>™</sup> technology, it achieves even higher sensitivity than the LCMS-8050.



High

Please consult your local sales representative for more information on which system is best suited for your application.

9

# **Optional Software**



## LC/MS/MS Rapid Toxicology Screening System Ver. 2

This system includes methods for simultaneously or individually analyzing 161 components, including abused drugs, psychiatric and neurological drugs, and hypnotics and sedatives.

## Predetermined Analytical Conditions

This system provides MRM transitions, LC separation conditions, retention times, spectral libraries, and report files, enabling rapid implementation of a screening method for many important toxicological substances.

### Screening for 161 Compounds

Developed in conjunction with front-line researchers in the toxicology field, the method analyzes several classes of substances. LCMS MRM optimization has been performed for common drugs of abuse, hypnotics, psychotropics, pharmaceuticals and other natural toxins. Shimadzu's Ultra-Fast MS technology enables triggered product ion scans, which provide MS/MS data during MRM acquisition.



# LC/MS/MS Forensic Toxicological Database

This database includes compound data necessary for analyzing abused drugs, psychiatric and neurological drugs, pharmaceuticals, pesticides, and natural toxins, and for other toxicological analysis.

### Simplified Acquisition and Identification

The Forensic Toxicology Database includes optimized LC-MS/MS data acquisition parameters and a library database to help clinical and forensic researchers build screening and quantitation methods quickly, simplifying method development.

## MRM & Spectral Library Database Contains Information on More Than 2,500 Compounds

The spectral library database is built using two separation conditions (ODS and Biphenyl). Both methods have information on clinical and forensic compounds of interest in routine analysis. The ODS method contains information on 1,250 compounds and the Biphenyl method contains 1,281 compounds. Compound datasheet includes: monoisotopic mass, RT, CAS number, formula and compound class. This package provides Synchronized Survey Scan<sup>™</sup> parameters (MRM parameters, MRM intensity threshold and triggered product ion scan parameters) optimized for screening analysis.

## Enhanced Identification by Merged Spectrum

Each certified reference material was acquired with three different collision energies to generate an information-rich merged-CE spectrum which can be used in library matching and compound verification. Matching with a merged-CE spectrum library can be a powerful tool to identify compounds with a library score.

# Consumables

### Specifications for Specialized Pretreatment Vials

Dedicated disposable pretreatment vials are used. Filter and collection vials are used as a pair.





Specialized	
Filter Vial	

Specialized Collection Vial

P/N	Description
241-16593-41	Vial Set (set of 100 pairs)
241-16593-42	Vial Set (set of 500 pairs)
241-16593-43	Vial Set (set of 2,000 pairs)
Vial Material	Polypropylene (PP)
Filter Material	PTFE, with 0.45 µm pore size
Reusability	None (disposable)

# Specifications of Applicable Sample and Reagent Vials

Sample Vials	Image: A problem of the point of the poin	
	1.5 mL Vial   P/N 228-16965-91   Vial (100 pieces in a bag)     P/N 241-18511-41   Septum sets for 1.5 mL vial (100 sets in a bag)     • Caps   • Septum	
Reagent Vials	6 mL Vial   P/N 241-16619-43   Vial (50 pieces in a bag)     P/N 241-18512-41   Septum sets for 6 mL vial (50 sets in a bag)     • Septum   • Filter paper     P/N 038-00199-33   Cap, 6 mL	
	12 mL Vial     P/N 223-81340-41   Vial     P/N 223-81468-02   Septum to Prevent Absorbing Humidity     P/N 038-00199-34   Cap, 12 mL	

11

►

CLAM, LCMS, Nexera, UFsweeper, UF-Qarray and Synchronized Servey Scan are trademarks of Shimadzu Corporation. Vacutainer is a registered trademark of Becton, Dickinson and Company. VENOJECT is a registered trademark of Terumo Corporation.



Shimadzu Corporation www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures. This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.