

Benchtop Linear MALDI-TOF Mass Spectrometer

MALDI-8020




MALDI-8020: A compact footprint without compromising on performance


The MALDI-8020 is the latest in a long line of MALDI-TOF products from Shimadzu. This benchtop, linear MALDI-TOF mass spectrometer delivers outstanding performance in a compact footprint, making it an ideal choice for today's increasingly demanding laboratories.

Key features:

- Linear mode (positive ion) MALDI-TOF
- 200 Hz solid-state laser, 355 nm
- Load-lock chamber for fast sample introduction
- UV laser-based source cleaning (patented)
- Small footprint/benchtop design
- Quiet operation (<55 dB)



Linear mode instruments have traditionally been the instruments of choice in MALDI-TOF-based quality control (QC) and profiling workflows, and the MALDI-8020 is no exception. Peptides, proteins, polymers or oligonucleotides - the instrument is capable of delivering the performance required for QC applications. Research laboratories will also benefit from the MS capabilities of the instrument through the rapid mass-measurement of samples.




Example of a 4.2 kDa peptide demonstrating an average resolution of >5000 FWHM. Resolution, calculated by the software, is shown in brackets.

Sample target solutions

Compatible with the FlexiMass™ series of microscope slide-format sample targets, these slides provide researchers with options depending on their experimental workflow.

The individually barcoded, single-use FlexiMass™-DS slides provide a convenient solution for more routine or defined workflows. Ready-to-use, these disposable targets eliminate the need for cleaning and the risk of carryover. Alternatively, the reusable stainless steel FlexiMass™-SR sample targets provide a cost-effective, longer-term solution to sample preparation.




MALDI-In-Source Decay (MALDI-ISD) spectrum of native Exendin-4 (top panel) showing matched fragment ions consistent with the amino acid sequence. Bottom panel: expanded region of the mass spectrum showing the isotopically resolved fragment ions.


TrueClean™ automated source cleaning

To maintain instrument performance over time, the MALDI-8020 features wide-bore ion optics - a feature used on our existing MALDI-TOF products - which minimize the risk of source contamination over time, providing a robust platform.

To maximize uptime, the system is equipped with TrueClean - an automated, rapid (<10 min) UV laser-based source cleaning which can be used to clean the extraction electrode *in-situ* without breaking instrument vacuum.




Polymer analysis: polycarbonate from a CD (green spectrum) and polycarbonate standard (orange spectrum). Inset: expanded view of the highlighted region.



Analysis of chlorhexidine gluconate in mouthwash solution. Sample was prepared with MALDI matrix without any clean-up or dilution. Inset: expanded view showing peak at 505.2 m/z.

Facilitates 21 CFR Part 11 compliance

Operating under the control of *MALDI Solutions™* software, the software features a centralized, secure Microsoft® SQL database which can be used to store everything from sample lists and acquisition Parameter Sets to acquired MALDI data. The system is managed by an Administrator and customizable user profiles provide control over access to the database and operation of the instrument. Along with a full audit trail, the software provides the tools to help achieve 21 CFR part 11 regulatory compliance.



Example of higher mass analysis using immunoglobulin A (IgA)


Quality as standard

As mass spectrometry moves into clinical laboratories, the MALDI-8020 is the perfect choice for researchers developing MALDI-based diagnostic methods[†]. The system is designed and manufactured to ISO 9001:2015 and ISO 13485:2016 quality management standards and is CE compliant.

Through the selection of components demonstrating improved longevity and performance, and by reducing the number of component parts, our engineers have created a reliable and robust platform capable of delivering outstanding performance in a small footprint.

The MALDI systems are maintained by trained service support engineers who form part of the Shimadzu global distribution and support network. Routine maintenance is simplified thanks to the instrument design, providing our engineers with easy access to commonly serviced parts.

[†] For Research Use Only. Not for use in diagnostic procedures.



The MALDI-8020 is designed and manufactured under the Kratos Analytical Ltd. quality management system and is CE compliant, CB certified and is tested to UL standards. The MALDI-8020 is a Class 1 laser product.



Kratos Analytical Ltd.
 Wharfside, Trafford Wharf Road, Manchester M17 1GP, UK
 Phone: +44 161 888 4400 Fax: +44 161 888 4402
 URL <http://www.shimadzu.com/an/>



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice.

Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.