



Screening System for Phthalate Esters





Making the Difficult Simple

The Py-Screener[™] system is designed to screen for phthalate esters in polymers. The use of phthalate esters in toys and food packaging is currently restricted. Moving forward, they are expected to be regulated as restricted substances under the Directive on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS (II) Directive) in Europe. The pyrolyzer GC/MS (Py-GC/MS) is used to selectively detect and quantify phthalate esters thermally extracted from samples. This screening system consists of a sampling toolkit, special standards, and special software and can be easily operated even by novices.



Easy to Operate Using Special Software

Using customized software, operations are easy, even for novices. To automatically start continuous analyses, just place the prepared standards and test samples in the autosampler, and enter the number of samples, the sample names, and their weights. Phthalate ester inspections compliant with the IEC62321-8 international analysis standard, the industry standard for RoHS inspections, can be performed easily by anyone. Shimadzu

offers both a simultaneous inspection method for phthalate esters and brominated flame retardants, which supports inspections for a wide range of regulated compounds, and a special high-speed inspection method for phthalate esters, which significantly reduces the inspection time.

8	Vial#	Sample Name	Sample Amt.
1	1	Blank_Cup	0.5
2	2	Phthalate_STD_Blank	0.51
3	3	Phthalate_STD_100	0.51
4	4	Phthalate_STD_1000	0.51
5	5	ERM-EC591	0.51
6	6	Test_Sample	0.51
7	7	Test_Sample	0.51
8	8	Test_Sample	0.51
9	9	Test_Sample	0.51
10	10	Test Sample	0.51

Easy to Operate Even for Novices Tabular Display of Concentrations and Criteria Clarifies the Results.

The concentrations of target components detected in continuous measurements are displayed in a table and color-coded using criteria based on concentration ranges. The results for continuously measured test samples can be checked at a glance. Also, the system is equipped with accuracy control functions in order to ensure the reliability of blank concentrations, instrument sensitivity, and other data, so even novices can feel confident that they are reporting reliable measurement results.



Py-Screener Software

 $\diamond - \times$


Data Processing

Maintenance

The special software displayed on the monitor helps you navigate the required procedures. Even novices can operate the system using the software.

Ample Maintenance Support Maintenance Navigation Supports Long-Term Operation with Periodic Replacement Kits


Using the Maintenance Navigator, the procedures appropriate for pyrolyzer and GC-MS maintenance can be performed easily and confidently. Procedures for locating and resolving leaks are included, and kits with common replacement parts simplify maintenance and troubleshooting procedures.



A Total Solution Proposed by Shimadzu for the RoHS (II) Directive

This system can be applied to polybrominated biphenyls (PBBs) and polybrominated diphenyl ethers (PBDEs), two types of brominated flame retardants already regulated under the RoHS (II) Directive.

In addition to this system, Shimadzu provides an Energy Dispersive X-ray (EDX) fluorescence spectrometer for inorganic compound screening and various other analytical systems for accurate quantitation. These systems provide a total solution for everything from screening to the accurate quantitation of substances already regulated under the RoHS (II) Directive, and substances for which regulation is anticipated. It is possible to create a report by integrating the inspection data obtained by EDX and GCMS, and since Shimadzu offers RoHS-compatible instruments, a comprehensive management system can be constructed.



Applicable Systems and Software

GC-MS	: GCMS-QP2020 NX, GCMS-QP2020, GCMS-QP2010 Ultra		
Pyrolyzer	: EGA/PY-3030D multi-shot pyrolyzer		
Autosampler	: AS-1020E auto-shot sampler		
GC/MS Workstation : GCMSsolution™ (Ver.4.30 or later, ver.4.50 or later for GCMS-QP2020 NX) + LabSolutions Insight™			
Py Workstation	: EGA-PY3030 program (Ver. 1.54 or later)		

Caution

1. Note that there are no guarantees regarding the accuracy of the information contained in the method files, or the usefulness of the information obtained from the results of their use. 2. In order to accurately identify the registered substances, perform the measurements using the system conditions in the method files contained in the product.

Py-Screener, GCMS-QP, GCMSsolution and LabSolutions Insight are trademarks of Shimadzu Corporation.



Shimadzu Corporation www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures. This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Source and your source of the Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "@". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.