

Ultra High Performance Liquid Chromatograph

# Nexera XS inert





SCL-40  
SYSTEM CONTROLLER

SHIMADZU

CONTROL ERROR

DGU-405  
DEGASSING UNIT

SHIMADZU

SPD-M40  
PHOTO DIODE ARRAY DETECTOR

SHIMADZU

0.3000ml MAX 105.0  
79.4MPa MIN 0.0  
PUMP | CONC | RUN | PURGE

LC-40D XSi  
SOLVENT DELIVERY MODULE

SHIMADZU

SET 25 ACTUAL 25.0 T. MAX 90  
OVEN | HEAT | RUN

Controls: CE, ^, <, >, v, del, Enter

CTO-40S  
COLUMN OVEN

# EXPERIENCE NEWFOUND CLARITY

For many biopharmaceuticals and medium-molecule drugs, it can be difficult to acquire reliable data using a typical HPLC system, due to adsorption caused by interactions with metal ions. That tendency is particularly true for molecules that contain phosphates, carbonyl groups, or other characteristic elements or functional groups, the adsorption of which can affect analysis in various ways, such as by causing poor peak shape or inadequate sensitivity. Therefore, to increase the sensitivity and accuracy of analysis, it is important to inhibit the adsorption of compounds to the system.

Nexera XS inert UHPLC systems achieve a 105 MPa pressure capacity while inhibiting adsorption to metals by using PEEK or ceramic materials in flow channels where samples flow. They also offer excellent resistance to corrosion and a broad range of pH levels, meaning they can be operated reliably using mobile phases that contain high haloid salt concentrations or have extreme pH levels.

## **Unconstrained Recovery and Sensitivity**

Reduces sample loss due to adsorption to metal and achieves excellent sensitivity.

## **Clear Resolution without Restrictions**

Improves peak shape and achieves excellent chromatographic separation.






## **Assured Reliability and Reproducibility**

Reliable data for metal-adsorbing compounds with high reproducibility.



## Addition of a New Bioinert Model to the Nexera Series

In addition to the standard Nexera models made with stainless steel for excellent pressure resistance, broad applicability, and robustness, a new bioinert model has been added to the Nexera series product line. The new Nexera XS inert model supports UHPLC analysis while also featuring PEEK and ceramic materials for flow channels exposed to sample flow. Similarly, the Nexera lite inert is an HPLC model with metal materials completely eliminated from flow channels.

	Bioinert HPLC/UHPLC	HPLC/UHPLC (stainless steel-based)
UHPLC <sup>*1</sup>	 <b>Nexera XS inert</b>	 <b>Nexera X3</b> <b>Nexera XS</b>
UHPLC-like <sup>*2</sup>		 <b>Nexera XR</b>
HPLC <sup>*3</sup>	 <b>Nexera lite inert</b>	 <b>Nexera lite</b>


\*1 ~130 MPa(X3), ~105 MPa(XS, XS inert)

\*2 ~70 MPa(XR)

\*3 ~44 MPa(lite), ~20 MPa(lite inert, for aqueous mobile phase with inert kit)


## Improved Peak Shape and Separation

The Nexera XS inert uses unique proprietary technology that minimizes exposure to metal ions on internal surfaces. That inhibits adsorption of target molecules within sample flow channels, ensures good peak shape, and provides excellent separation.



## Adsorption Inhibited and Detection Sensitivity Improved

The Nexera XS inert prevents target component recovery rates from decreasing by inhibiting adsorption caused by interactions between metal ions and proteins, nucleic acids, or other components. As a result, the Nexera XS inert system can detect phosphorylated peptides and other molecules with a strong tendency to adsorb to metals.




Example of Phosphopeptide Analysis

# Assured Reliability and Reproducibility


## Consistent Analysis Reproducibility

In order to reduce the effect of adsorption of metal-sensitive compounds, a system is often passivated by repeatedly injecting samples that contain the target compounds before starting data acquisition. However, this approach not only wastes valuable samples and time, but can make it extremely difficult to acquire reliable data due to changes in the state of passivation during continuous analysis. Nexera XS inert systems eliminate the need for the preliminary passivation and provide highly reliable data from the first injection and throughout the analytical session.



## Improved Quantitative Analysis Performance



Adsorption of target compounds tends to be especially apparent in low-concentration regions, which can worsen calibration curve linearity, as well as detection and quantitation limits. Adsorption can also cause carryover. The Nexera XS inert achieves a wider dynamic range and higher quantitative accuracy in low-concentration regions by minimizing adsorption from interactions with metals.



Accuracy and Contribution Rate ( $R^2$ ) at Each Concentration

Spiked Conc. (ng/mL)	Nexera XS inert		Standard UHPLC (stainless steel-based)	
	Measured Conc. (ng/mL)	Accuracy (%)	Measured Conc. (ng/mL)	Accuracy (%)
0.5	0.57	113.5	2.28	455.7
1	0.93	93.0	-1.04	-104.4
5	5.42	108.0	2.43	48.5
10	9.13	91.3	5.62	56.2
50	50.28	100.6	26.63	53.3
100	92.77	92.8	76.39	76.4
500	497.04	99.4	588.74	117.7
1000	1010.36	101.0	965.46	96.5
Contribution Rate ( $R^2$ )	0.9996		0.9721	

Carryover is a phenomenon where sample substances remaining inside the system or column after analysis, due to adsorption or other reasons, affect the next analysis. Therefore, carryover should be inhibited particularly for quantitative analysis of trace quantities, because of the large detrimental effects it can have on analytical results. Using a Nexera XS inert system can practically eliminate carryover by inhibiting adsorption to metal in the system.



Carryover test of oligonucleic acid

# Full-range of Bioinert UHPLC Capabilities

The Nexera XS inert system eliminates the risk of sample adsorption or surface corrosion, while still providing all the exceptional features of the Nexera series, making it the perfect solution for a wide variety of applications.



The Nexera XS inert system combines the inertness of the sample flow path with excellent resistance to corrosion. It delivers ideal chromatographic separation performance for a wide range of applications and excellent durability, even in the case of mobile phases with high salt concentrations or extreme pH values.

## System Controller **SCL-40, CBM-40/40lite**

Used to coordinate actions of the overall system, this controller offers intuitive operation that minimizes the stress involved in operating the system, from startup to shutdown.

## Detector **SPD-40/40V/M40**

Inert-type cells for UHPLC analysis eliminate the risk of adsorption on the detector. A low-diffusion cell, which has 5 mm of optical path length, can also be selected.

## Solvent Delivery Unit **LC-40D XSi**

Designed with corrosion-resistant non-stainless steel materials, it offers rugged, low-pulsation performance, advanced AI features and solvent blending capabilities. (optional).

## Autosampler **SIL-40C XSi**

This high-performance autosampler features non-metal materials for all surfaces that contact liquids. That inhibits metal-adsorption of biomolecules.

## UHPLC Inert Switching Valve **FCV-0206H2i/FCV-0607H2i**

Designed with adsorption-inhibiting materials for all wetted surfaces.

## Online pH Monitor **pHM-40 (optional)**

This monitors the pH level of mobile phases in real time and records the data in a data file together with chromatograms and instrument information.






## Inhibits Metal Adsorption and Achieves 105 MPa Pressure Capacity

To reduce the various chromatogram problems that metal adsorption can cause, the Nexera XS inert features sample flow channels made with PEEK material. Though it is generally difficult to use PEEK plastic for high-pressure conditions, the PEEK-lined tubing (developed independently by Shimadzu) used in the Nexera XS inert offers a maximum pressure capacity of 105 MPa\*. Furthermore, a proprietary butt-sealed tubing configuration eliminates the need for ferrules and results in tubing joints with no dead volume. In addition, a proprietary ratcheted fitting design allows column connections to be hand-tightened even at high pressures.

\* Also supports mobile phases spiked with HFIP (max. 5 %).



Zero dead volume with  
butt-sealed tubing configuration




Columns connected easily with  
ratcheted fittings

## High Corrosion Resistance and Broad pH Resistance

Nexera XS inert systems also offer superior corrosion resistance. To support mobile phases that contain high haloid salt concentrations for analyzing proteins, for example, corrosion-resistant metal is used for mobile phase contact surfaces upstream from the sample injection unit, such as in the solvent delivery pump and mixer, whereas PEEK plastic and ceramic materials are used for sample contact surfaces. In addition, the superior pH resistance of PEEK plastic ensures the system can be used to reliably acquire data for long periods even when using high-pH\* mobile phases, such as for oligonucleic acid analysis.

\* Supports pH levels from 1 to 14.



Example of Analyzing mAb Aggregates Using a Mobile Phase with a High Haloid Salt Concentration

## Improves Efficiency of High-Sensitivity Analysis and Characterization of Oligonucleotide Acid

### LabSolutions Insight™ Biologics

LabSolutions Insight Biologics is a dedicated software platform for oligonucleotide characterization using the LCMS-9030 or LCMS-9050 quadrupole time-of-flight type (Q-TOF) mass spectrometer. Main product and impurity identification includes several core editors for sequence, nucleotides, linkers, ribose and base modifications. Together with processing and integration for target modifications, Insight Biologics is a comprehensive workspace for data review, processing, and reporting. Because the Nexera XS inert system prevents the adsorption of target components to metal surfaces, combining it with LabSolutions Insight Biologics creates an ideal setup for oligonucleotide characterization.



## A Superior User Experience

### UX 1


#### Depiction of the Structural Formula

In the window for setting the oligonucleotide sequences, the structural formula of the sequence that was entered is displayed in real time, enabling quick and easy verification of the information. Also, the nucleobases used in the sequence, as well as the backbone linker, ribose, and base modifications, can be added and edited.


##### Input oligonucleotide sequences

#	Name	Base	Base Modification	Linker	Ribose	Formula	Mono-isotopic
1	Cd	Cytosine	None	----	Deoxy	C <sub>9</sub> H <sub>11</sub> N <sub>3</sub> O <sub>2</sub>	193.08513
2	sTd	Thymine	None	Phosphorothioate	Deoxy	C <sub>10</sub> H <sub>13</sub> N <sub>2</sub> O <sub>6</sub> P S	320.02319
3	sGd	Guanine	None	Phosphorothioate	Deoxy	C <sub>10</sub> H <sub>12</sub> N <sub>5</sub> O <sub>5</sub> P S	345.02968
4	sCd	Cytosine	None	Phosphorothioate	Deoxy	C <sub>9</sub> H <sub>12</sub> N <sub>3</sub> O <sub>5</sub> P S	305.02353
5	sTd	Thymine	None	Phosphorothioate	Deoxy	C <sub>10</sub> H <sub>13</sub> N <sub>2</sub> O <sub>6</sub> P S	320.02319
6	sAd	Adenine	None	Phosphorothioate	Deoxy	C <sub>10</sub> H <sub>12</sub> N <sub>5</sub> O <sub>4</sub> P S	329.03476
7	sGd	Guanine	None	Phosphorothioate	Deoxy	C <sub>10</sub> H <sub>12</sub> N <sub>5</sub> O <sub>5</sub> P S	345.02968

The structural formula of sequences is shown in real time



Base modifications can be added and edited




### UX 2

#### Display of Fragment Coverage


The software includes a coverage display which indicates fragment spectral assignments. The coverage display switches to match the items to be checked. Reports can also be output.

Checking ion intensity and coverage




Switchable coverage display

Checking the fragment sequence




Modification positions are also clearly identified



### UX 3

#### Component Chromatogram Display

The impurity peaks are displayed as a component chromatogram. The UV and MS chromatograms can be checked simultaneously.



## Achieving Maximum Inhibition of Metal Adsorption in Combination with Nexera XS inert System

### Shim-pack Scepter Claris

The Shim-pack Scepter Claris is a new metal-free column with a newly developed bioinert-conditioned column body. It is packed with a Scepter series stationary phase that offers superior durability and chromatography performance.

- Bioinert coating is applied to the column body and stainless steel frit
- Ideal for analysis of metal-coordinating and hydrophobically adsorbing compounds such as nucleic acids, proteins, and lipids
- Outstanding pH and lifetime stability due to Scepter organic silica hybrid packing

Shim-pack Scepter	Reversed Phase				
	C18-120	C18-300	HD-C18-80	C8-120	C4-300
Ligand Type	Trifunctional C18	Trifunctional C18	Trifunctional C18	Trifunctional C8	Trifunctional C4
	General Purpose	General Purpose for Large Molecules	High Density Type for Increased Retention		
Particle	Organic Silica Hybrid				
Particle Size	1.9 µm, 3 µm, 5 µm				
Pore Size	12 nm (120Å)	30 nm (300Å)	8 nm (80Å)	12 nm	30 nm
End Capping	Proprietary				
pH Range	1 - 12				1 - 10
100% Aqueous Condition	YES	YES	x	x	YES
USP Classification	L1	L1	L1	L7	L26


  

Shim-pack Scepter	Reversed Phase		HILIC
	Phenyl	PFPP	Diol-HILIC
Ligand Type	Trifunctional Phenylbutyl	Trifunctional Pentafluorophenylpropyl	Trifunctional Dihydroxypropyl
Particle	Organic Silica Hybrid		
Particle Size	1.9 µm, 3 µm, 5 µm		
Pore Size	12 nm (120Å)		
End Capping	Proprietary	None	
pH Range	1 - 10	1 - 8	2 - 10
100% Aqueous Condition	YES	YES	—
USP Classification	L11	L43	L20

### Superior Sensitivity and Separation Performance in Nucleic Acid Analysis

Shim-pack Scepter Claris C18-120 (bioinert coating) and Scepter C18-20 (stainless steel body) were compared in this analysis of Oligodeoxythymidylic acid.

Results from Scepter C18-120 show low peak intensity, suggesting adsorption on metal surfaces. In contrast, Scepter Claris C18-120 produced a sharp, high-intensity peak.



## Ion Exchange Chromatography Columns for Protein and Nucleic Acid Analysis

### Shim-pack Bio IEX Columns

The Shim-pack Bio IEX columns are available in the following chemistries.



- Anion exchange quaternary ammonium (Q)
- Cation exchange sulfopropyl (SP)

The following porous and non-porous particles, both of which use hydrophilic polymer packing, are available.

- Hydrophilic porous polymer (Q and SP columns)
- Hydrophilic non-porous polymer (Q-NP and SP-NP columns)

The porous packing materials offer highly efficient and excellent binding capacity, and the non-porous packing materials offer high recovery rates and excellent resolution.

	Q-NP	SP-NP	Q	SP
Particle Material	Hydrophilic non-porous polymer		Hydrophilic porous polymer	
Particle Size	3 $\mu\text{m}$ , 5 $\mu\text{m}$		5 $\mu\text{m}$	
Ion Exchange Group	- $\text{CH}_2\text{N}^+(\text{CH}_3)_3$	- $(\text{CH}_2)_3\text{SO}_3^-$	- $\text{CH}_2\text{N}^+(\text{CH}_3)_3$	- $(\text{CH}_2)_3\text{SO}_3^-$
Operating pH Range	2 - 12			
Operating Temp. Range	4-60°C			
Column Material	PEEK			



Application




## Size Exclusion Chromatography Columns for Mid- and Large-Size Molecules

### Shim-pack Bio Diol Columns

Shim-pack Bio Diol columns are size exclusion chromatography columns with packing particles of various pore sizes for analyzing aggregates and fragments of antibodies, oligonucleotides, carbohydrates, and other analytes. The lineup covers a wide range of use, from traditional analysis to laboratory-scale purification. In addition, the 2  $\mu\text{m}$ -particle columns for high-speed analysis enable rapid characterization.

	Diol-60	Diol-120	Diol-200	Diol-300
Particle Material	Silica gel			
Functional Group	Dihydroxypropyl (Diol)			
Particle Size	3 $\mu\text{m}$ , 5 $\mu\text{m}$		2 $\mu\text{m}$ , 3 $\mu\text{m}$ , 5 $\mu\text{m}$	
Pore Size	6 nm (60Å)	12 nm (120Å)	20 nm (200Å)	30 nm (300Å)
Operating pH Range	5.0 - 7.5			
Molecular Weight Range	10,000 or less	1,000 - 100,000	5,000 - 300,000	20,000 - 1,000,000
Maximum pressure	2 $\mu\text{m}$ : 45 MPa (Usually at 30 MPa or less)		3, 5 $\mu\text{m}$ : 20 MPa	




Application



## Improving the Efficiency of Process Steps in Analytical Method Development

### LabSolutions™ MD


When evaluating manufacturing processes or testing the stability of pharmaceuticals, separating principal components and impurities in chromatograms is very important for characterization analysis and quality control, but establishing the appropriate separation parameters requires a significant amount of time and effort. LabSolutions MD enables efficient development of analysis methods using an Analytical Quality by Design (AQbD) approach based on scientific evidence and risk. All steps involved in the screening, optimization, and validation phases of the method development process can be completed using LabSolutions MD, from analyzing samples by automatically switching between respective mobile phases and columns to using the analytical results to build a design space, deciding the optimal analytical condition settings, and then evaluating robustness (validation).



#### Example 1:

### Nexera XS inert with LabSolutions MD Provides Optimized Separation Parameters

When developing methods for testing substances related to oligonucleic acids, separation patterns can differ depending on the chain length and base composition, the presence/absence of modifications, and other factors. Therefore, separation must be optimized for each target sequence. Good peak shape and separation can be achieved by using a Nexera XS inert system with a Shim-pack Scepter Claris column to inhibit oligonucleic acid adsorption to metals. Using LabSolutions MD to search for separation parameters enabled a comprehensive evaluation of a variety of LC parameter settings, such as the mobile phase, oven temperature, and gradient program.



MPA Nick Name	MPB Nick Name	Evaluation Value
100mM HFIP 10mM TEA	ACN 50%_MeOH 50%	54.074
100mM HFIP 15mM TEA	ACN 50%_MeOH 50%	53.77
100mM HFIP 5mM TEA	ACN 50%_MeOH 50%	52.477
100mM HFIP 20mM TEA	ACN 50%_MeOH 50%	51.919
200mM HFIP 20mM TEA	ACN 50%_MeOH 50%	47.016
200mM HFIP 15mM TEA	ACN 50%_MeOH 50%	46.926
200mM HFIP 10mM TEA	ACN 50%_MeOH 50%	46.836
200mM HFIP 5mM TEA	ACN 50%_MeOH 50%	45.719
100mM HFIP 10mM TEA	ACN 100%_MeOH 0%	38.822
200mM HFIP 10mM TEA	ACN 100%_MeOH 0%	37.732

Parameter Settings Ranked by Evaluation Score  
(Top 10 settings listed in descending order of evaluation scores)


Application



#### Example 2:

### Corrosion-Resistant System Assists Method Development for Analyzing Antibody Drugs

Antibody drugs made with monoclonal antibodies (mAb) can form dimers, multimers, or other aggregates depending on the manufacturing processes or storage conditions used. Therefore, it is important to monitor the aggregates by size exclusion chromatography during manufacturing processes. Using a corrosion-resistant Nexera XS inert system for that analysis enables reliable data acquisition even for mobile phases with high salt concentrations. The visualization of the design space with LabSolutions MD streamlines the process of identifying the analytical conditions, resulting in ideal peak separation of aggregates / monomers / fragments.



Application



Nexera, Shim-pack, Analytical Intelligence logo, i-Peak Finder and LabSolutions Insight are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation

[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.