

Comprehensive Two-Dimensional Liquid Chromatograph

# Nexera-e



# Nexera™-e


Comprehensive Two-Dimensional Liquid Chromatograph

## Excellent separation method for complex samples

Comprehensive 2D-LC methodology is a paradigm shift in liquid chromatography separation. By combining two independent separation modes orthogonally in combination with a dual-loop/dual-valve alternate switching design, the highest possible peak capacity is achieved.

The new Nexera-e is the ideal solution for the separation and characterization of most challenging and complex samples. It is particularly useful for the analysis of pharmaceutical impurities, proteolytic digests, food extracts, natural products and synthetic polymers.






Comprehensive Two-Dimensional Liquid Chromatograph



# Separation technology for analyzing complex components


## ■ Maximizing the full potential of two separation modes to perform comprehensive two-dimensional separation

Comprehensive 2D-LC is an analytical methodology combining two independent separation modes orthogonally, greatly increasing separation efficiency. The combination of different modes enables the separation of peaks that are difficult to separate using conventional LC, providing excellent results for the analysis of complex sample matrices.



Differences between a conventional 1D-LC and a comprehensive 2D-LC


In comprehensive 2D-LC, all eluates from the first dimension can also be analyzed in the second dimension. The eluates from the first dimension are injected on to the column in the second dimension using two sampling loops that are alternately repeated by a switching valve.



Flow line and mechanism of the comprehensive 2D-LC system

## Improved resolving power compared to conventional 1D-LC

How many peaks are co-eluted under one peak? The more complex the sample, along with a high percentage of similar compounds, there exists a high probability of co-eluting compounds found under single peaks. Even though the co-eluting compounds may not be separated by mono-dimensional separation, an orthogonal two-dimensional system will provide the best possible result.




Chromatogram of a traditional Chinese medicine sample using Nexera-e

## Simultaneous detection in positive and negative mode

Shimadzu triple quadrupole mass spectrometers (TQMS) have the world-class polarity reversal speed (5 ms) and data capture speed (1.5 ms). In comprehensive 2D-LC analysis, the analysis time of second dimension is short, therefore the high-speed performance of the LC-MS affects the quality of the data. Thanks to the high-speed analytical performance of the Shimadzu TQMS, it was possible to obtain stable data even for high-speed analysis of multiple component in comprehensive 2D-LC.

The figures below show contour plotting data obtained from simultaneous analysis of 10 polyphenol components using Nexera-e with shimadzu TQMS. Although some components are detected in positive mode and others in negative mode, it is possible to acquire both data in a single analysis using the Shimadzu TQMS.



Simultaneous analysis of polyphenols

(1) gallic acid, (2) catechin, (3) hesperidin, (4) quercetin-3-o-glucoside, (5) isorhamnetin-3-o-glucoside, (6) sinapic acid, (7) quercetin, (8) naringenin, (9) kaempferol, (10) apigenin

# Nexera Series




Comprehensive 2D-LC by the Shimadzu Nexera series enables comprehensive data collection by analyzing at an extremely low flow rate in the first dimension and ultra-high speed in the second dimension.

The Nexera is equipped with the "Auto-Diagnostic" function, which monitors baseline changes and pressure fluctuations to check for abnormalities. When it detects an unusual fluctuation, it can automatically pause the analysis.

This feature prevents the loss of valuable samples due to abnormal pumping.



Nexera-e (PDA Model)



Auto-Diagnostics Algorithm

## Auto-Diagnostics for Trapped Bubbles (Patent Pending)

If an air bubble becomes trapped in the pump head, it causes an abrupt drop in pressure, after which periodic pressure changes (pulsations) will occur. This can be detected using Shimadzu's proprietary bubble detection algorithm which assesses these distinctive pressure fluctuations caused by trapped bubbles to distinguish them from any expected changes in pressure.

# LCMS-8060NX: suitable for analyzing complex samples


The triple quadrupole mass spectrometer LCMS-8060NX is highly robust due to its superior ionization and ion transport sections. This enables stable analysis of complex samples that tend to contaminate instruments.

## ■ IonFocus™ unit

A newly-developed ESI probe with focus electrodes introduces ions into the mass spectrometer more efficiently, while expelling contaminants to reduce noise and provide more stable data. In addition, an improved heat-assist design promotes the ionization of a wide range of compounds. (Patented technology.)




Simple design without cables or tubes. Remove or replace easily with the "one-touch" lever.



Heating gas with improved heat transfer efficiency

Contaminants

Ions are introduced more efficiently



## ■ UF-Qarray™ II, UF-Lens™ II


Re-engineered ion guide improves robustness without compromising on sensitivity. Maintenance of UF-Qarray II and UF-Lens II can be performed easily without tools.



# Easy setting of 1D and 2D analysis parameters

Comprehensive 2D-LC requires the setting of analytical parameters for the first and second dimensions, respectively. However, setting up the first dimension for ultra-low flow rate analysis and the second dimension for ultra-high-speed analysis can be complicated.


Nexera-e has a dedicated support function to facilitate these settings. The analytical conditions for the first dimension can be set up in the same way as for normal analytical LC, and the gradient settings for the second dimension are automatically generated based on the modulation time.



## ChromSquare, the data analysis software for 2D-LC

### 2D qualitative and quantitative analysis using Contour Graphics

Acquired data is converted to two-dimensional contour plotting using ChromSquare, the software for comprehensive 2D-LC analysis. A peak on a chromatogram is recognized as a spot on the contour plot. The qualitative and quantitative data processing are performed for the target spot.



Two-dimensional chromatogram data processing screen

ChromSquare is a product of Chromaleont S.r.l., Italy.





**ANALYTICAL  
INTELLIGENCE**

Automated support functions utilizing digital technology, such as M2M, IoT, and Artificial Intelligence (AI), that enable higher productivity and maximum reliability. Allows a system to monitor and diagnose itself, handle any issues during data acquisition without user input, and automatically behave as if it were operated by an expert. Supports the acquisition of high quality, reproducible data regardless of an operator's skill level for both routine and demanding applications.

Nexera, Analytical Intelligence logo, IonFocus, UF-Qarray and UF-Lens are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.  
ChromSquare is a trademark of Chromalont S.r.l.



**Shimadzu Corporation**

**[www.shimadzu.com/an/](http://www.shimadzu.com/an/)**

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.