

### Pittcon 2016 1780-4

William Hedgepeth, Ken Tanaka Shimadzu Scientific Instruments, Inc., Columbia, Maryland

PO-CON1613E




## Introduction

Current analytical methodology typically consists of a separate off-line sample preparation technique that is followed by a chromatographic analysis. Sample preparation and manual transfer to the analytical instrument often consumes a majority of the analyst's time and effort. Recently, an innovative new concept was introduced that greatly reduces sample preparation times and the variability associated with manual procedures. This new technique automates the sample preparation and analysis of samples by supercritical fluid extraction of compounds from the sample matrix, which are then transported to the analytical column for analysis without any human intervention. A number of applications that include food, environmental, and pharmaceutical areas will be shown that show the flexibility of this technique.



## SFC Review

- Supercritical CO<sub>2</sub> is a fluid state of carbon dioxide where it is held at or above its critical temperature (31.1 °C) and critical pressure (1,070 psi)
- Rule-of-thumb: any molecule that dissolves in methanol (or less polar solvents) is compatible with SFC
- CO<sub>2</sub> at its critical point is non-polar, solvent strength is increased by using a polar co-solvent
- The benefits of supercritical fluids are still retained with the addition of co-solvents that greatly expand the application range. (subcritical fluid chromatography)



#### 

#### A New Standard in Analytical Workflow Design


## SFC Advantages vs. LC and GC

Lower viscosity of the mobile phase

- faster analysis (5 to 10 x faster)
- less pressure drop across the column
- better sensitivity

Greater separation efficiency than that of HPLC (chiral) "Green" technique – reduced organic solvent usage/waste Little residual solvent obtained in preparative chromatography Analysis of non-volatile, polar or adsorptive solutes without derivatization Analysis of thermally labile compounds Preparative chromatography

### System Configuration – SFE-SFC-MS/MS Online SFE – SFC – MS System




extraction vessel

## Splitless BPR Improves Sensitivity

#### • High sensitivity detection

- Improved sensitivity due to low dead volume BPR
- Splitless injection into the MS without band broadening



# Streamlined, On-line Sample Prep

• SFE (supercritical fluid extraction) sample preparation



Patented absorbent effective for dehydration of samples with high water content



Up to 48 samples can be extracted and analyzed in an automated workflow using the Rack Changer.

approx. 30 min time saving / sample




Conventional sample preparation method (QuEChERS)


| Homogenize Add solvent Stir Add reagent Stir Centrifuge Transfer Add reagent Stir Centrifuge | • |
|--|---|
|--|---|

## Preserves Labile Compounds

### Analysis of coenzyme Q<sub>10</sub> in supplements

• Conventional Method: Coenzyme  $Q_{10}$  is prone to oxidation






• Nexera UC: No oxidation with SFE

Extraction with 5 % MeOH in CO<sub>2</sub> for 4 min

Gradient: 5 %B (0 min – 4 min for SFE) 5 % - 50 %B in 5 min



## DBS Analysis


### Analysis of biomarkers from dried blood spot (DBS)

Conventional Method: time consuming

| Spot the sample | Cut DBS | Add extraction<br>solvent | Stir / Shake | Soak for 1 h | Filtrate | Evaporate /<br>reconstitute |  |
|-----------------|---------|---------------------------|--------------|--------------|----------|-----------------------------|--|

• Nexera UC:






Extraction and analysis of 1 ppm phospholipid spiked into plasma on DBS


# Polymer Additives Analysis

### Extraction of polymer additives

- Conventional Method: 10 20 h soxhlet extraction
- Nexera UC: 7 min CO<sub>2</sub> extraction
  - Crushed polymer samples are placed in the extraction vessel
  - After 7 min CO<sub>2</sub> extraction, the sample is ready for SFC analysis







#### 


A New Standard in Analytical Workflow Design

# Cleaning Analyses


### **Cleaning validation**

- Conventional Method:
  - Sampling Solvent extraction Concentration HPLC analysis
  - Organic solvents can't be used for TOC analysis
- Nexera UC:
  - Swab is enclosed in the extraction vessel for SFE





# Analysis of API in a Tablet



| Sample        | Area      | Area converted to 2mg of sample | Recovery (%) |
|---------------|-----------|---------------------------------|--------------|
| Tablet 2.27mg | 4,508,728 | 3972447.577                     | 100.1        |
| Tablet 2.03mg | 3,954,148 | 3895712.315                     | 98.2         |
| Tablet 2.10mg | 4,378,007 | 4169530.476                     | 105.1        |
| Tablet 1.96mg | 3,995,147 | 4076680.612                     | 102.7        |
| Tablet 2.25mg | 4,473,899 | 3976799.111                     | 100.2        |

Approx. 2mg of crushed tablet sample was weighed and transferred to the vessel for analysis. Recovery was calculated by comparing to the standard that is equivalent to the amount of API in 2mg of the tablet.



### Summary

- A new automated online sample preparation/chromatography system that uses supercritical fluids was recently introduced.
- The system can greatly reduce sample preparation times when compared to conventional methods like QuEChERS.
- Automated supercritical fluid extraction (SFE) can reduce solvent waste while improving reproducibility of results.
- SFE was shown to be a viable technique for a number of applications.

First Edition: March, 2016



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names of the names of the

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Shimadzu Corporation www.shimadzu.com/an/