() SHIMADZU

Simultaneous determination of four immunosuppressants in whole blood – A fast method on a Shimadzu tandem LCMS-8030 mass analyser

Katharina Rosing²; Manfred Fobker²; Frank Kannenberg²; Joerg Stypmann³; Klaus Dr. Bolliq¹; Sven Vedder⁴; Roch Nofer²

1Shimadzu Germany GmbH, Duisburg, Germany; 2Center for Laboratory Medicine, University Hospital Muenster, Germany; 3Department for cardiology and angiology, University Hospital Muenster, Germany; 4Shimadzu Europe GmbH, Duisburg, Germany

1: Introduction

The amount of different kind of organ transplantations is increasing steadily since the improvement of this technology in the past century. The European organization Eurotransplant counts several thousand surgical interventions in Europe every year. To avoid rejection of the transplanted organ by the immune system, a complex therapy of immunosuppressants in combination with other drugs is necessary. A therapeutic concentration range in blood minimizes the risk of unwanted side effects. Immunoassays which are an alternative method for monitoring immunosuppressants, often lead to insufficient results due to their cross reactivity to drug metabolites. In this study we established a fast, kit-based screening method for routine use of immunosuppressants based on LCMS/MS technology.

2: Materials and Methods

The four immunosuppressants Cyclosporine A. Tacrolimus, Sirolimus and Everolimus are often used in therapy followed by autologous organ transplantation. These drugs were analyzed in whole blood samples using a commercially available kit supplied by Chromsystems®2, München, Germany. The samples are first prepared by a short manual protein precipitation, injected into the Shimadzu LCMS-8030 triple-guad system and afterwards analyzed by LabSolution software. An online trapping system separates the analytes from

LC/MS/MS:

Key feature of LCMS-8030 triple quadrupole mass spectrometer

matrix signals and eliminates unwanted effects like ion suppression.

- Ultra fast polarity switching of 15msec
- · Ultra fast scan speed of up to 15,000 u/sec
- UFsweeper™ technology dramatically minimizes cross talk
- · Excellent linearity with wide dynamic range



Figure 1: LCMS-8030 triple quadrupole mass spectrometer

3: Results

3-1 Method Validation Linearity

Linearity was evaluated by analysis of everolimus standards over the concentration range of the calibration curve (Table 1). Correlation coefficient was 0,99956 for everolimus, 0,99916 for tacrolimus, 0,99944 for sirolimus and 0,99973 for cyclosporine a respectively.

Range (µg/l)	Lower limit of quantification (µg/l)
25-900	3.9
2-40	0.395
2-40	1.003
2-50	0.997
	Range (μg/l) 25-900 2-40 2-40 2-50

Table 1. Range and LOQ.

Precision Interday- and Intraday-Precision were carried out by evaluation of quality control samples prepared from the same specimen(n=5) and measured at different times of the day(intraday) or different days (interday) with the same instrument parameters. Inter- and intraday precision were assessed by daily routine samples of hospital patients in replicates prepared for each method. Interday-precison for low-quality-control-samples are shown in table 2. Interday precision was also obtained by quality control samples of 4 different levels for all four immuno-suppressants supplied by Chromsystems. Stability of these control samples were monitored by measurement of the same batch/lot at different times during the course of 14 days (Table 3).

Drug	CV
CsA	1.4%
Everolimus	6.4%
Tacrolimus	5.2%
Sirolimus	9.4 %

Table 2. CV for interday precision of QK1.





Accuracy

Accuracy limits were set within the following criteria: 90-110% and imprecision CV of less than 10%. LC-MS/MS within run accuracy were determined by analyzing 4 quality controls at different levels (therapeutic range). The accuracy biases of the quality controls are shown in Table 4:

Drug	AB	AB	AB	AB
	QC1	QC2	QC3	QC4
CsA	-0,5%	2.13%	1.8%	6.8%
Everolimus	3.01%	2.28%	2.33%	-1.11%
Tacrolimus	-5.62%	-2.72%	4.1%	9.3%
Sirolimus	-7.32%	2.77%	-0.37	4.43%

3-2 Passing-Bablok analysis Method comparison

Comparison of methods is an important process in validation of a new analytical measurement brought into laboratory routine. In this case it was evaluated using Passing-Bablok-analysis to confirm a distinct relation³ between either HPLC-MS vs. HPLC-MS/MS (everolimus) or Immunoassay vs. HPLC-MS/MS (csA, tacrolimus, sirolimus). All four regression curves show good correlation. Slope: 1,02 (sirolimus), 0,99 (tacrolimus), 0,97 (everolimus) and 1,03 for csA. 25-69 samples were examined under the same conditions over a course of 2 weeks to give a spectrum of repeatability under the given circumstances.

Drug	Method comparison	
CsA	Immunoassay	
Everolimus	HPLC-MS	
Tacrolimus	Immunoassay	
Sirolimus	Immunoassay	

regression Bablok I Passing | 95% CI



Пкм

WESTERINGON

Figure 2. Passing-Bablok-Regressions

4: Conclusion

Quantitative Assessment of immunosuppressant drug levels is, though extremely important to patients needs and sufficient organ survival in transplantation medicine, a challenging part of clinical laboratory procedure⁴. Although there have been many approaches to simplify HPLC-MS by immunoassay methods it is still the most reliable method and gold standard procedure in modern clinical monitoring of immuno-suppressant therapy.

5: References

1.Vogeser M, Spöhrer U. Automated processing of whole blood samples for the determination of immunosuppressants by liquid chromatography tandem-mass spectrometry. Clin. Chem. Lab. Med. 2006;44(9):1126-1130. 2.Chromsystems application manual, MassTox® Immunsuppressiva, 2010 3.Passing H. Bablok, A new biometrical procedure for testing the equality of measurements from two different analytical methods. Application of linear regression procedures for method comparison studies in clinical chemistry, Part I. J. Clin. Chem. Clin. Biochem. 1983;21(11):709-720. 4. Vogeser M. Seger C. A decade of HPLC-MS/MS in the routine clinical laboratory--goals for further developments. Clin. Biochem. 2008;41(9):649-662.