

Application of a Sensitive Liquid Chromatography-Tandem Mass Spectrometric Method to Pharmacokinetic Study of Telbivudine in Humans

ASMS 2014 WP 629

Bicui Chen¹, Bin Wang¹, Xiaojin Shi¹, Yuling Song²,
Jinting Yao², Taohong Huang², Shin-ichi Kawano²,
Yuki Hashi²

1 Pharmacy Department, Huashan Hospital,
Fudan University,

2 Shimadzu Global COE, Shimadzu (China) Co., Ltd.

Application of a Sensitive Liquid Chromatography-Tandem Mass Spectrometric Method to Pharmacokinetic Study of Telbivudine in Humans

Introduction

Telbivudine is a synthetic L-nucleoside analogue, which is phosphorylated to its active metabolite, 5'-triphosphate, by cellular kinases. The telbivudine 5'-triphosphate inhibits HBV DNA polymerase (a reverse transcriptase) by competing with the natural substrate, dTTP. Incorporation

of 5'-triphosphorylated-telbivudine into viral DNA obligates DNA chain termination, resulting in inhibition of HBV replication. The objectives of the current studies were to develop a selective and sensitive LC-MS/MS method to determine of telbivudine in human plasma.

Method

Sample Preparation

- (1) Add 100 μ L of plasma into the polypropylene tube, add 40 μ L of internal standard working solution (33 μ g/mL, with thymidine phosphorylase) to all other tubes.
- (2) Incubate the tubes for 1 h at 37 °C in dark.
- (3) Add 200 μ L of acetonitrile to all tubes, seal and vortex for 1 minutes.
- (4) Centrifuge the tubes for 5 minutes at 13000 rpm.
- (5) Transfer 200 μ L supernatant to a clean glass bottle and inject into the HPLC-MS/MS system.

LC-MS/MS Analysis

The analysis was performed on a Shimadzu Nexera UHPLC instrument (Kyoto, Japan) equipped with LC-30AD pumps, CTO-30A column oven, DGU-30A₅ on-line degasser, and SIL-30AC autosampler. The separation was carried out on GL Sciences InertSustain C18 column (3.0 mmI.D. x 100

mmL.) with the column temperature at 40 °C. A triple quadruple mass spectrometer (Shimadzu LCMS-8050, Kyoto, Japan) was connected to the UHPLC instrument via an ESI interface.

Analytical Conditions

HPLC (Nexera UHPLC system)

Column	: InertSustain (3.0 mmI.D. x 100 mmL., 2 μ m, GL Sciences)
Mobile Phase A	: water with 0.1% formic acid
Mobile Phase B	: acetonitrile
Gradient Program	: as shown in Table 1
Flow Rate	: 0.4 mL/min
Column Temperature	: 40 °C
Injection Volume	: 2 μ L

Table 1 Time Program

Time (min)	Module	Command	Value
0.00	Pumps	Pump B Conc.	5
4.00	Pumps	Pump B Conc.	80
4.10	Pumps	Pump B Conc.	5
6.00	Controller	Stop	

Application of a Sensitive Liquid Chromatography-Tandem Mass Spectrometric Method to Pharmacokinetic Study of Telbivudine in Humans

MS (LCMS-8050 triple quadrupole mass spectrometer)

Ionization	: ESI
Polarity	: Positive
Ionization Voltage	: +0.5 kV (ESI-Positive mode)
Nebulizing Gas Flow	: 3.0 L/min
Heating Gas Flow	: 8.0 L/min
Drying Gas Flow	: 12.0 L/min
Interface Temperature	: 250 °C
Heat Block Temperature	: 300 °C
DL Temperature	: 350 °C
Mode	: MRM

Table 2 MRM Parameters

Compound	Precursor <i>m/z</i>	Product <i>m/z</i>	Dwell Time (ms)	Q1 Pre Bias (V)	CE (V)	Q3 Pre Bias (V)
Telbivudine	243.10	127.10	100	-26	-10	-13
Telbivudine-D3	246.10	130.10	100	-16	-9	-25

Results and Discussion

Human plasma samples containing telbivudine ranging from 1.0 to 10000 ng/mL were prepared and extracted by protein precipitation and the final extracts were analyzed by LC-MS/MS. MRM chromatograms of telbivudine (1 ng/mL) and deuterated internal standard are presented in Fig. 1 (blank) and Fig. 2 (spiked). The linear regression for telbivudine was found to be >0.9999. The calibration curve with human plasma as the matrix were shown in Fig. 3. Excellent precision and accuracy were maintained for four orders of magnitude, demonstrating a linear dynamic range suitable for real-world applications. LLOQ for telbivudine was 1.0 ng/mL, which met the criteria for bias (%) and precision within ±15% both within run and between run. The

intra-day and inter-day precision and accuracy of the assay were investigated by analyzing QC samples. Intra-day precision (%RSD) at three concentration levels (3, 30, and 8000 ng/mL) were below 2.5% and inter-day precision (%RSD) was below 2.5%. The recoveries of telbivudine were 100.6±2.5 %, 104.5±1.5% and 104.3±1.6% at three concentration levels, respectively. The stability data showed that the processed samples were stable at the room temperature for 8 h, and there was no significant degradation during the three freeze/thaw cycles at -20 °C. The reinjection reproducibility results indicated that the extracted samples could be stable for 72 h at 10 °C.

Application of a Sensitive Liquid Chromatography-Tandem Mass Spectrometric Method to Pharmacokinetic Study of Telbivudine in Humans




Figure 1 Representative MRM chromatograms of blank human plasma (left: transition for telbivudine, right: transition for internal standard)




Figure 2 Representative MRM chromatograms of telbivudine (left, 1 ng/mL) and internal standard (right) in human plasma




Figure 3 Calibration curve of telbivudine in human plasma

Application of a Sensitive Liquid Chromatography-Tandem Mass Spectrometric Method to Pharmacokinetic Study of Telbivudine in Humans

Compound	Calibration Curve	Linear Range (ng/mL)	Accuracy (%)	r
Telbivudine	$Y = (2.77 \times 10^{-4})X + (3.39 \times 10^{-5})$	1~10000	93.1~116.6%	0.9998

Table 3 Accuracy and precision for the analysis of amlodipine in human plasma (in pre-study validation, n=3 days, six replicates per day)

Added Conc. (ng/mL)	Intra-day Precision (%RSD)	Inter-day Precision (%RSD)	Accuracy (%)
3	2.18	2.11	107.7~114.4
400	1.52	1.58	91.6~95.9
8000	1.76	1.68	95.4~101.3

Table 4 Recovery for QC samples (n=6)

QC Level	Concentration (ng/mL)	Recovery (%)
LQC	3	100.6
MQC	400	104.5
HQC	8000	104.3

Table 5 Matrix effect for QC samples (n=6)

QC Level	Added Conc. (ng/mL)	Matrix Factor	IS-Normalized Matrix Factor
LQC	3	82.3%	99.0%
MQC	400	81.7%	101.0%
HQC	8000	90.8%	101.5%




Figure 4 Representative MRM chromatograms of real-world sample

Application of a Sensitive Liquid Chromatography-Tandem Mass Spectrometric Method to Pharmacokinetic Study of Telbivudine in Humans

Conclusion

Results of parameters for method validation such as dynamic range, linearity, LLOQ, intra-day precision, inter-day precision, recoveries, and matrix effect factors were excellent. The sensitive LC-MS/MS technique provides a powerful tool for the high-throughput and highly selective analysis of telbivudine in clinical trial study.