

Rapid Quantification of Perfluorinated Compounds in Drinking and Surface Water Using LC-MS/MS



Jeremy Post, Christopher Gilles, William Lipps Shimadzu Scientific Instruments, Columbia, MD

PO-CON1566E

Rapid Quantification of Perfluorinated Compounds in Drinking and Surface Water Using LC-MS/MS

Novel Aspect

Fast, low chemical background, high sensitivity measurement of perfluorinated compounds in water.

Introduction

Perfluorinated compounds (PFCs) are synthetic molecules with carbon chains to which fluorine is bound by high strength carbon-fluorine covalent bonds. Useful due to their surfactant and hydrophobic physicochemical properties these materials are used broadly since the 1950s in industrial applications such as the conveyance of stain and stick resistance to a wide variety of materials and products including clothing, cookware, carpet, and food containers. These molecules, now classified as persistent organic pollutants, are highly stable and resistant to biological degradation, therefore they persist in the environment, bio-accumulate, and due to their ready transport by water are found even in remote regions of Earth. As contaminants of emerging concern, research on PFCs is ongoing to determine the impacts of these materials on human health and the environment. A CDC study detected four of these analytes in >98% of serum samples representative of the general US population above age twelve. Low level analysis of these compounds is challenging due to their ubiquitous presence in the laboratory and analytical instrumentation components.



Figure 1: Central Maryland surface water source locations



Figure 2: Direct method surface water sampling collection technique

Rapid Quantification of Perfluorinated Compounds in Drinking and Surface Water Using LC-MS/MS

Method

Surface water samples were collected and stored according to document SOP #EH-01 Surface Water Collection using the Direct Method. As the mobile phase transfer lines, degasser membrane, pump inlet block, and sample vial septa all utilize PTFE, made from a prominent PFC, perfluorooctanoic acid (PFOA), these components were either exchanged or bypassed to eliminate all sources of extraneous analyte from the instrumentation. The calibration curve solutions were prepared in glass vials using glass syringes with a metal plunger. PFC standards were purchased from Wellington Laboratories. A Shimadzu LCMS-8050 was used for analyte separation and mass spectrometry. Analytes were ionized using negative mode electrospray ionization. All samples were passed through 0.2 µm glass fiber filter prior to analysis. The method required only 3.15 minutes. Table 1 provides optimized transitions and retention times for the PFCs measured in this method.

Chromatography Parameters					
Flow Rate	: 0.35 mL/min				
Column oven	: 35 °C				
Injection volume	: 20 µL				
Mobile Phase A	: 10 mM ammonium acetate in water				
Mobile Phase B	: 10 mM ammonium acetate in 80:20 MeOH:ACN				
Column	: Supelco Titan C18 (50 x 2.1mm 1.9 μm)				
Gradient	: Linear 45% - 97% B 0-2.6 minutes				
Mass Spectrometry Parameters					
Mass Spectrometry Faramet					
 Nebulizing Gas	: 2.2 L/min				
Nebulizing Gas Heating Gas	: 2.2 L/min : 8 L/min				
Nebulizing Gas Heating Gas Drying Gas	: 2.2 L/min : 8 L/min : 10 L/min				
Nebulizing Gas Heating Gas Drying Gas Interface Temperature	: 2.2 L/min : 8 L/min : 10 L/min : 230 °C				
Nebulizing Gas Heating Gas Drying Gas Interface Temperature Desolvation Line Temperature	: 2.2 L/min : 8 L/min : 10 L/min : 230 °C : 95 °C				
Nebulizing Gas Heating Gas Drying Gas Interface Temperature Desolvation Line Temperature Heat Block Temperature	: 2.2 L/min : 8 L/min : 10 L/min : 230 °C : 95 °C : 225 °C				
Nebulizing Gas Heating Gas Drying Gas Interface Temperature Desolvation Line Temperature Heat Block Temperature Column Oven	: 2.2 L/min : 8 L/min : 10 L/min : 230 °C : 95 °C : 225 °C : 35 °C				

Table 1: Method information for the perfluorinated compounds measured in surface and drinking water

Name	Abbreviation	Structure	Retention Time	Transitions
Perfluoro-n-hexanoic acid	PFHxA		1.28	313>119 313>269
Perfluoro-n-heptanoic acid	PFHpA		1.63	363>169 363>319
Perfluoro-n-octanoic acid	PFOA	E C C C C C C C C C C C C C C C C C C C	1.91	413>169 413>369
Perfluoro-n-nonanoic acid	PFNA	, , , , , , , , , , , , , , , , , , , ,	2.15	463>219 463>419
Perfluoro-n-decanoic acid	PFDA		2.36	513>269 513>469
K Perfluoro-n-butanesulfonate	L-PFBS	F F F F F S03K*	1.02	299>80 299>99
Na Perfluoro-n-hexanesulfonate	L-PFHxS		1.70	399>80 399>99
Na Perfluoro-n-octanesulfonate	L-PFOS		2.18	499>80 499>99

Excellence in Science

Rapid Quantification of Perfluorinated Compounds in Drinking and Surface Water Using LC-MS/MS

Results



Figure 3: PFC calibration curves demonstrating excellent linearity at low analyte concentration matched with spectra showing the analyte LOQ and a chromatogram with 60 fg on column.

Excellence in Science

Rapid Quantification of Perfluorinated Compounds in Drinking and Surface Water Using LC-MS/MS



Figure 4: Representative PFC LC Chromatograms from 20uL injections of surface and drinking water.

Table 2:Calculated amounts of PFC materials in the surface and drinking water samples.
The measured amounts were averaged (n=4) and reduced by the amount of
each analyte detected in the Lab Blank (18MOhm purified water).
These values were then divided by the injection volume to determine the final
concentration in ng/L.

	Lab Blank	Centennial Lake	Liberty Lake	Little Patuxent	Middle Patuxent	Drinking Water
				River	River	
PFHxA	0.34	1.19	0.59	0.68	1.17	0.51
PFHpA	0.1	0.36	0.33	0.33	0.33	0.16
PFOA	0.79	1.09	0.92	1.04	1.05	0.37
PFNA	0.22	0.24	0.24	0.25	0.26	0
PFDA	0.71	0.43	0.16	0.07	0.08	0
L-PFBS	0	0	0	0	0	0
L-PFHxS	0	0	0	0.06	0.06	0.05
PFOS	0	0	0.13	0.19	0.1	0



Rapid Quantification of Perfluorinated Compounds in Drinking and Surface Water Using LC-MS/MS

Discussion

Creating an environment suitable for measuring PFCs required several adjustments to the standard LCMS workflow in both the way samples were prepared and data was acquired. Iterative modifications started with tubing replacements and LC hardware reconfiguration, which was followed by using glass and metal only for sample preparation. Even a seemingly minor adjustment such as inverting the autosampler vial septa proved important to reducing the PFC chemical background. Following these modifications the background reduced significantly making possible the acquisition of calibration curves measuring down to 2.5 fg on column. The low standard %RSD value for PFOA was 4.1 and for PFOS, 5.6. Seven of the eight compounds were detected in surface waters across central Maryland whereas four were found in drinking water, though in reduced amounts relative to surface water.

Conclusions

Though nearly ubiquitous, perfluorinated compounds can be effectively and quickly measured directly from surface and drinking water using a modified configuration of the LCMS-8050. This method can be used for routine monitoring of PFCs in water samples requiring only minimal sample preparation.

First Edition: June, 2015



Shimadzu Corporation www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2015