



Motoshi Sakakura<sup>1</sup>, Teruhisa Shiota<sup>1</sup>, Jun Watanabe<sup>2</sup> <sup>1</sup> AMR, Inc., Tokyo, JAPAN;

<sup>2</sup> Shimadzu Corporation, Kyoto, JAPAN



PO-CON1641E



## Introduction

Acylglycerols included in a food play the role important to the food function such as flavor and nutrition. There were some problems such as complex preparation steps and carry over to analytical instruments with an analysis of acylglycerols, and it was difficult to analyze a great deal of samples quickly.



In this study, we analyzed food samples including acylglyserols directly with DART (Direct Analysis in Real Time) and mass spectrometer with less sample pretreatment and tried to detect the characteristic ingredients in samples and distinguish each sample.



## Methods and Materials

Mass spectrometer LCMS-2020 (Shimadzu Corporation, Kyoto, Japan) equipped with real time direct analysis ion source DART-OS (IonSense, Inc., Saugus, MA, USA) was used for this study. Marketed soybean milk, butter and margarine, etc. were used for food samples including acylglycerols. Samples were analyzed by making a small amount samples adherent onto a glass stick and holding it up to the DART gas stream about 10 seconds. Further, without pretreating a sample, it was set to the analysis. LCMS-2020 can achieve the polarity switching time of 15msec and the scanning speed of up to 15000u/sec, therefore the loop time can be set at less than 1 second despite the relatively large scanning range of 10-1000u.





Figure 1 DART-MS System; LCMS-2020 equipped with DART-OS



# Result

Various samples were analyzed by DART-MS. Analyzing time of each sample was about 10 seconds and analyzing interval of each sample was less than 1 minute.

### **Positive Mass Spectum**

From the results of food samples of plant origin like soybean milk and margarine, several signals of triacylglycerides (m/z 900 near here) and diacylglycerides origin mainly comprised of oleic acid and the palmitic acid were detected in the positive ion spectrum. In the case of measuring food samples of animal origin like milk and butter, many signals of triacylglycerides and diacylglycerides which comprised of middle chain fatty acid as well as long chain fatty acid were detected, and the spectrum pattern was indicated different from vegetable fats and oils.



Figure 2 Positive mass spectra of samples analyzed by DART-MS



The marketed food sample which was margarine added butter was analyzed. It had just mixed mass spectrum in which some glyceride signals mainly comprised of oleic acid and the palmitic acid were detected and many glyceride signals which comprised of middle and long chain fatty acids were detected.



Figure 3 Positive mass spectum of margarine added butter analyzed by DART-MS

#### **Negative Mass Spectrum**

The negative mass spectra were indicated (figure 4). A signal was detected by m/z 215 in all spectra except for soybean milk both. It was inferred that these signals were chloric adduct ion [M+CI]- of monosaccharide (C6H12O6).



Figure 4 Positive mass spectra of samples analyzed by DART-MS

Total ion current chromatograms (TIC) and extracted ion current chromatograms (XIC) of m/z 215 in negative mode when each sample being analyzed were shown in figure 5. XIC signal of monosaccharide origin wasn't detected in soybean milk. It was detected weakly in butter and hard in milk and margarine.



Figure 5 TIC chromatogram and XIC chromatogram (m/z 215)

High throughput lipid analysis in foods containing oils and fats was achieved by a mass spectrometer with high speed polarity switching and high speed scanning integrated with DART ion source which could real-time-analyze and it was possible to profile plant acylglycerol and animal acylglycerol.

## Conclusions

#### Various food samples were analyzed by DART-MS;

High throughput lipid analysis in foods containing oils and fats was achieved by a mass spectrometer with high speed polarity switching and high speed scanning.

#### From the results of food samples;

[Food samples of plant origin] several signals of triacylglycerides and diacylglycerides origin mainly comprised of oleic acid and the palmitic acid were detected.

[Food samples of animal origin] many signals of triacylglycerides and diacylglycerides which comprised of middle and long chain fatty acids were detected.



### References

- 1. A. McANOY et.al, J. Am. Mass Spectrom., 2005, 16, 1498-1509.
- 2. B. Winther et.al, Lipids, 2011, 46, 25-36.
- 3. V. A.Tyurin et.al, Journal of Neurochemistry, 2008, 107, 1614-1633.
- 4. H. Mikuma, A. Kaneko, BUNSEKI KAGAKU, 2010, 59(5), 399-404.

Disclaimer: The products and applications in this presentations are intended for Research Use Only (RUO). Not for use in diagnostic procedures.





Shimadzu Corporation

www.shimadzu.com/an/

#### For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "@". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names of ther the networks and trade names of the na

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.