🕀 SHIMADZU

Fast polarity switching and MRM triggered automatic MS/MS applied to benzodiazepines and their metabolites in clinical and forensic analysis

ASMS 2011 WP171

Toshikazu Minohata¹, Natsuyo Asano¹, Kiyomi Arakawa¹, Jun Watanabe¹, Junko Iida¹, Hitoshi Tsuchihashi², Koichi Suzuki², Kei Zaitsu³, Noriaki Shima³, Munehiro Katagi³

- ¹ Shimadzu Corporation, Kyoto, JAPAN
- ² Dep. of Legal Medicine, Osaka Medical Collage, Takatsuki, JAPAN
- ³ Forensic Science Lab, Osaka Prefectural Police, Osaka, JAPAN

Fast polarity switching and MRM triggered automatic MS/MS applied to benzodiazepines and their metabolites in clinical and forensic analysis

Introduction

Benzodiazepines are now among the most commonly-prescribed drugs, which increases their potential for addiction and abuse, and they are often found in combination with other drugs in drug-related fatalities or drug facilitated sexual assault cases. Using fast polarity switching and MRM triggered automatic MS/MS, a new method package was developed for the simultaneous screening and quantification of benzodiazepines and benzodiazepine-like substances, which are available in Japan and that are relevant in clinical and forensic cases. This method, and associated MS/MS data base, will be applied to both low dose and high dose benzodiazepines abuse in forensic analysis.

Methods

A standard sample solution was prepared by mixing 35 benzodiazepines, Zopiclone, Zolpidem and their metabolites in water (Table 1). Human urine samples (1mL) were initially adjusted to pH 5 with acetic acid and extracted into

chloroform / isopropanol (3:1) by vortex mixing. After centrifugation, the organic phase was transferred to a clean tube and evaporated to dryness by a gentle stream of N2. The residue was dissolved in 500uL of water/methanol (9:1).

	Table 1. List of Compound	ds	Tabl	e 2. Analytical conditions
Hypnotic Drugs	M	letabolites	Accelerated Nexera LC	method
Alprazolam	a-Hydroxyalprazolam		Column	Shim-pack XR-ODSⅢ
Bromazepam	3-OH ABBP			(2.0 mml.D. x 50 mmL., 1.6 um)
Brotizolam	a-Hydroxybrotizolam	4-Hhydroxybrotizolam	Mobile phase A	10 mM Ammonium formate
Chlordiazepoxide	Desmethyldiazepam	Oxazepam	Mobile phase B	Methanol
Clobazam	N-desmethylclobazam		Gradient program	5 %B (0 min) - 95 %B (5-6 min):
Clonazepam	7-Aminoclonazepam	7-Acetamidoclonazepam		5 %B (6 01 - 10 min)
Clorazepic acid	Desmethyldiazepam	Oxazepam	Flow rate	0 3ml /min
Clotiazepam	Desmethylclotiazepam	Clotiazepam Y-10247		40 °C
Cloxazolam	Delorazepam			40 C
Delorazepam	-	_	Injection volume	5 UL
Diazepam	Desmethyldiazepam	Oxazepam	Cycle time	10 mins
Estazolam	Estazolam M-II	Estazolam M-IV	Conventional LC condi	tions
Ethyl loflazepate	Desalkylflurazepam		Column	Shim-pack FC-ODS
Etizolam	8-Hydroxyetizolam	a-Hydroxyetizolam		(2.0 mml.D. x 150 mmL., 3 um)
Fludiazepam	Desalkylflurazepam		Gradient	5 %B (0 min) - 95 %B (15-20 min); 5 %B
Flunitrazepam	3-Hydroxyflunitrazepam	7-Aminoflunitrazepam		(20.01 - 30 min)
	7-Acetamidoflunitrazepam		Cycle time	30 mins
Flurazepam	1-Ethanolflurazepam	Desalkylflurazepam		
Flutazolam	1-Ethanolflurazepam		LCIVIS-8030 Iriple quad	
Flutoprazepam	Desalkylflurazepam		Ionization	ESI
Haloxazolam	Haloxazolam KAZ-609		Polarity	Simultaneous positive/negative
Lorazepam	-			switching (15msecs)
Lormetazepam	Lorazepam	0	Gas flows	Nebulising 1.5L/min; Drying 10L/min
Medazepam	Desmethyldiazepam	Oxazepam	Source temperatures	Capillary 250 °C; Block 400 °C
Mexazolam	Delorazepam			
IVIIdazolam	a-Hydroxymidazolam	7 A anto an international and	Sustan design options	
Nimetazepam	7-Aminonimetazepam	7-Acetamidonimetazepam	System design options	Development of the table of the second
Nitrazepam	7-Adminonitrazepam		Nexera UPLC	Pressure range up to 130 MPa high-speed
Oxazepam	- Desmethuldiazonam	Overenem		injection, overlap injection, near-zero carryover,
Dxazolalli	2 Hudrow prozonom	Oxazepani		precise solvent delivery and excellent
Plazepalli	5-Hydroxyprazepam	Oxazepani Ovazepam M 4		reproducibility
Pilmazafono		Quazepain M-4	LCMS-8030	Ultra fast polarity switching of 15
Tomazonam				msec & ultra fast scan speed of up to
Tofisonam				15,000 u/sec
Triazolam	- a-Hydroxytriazolam	4-Hydroxytriazolam		UFsweeper [®] technology dramatically
	a-riyuloxythazolam			minimizes cross talk and delivers
Zoniclone	Zopiclone-N-oxide	N-desmetylzoniclone		excellent linearity with a wide
Zolnidem	Zolpidem-COOH [main]	Zolpidem-COOH [minor]		dynamic range

🕀 SHIMADZU

Fast polarity switching and MRM triggered automatic MS/MS applied to benzodiazepines and their metabolites in clinical and forensic analysis



Results

Benzodiazepine analysis has been accelerated by bringing together high resolution LC separations and high speed data acquisition LC/MS/MS system with fast polarity switching.



The method uses Synchronized Survey Scan[®] (in this mode MS/MS scanning is triggered by MRM signals) generating a full-product ion mass spectrum and MRM data.

Figure 1.

Accelerated benzodiazepine analysis

MRM chromatograms obtained from analysis of a standard solution of 35 benzodiazepines, Zopiclone, Zolpidem and their metabolites (each 100 ng/mL). Using high resolution LC conditions and high speed MS/MS acquisitions the sample analysis time is reduced to 10 minutes compared to 30 minutes for the standard method.

It is also important to note the near zero carry over on sample injection plays a key part in minimizing errors and helping to improve sensitivity and precision.

Figure 2.

Standard benzodiazepine analysis

MRM chromatograms obtained from analysis of a standard solution of 35 benzodiazepines, Zopiclone, Zolpidem and their metabolites (each 100 ng/mL) using a low resolution LC separation Shim-pack FC-ODS. The analysis time was 30 minutes.

🕀 SHIMADZU

Fast polarity switching and MRM triggered automatic MS/MS applied to benzodiazepines and their metabolites in clinical and forensic analysis

MS/MS library matching

Using a polarity switching speed of 15msec and a scan speed of 15,000u/sec, product ion spectra were generated in both positive and negative ionization which could be matched against a user library of 70 compounds an

Туре	Event#	+/-	Compound Name n/z	Time (2.898 min - 5.711 min) 🔺	-
MRM	19	+	66_N-desmetylzopiclone 375.		
- Product Ion Scan	20	+	66_N-desmetylzopiclone 100.	and the second se	_
MRM	21	+	57_7-acetamidonimetazepam		
- Product Ion Scan	22	+	57_7-acetamidonimetazepam		
MRM	143	-	57_7-acetamidonimetazepam		
- Product Ion Scan	144	-	57_7-acetamidonimetazepam		
MRM	23	+	63_4-hydroxytriazolam 359.0		
- Product Ion Scan	24	+	63_4-hydroxytriazolam 100.0		
MRM	25	+	44_clotiazepam Y-10241321.		
- Product Ion Scan	26	+	44_clotiazepam Y-10241100.		-

Figure 3.

Method set-up for accelerated benzodiazepine analysis using overlapping MRM/product ion scan acquisitions in both positive and negative mode (17 of 70 compounds were set both positive and negative mode). automated aid to screening and compound identification. Fast polarity switching helps to provide information rich product ion spectra resulting in better detection and identification for each benzodiazepine.

Ch1 375.20 245.10 250 -14.0 -21.0 -20.0 Ch3 Ch4 -21.0 -20.0 -20.0 -20.0 -20.0 Ch4 Ch4 -20.0 -20.0 -20.0 -20.0 -20.0 -20.0 Ch4 Ch4 -20.0 -2	Pre Bias(V) 🔺	Q3 Pre	CE	as(V)	Q1 Pre Bias	(msec)	Time	m/z Dwell	Product	recursor m/z	Ch
Ch2 Ch3 Ch4 Egent Time: 0.028 sec QQ Resolution: Link ¥ Queen the vert: Product Ion Scan ¥ Queen the vert: V Queen the vert: V	5 –	-20.0	-21.0		-14.0			25.0	245.10	20	Ch1
Crist Eyent Time: 0.028 sec Q1 Resolution: Unit Advance Q3 Resolution: Unit Advance Q4 Resolution: Unit Advance Q5 Resolution:											Ch2
Eyent Time: 0.028 sec Q1 Resolution: Unit Q3 Resolution: Unit Q3 Resolution: Unit											Ch3
Eyent Time: 0.028 sec Q1 Resolution: Unit Imit Imit <t< td=""><td>ЪĒ</td><th>1</th><th></th><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td>() ↓ </td></t<>	ЪĒ	1									() ↓
Start m/z: 50.00 Epd m/z: 400.00 Precursor Ion m/z: 100.00 Collision Energy: -15.0 v	ed Settings	Advanced	it 👤 it 👤 ve <u>y</u> Event Set ct Ion Scan	on: Un on: Un Sur Produ	Q1 Resdution Q3 Resdution urvey Evert	sec Use as Si <u>D</u> eper	0.028	E <u>v</u> ent Time:			
Precursgr Ion m/z: 00.00 Collision Energy: 15.0 v							400.00	E <u>n</u> d m/z:	I	: m/z: 50.00	5
						v	-15.0	Collision Energy:		n m/z: 100.00	Precurso
Scan Speed: 15000 u/sec Eyent Time: 0.030 sec Q1 Resolution: Unit Advanced	ed Settings	Advanced S		n: Unit	Q1 Resolution	sec	0.030	Event Time:	u/sec	geed: 15000	Sca
Q3 Resolution: Unit			-	n: Linit	O3 Resolution						

SHIMADZU

Fast polarity switching and MRM triggered automatic MS/MS applied to benzodiazepines and their metabolites in clinical and forensic analysis



Figure 4. MRM chromatograms of 4 compounds (each 100 ng/mL) spiked into urine and analyzed by Nexera using Shim-pack XR-ODS III coupled to LCMS-8030 after sample preparation. As the LC/MS/MS system has a high speed of data acquisition, the assay generates both MRM and Product Ion Scan (MS/MS) spectra resulting in quantitative data and library searching/product matching to help product confirmation. Fast polarity switching helps to provide information rich product ion spectra resulting in better detection and identification for each benzodiazepine.

Conclusion

- With Nexera using Shim-pack XR-ODS III coupled to LCMS-8030 provides significant advantage over other method for benzodiazepine analysis: fast analysis time (all compounds are eluted in a retention window less than 3 minutes) and detected in both positive and negative ion using a single analytical run.
- To help forensic chemists this high speed MRM triggered automatic MS/MS and a new method package and database enabled simultaneous screening and quantification of benzodiazepines and their metabolites.
- This method will be applied to the forensic analysis of urine samples taken to confirm administration in cases of benzodiazepine abuse.

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. We maintain a global network of sales, service, technical support and applications centers on six continents, and have established long-term relationships with a host of highly trained distributors located in over 100 countries. For information about Shimadzu, and to contact your local office, please visit our Web site at **www.shimadzu.com**



SHIMADZU CORPORATION. International Marketing Division 3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan Phone: 81(3)3219-5641 Fax. 81(3)3219-5710 URL http://www.shimadzu.com

For Research Use Only. Not for use in diagnostic procedures. Shimadzu Corporation ("Shimadzu") reserves all rights including copyright in this publication. The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to, or arising out of the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.