

Mice fecal metabolomics by LC-MS/MS: Comparison between young mice and old mice

ASMS 2017 MP 455

Takanari Hattori^{1,2}, Shuichi Kawana¹, Yoshihiro Hayakawa¹,
Eiichiro Fukusaki³, Mitsuharu Matsumoto⁴

1 Shimadzu Corporation, Kyoto, Japan.

2 Osaka University Shimadzu Analytical Innovation
Research Laboratory, Osaka, Japan.

3 Osaka University, Osaka, Japan.

4 Kyodo Milk Industry Co. Ltd., Tokyo, Japan

Mice fecal metabolomics by LC-MS/MS: Comparison between young mice and old mice

Introduction

The intestinal microbiome plays an important role in health and/or disease because it influences pathological and normal homeostatic functions. Low-molecular-weight metabolites produced by intestinal microbiome are absorbed constantly from the intestinal lumen and carried to systemic circulation; they play a direct role in health and/or disease. There are limited reports concerning the function of metabolites produced by intestinal microbiome.

Furthermore, these studies are targeting specific metabolites such as short chain fatty acids but not global metabolites (metabolome). For clarifying the relationship between health and/or disease and metabolome produced by intestinal bacteria, only free bacterial metabolites in the intestinal luminal content (or feces) should be analyzed. In this study, metabolites in feces of young and old mice were analyzed by LC-MS/MS.

Methods and Materials

Sample Preparation

Mice feces were diluted with phosphate-buffered saline and extracted by intense mixing. The upper portion was centrifuged and the supernatant was centrifugally filtered.

Figure 1 shows the detailed procedures.




Figure 1 Process flow of sample preparation

Mice fecal metabolomics by LC-MS/MS: Comparison between young mice and old mice

Experimental

Two LC-MS/MS methods (ion-pairing and non-ion-pairing method) were used for analysis. Nexera X2 system coupled with a LCMS-8040 or a LCMS-8050 (Shimadzu Corporation, Japan) was used. In ion-pairing method, separation was achieved on Mastro C18 column (Shimadzu GLC Ltd., Japan). The mobile phase consisted of 15 mM acetic acid and 10 mM tributylamine in water

(A) and methanol (B) under gradient elution. In non-ion-pairing method, separation was achieved on Discovery HS F5-3 column (Sigma-Aldrich, U.S.A.). The mobile phase consisted of 0.1% formic acid in water (A) and 0.1% formic acid in acetonitrile (B) under gradient elution.

Analytical Condition (Ion-pairing method)

UHPLC (Nexera X2 system)	
Column	: Mastro C18 (150 mmL.x2.0 mmI.D., 2.0 µm)
Mobile phase A	: 15 mM Acetate, 10 mM Tributylamine/water
B	: Methanol
Flow rate	: 0.3 mL/min
Injection vol.	: 3 µL
Column temp.	: 40°C
MS (LCMS-8040)	
Ionization	: Negative, MRM mode
DL temp.	: 250°C
HB temp	: 400°C
Drying gas	: 10 L/min
Nebulizing gas	: 2.0 L/min



Analytical Condition (Non-ion-pairing method)

UHPLC (Nexera X2 system)	
Column	: Discovery HS F5 (150 mmL.x2.1 mmI.D., 3.0 µm)
Mobile phase A	: 0.1% Formate/water
B	: 0.1% Formate/acetonitrile
Flow rate	: 0.25 mL/min
Injection vol.	: 3 µL
Column temp.	: 40°C
MS (LCMS-8050)	
Ionization	: Positive/Negative, MRM mode
DL temp.	: 250°C
HB temp.	: 400°C
Interface temp.	: 400°C
Nebulizing gas	: 2.0 L/min
Drying gas	: 10 L/min
Heating gas	: 10 L/min



Mice fecal metabolomics by LC-MS/MS: Comparison between young mice and old mice

Result

Analysis of Metabolites in Mice Fecal by LC-MS/MS

We evaluated 55 and 97 of metabolites in ion-pairing method and non-ion-pairing method respectively. In ion-pairing method, 17 of metabolites were detected in the extract of mice feces (peak area RSD <20%). Main compounds of 17 metabolites were amino acids. In non-ion-pairing method, 75 of metabolites were detected in the 10-fold diluted extract of mice feces

(peak area RSD <20%). Main compounds of 75 metabolites were amino acids and nucleic acid-related substances. Figure 2 shows the MRM chromatogram of extract of mice feces obtained using LC-MS/MS. The list of detected compounds and their repeatability is shown in Table 1.




Figure 2 MRM chromatograms of extract of mice feces

Mice fecal metabolomics by LC-MS/MS: Comparison between young mice and old mice

Table 1 Detected compounds and their repeatability

Category	Compounds	Repeatability (%RSD, n=6)		Category	Compounds	Repeatability (%RSD, n=6)		
		R.T.	Peak area			R.T.	Peak area	
Amino acid	4-Hydroxyproline	0.11	3.8	Nucleosides and Nucleotides	Guanine	0.17	5.8	
	Alanine	0.10	3.0		Guanosine	0.32	2.1	
	Arginine	0.07	2.6		Guanosine monophosphate	0.16	6.4	
	Asparagine	0.10	1.4		Inosine	0.33	1.5	
	Aspartic acid	0.07	1.7		Thymidine	0.26	10.1	
	Asymmetric dimethylarginine	0.18	1.5		Thymidine monophosphate	0.18	13.1	
	Citrulline	0.09	2.1		Thymine	0.05	6.9	
	Cystine	0.16	7.2		Uracil	0.12	3.1	
	Dimethylglycine	0.07	5.5		Uridine	0.07	3.7	
	Glutamic acid	0.09	1.6		Uridine diphosphate	0.03	8.1	
	Glutamine	0.10	1.0	Xanthine	0.15	7.6		
	Glycine	0.08	4.0	Organic acid	4-Aminobutyric acid	0.06	3.9	
	Histidine	0.06	4.7		Argininosuccinic acid	0.14	14.1	
	Isoleucine	0.13	1.7		Cholic acid	0.02	2.6	
	Leucine	0.10	2.4		Creatine	0.28	2.5	
	Lysine	0.08	0.7		Lactic acid	0.04	5.0	
	Methionine sulfoxide	0.08	2.1		Nicotinic acid	0.10	3.4	
	Ornithine	0.06	2.8		Pantothenic acid	0.16	2.6	
	Phenylalanine	0.07	1.0		Taurocholic acid	0.04	4.1	
	Proline	0.04	2.3		Uric acid	0.15	7.0	
	Serine	0.09	1.9		TCA cycle	2-Ketoglutaric acid	0.07	2.9
	Symmetric dimethylarginine	0.22	4.3	Aconitic acid		0.14	2.4	
	Threonine	0.09	2.4	Citric acid		0.13	6.4	
	Tryptophan	0.04	1.2	Isocitric acid		0.12	5.0	
	Tyrosine	0.13	0.9	Malic acid		0.10	2.8	
	Valine	0.05	4.8	Succinic acid		0.06	2.2	
	Nucleosides and Nucleotides	Adenine	0.20	2.1	Methylation and Transsulfuration cycle	Cystathionine	0.14	5.1
		Adenosine	0.23	9.6		Cysteine	0.20	7.1
		Adenosine 3',5'-cyclic monophosphate	0.42	14.8		Methionine	0.06	2.5
		Adenosine monophosphate	0.14	3.6	5-Adenosylhomocysteine	0.12	8.9	
Cytidine		0.23	1.2	Others	2-Aminobutyric acid	0.24	4.8	
Cytidine monophosphate		0.07	1.3		6-Phosphogluconic acid	0.04	5.6	
Cytosine		0.12	9.2		Acetylcarnitine	0.15	3.0	
					Allantoin	0.11	4.2	
			Carnitine		0.12	6.9		
			Carnosine		0.16	2.0		
			Choline		0.05	1.6		
			Creatinine		0.06	3.0		
			Dopamine		0.26	6.4		
			FMN		0.11	8.0		
			Histamine		0.22	2.6		
			Hypoxanthine		0.04	1.4		
			Norepinephrine	0.54	1.7			
			Serotonin	0.06	2.4			

* Red: the repeatability of compounds using ion-pairing method
Black: the repeatability of compounds using non-ion-pairing method

Mice fecal metabolomics by LC-MS/MS: Comparison between young mice and old mice

Comparison between Young Mice and Old Mice

As an application to fecal metabolomics, feces derived from mice that the age was different (10-week-old and 70-week-old) were analyzed by non-ion-pairing method. As results of analyses for the extracts of the mice feces, 66 of metabolites were detected. By Traverse MS software, principal component analyses (PCA) were performed. As shown in Figure 3, two mice feces were successfully classified. The concentrations of 21 in old mice (70-week-old) feces were obviously lower ($p < 0.01$) than those in young mice (10-week-old) feces. The lower concentration of fecal cholic acid in old mice feces indicates that the secretion of bile acid was decreased by the ageing. On the other hand, there is a likelihood that

the number of bacteria that deconjugate bile acid (taurocholic acid) to taurine and cholic acid, is more in young mice than old mice. In contrast, the concentrations of 26 in old mice feces were obviously higher ($p < 0.01$) than those in young mice feces. Intestinal luminal arginine is absorbed and converted to polyamines that are essential substances for the proliferation of somatic cells, through ornithine. The lower concentration of arginine in young mice feces indicates that colonocytes of young mice absorb luminal arginine actively to grow compared with old mice. Table 2 shows the list of compounds that showed significant differences between young mice and old mice.




Figure 3 Score plot and loading plot

Mice fecal metabolomics by LC-MS/MS: Comparison between young mice and old mice

Table 2 Compounds that showed significant differences between young mice and old mice

Compounds			
Young mice > Old mice		Old mice > Young mice	
Adenine	Proline	4-Aminobutyric acid	Guanosine monophosphate
Adenosine	Succinic acid	Adenosine monophosphate	Histamine
Aspartic acid	Thymidine	Arginine	Histidine
Carnitine	Thymidine monophosphate	Asymmetric dimethylarginine	Hypoxanthine
Cholic acid	Thymine	Citrulline	Lactic acid
Choline	Uracil	Cystine	Lysine
Cystathionine	Uric acid	Cytidine	Methionine
Cytosine	Xanthine	Cytidine monophosphate	Methionine sulfoxide
Fumaric acid		Dimethylglycine	Nicotinic acid
Glutamic acid		FMN	S-Adenosylmethionine
Glycine		Glutamine	Serotonin
Homosysteine		Guanine	Symmetric dimethylarginine
Isoleucine		Guanosine	Threonine

Conclusions

- We analyzed metabolites in mice fecal by LC-MS/MS (ion-pairing method and non-ion-pairing method) and 17 and 75 of metabolites were detected, respectively.
- As results of analyses by non-ion-pairing method and principal component analysis (PCA), two types of mice fecal (young and old) was successfully categorized.
- These analytical methods could be a very powerful tool in fecal metabolomics.

Disclaimer: The products and applications in this presentation are intended for Research Use Only (RUO).
Not for use in diagnostic procedures. Not available in China.

First Edition: June, 2017