Application News

No. 01-00949-ENK

GC-MS GCMS-QP2020NX

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■ Reference

- 1) IEC 62321-8:2017, *Determination of certain substances in electrotechnical products -Part 8: Phthalates in polymers by gas chromatography-mass spectrometry(GC-MS),gas chromatography mass spectrometry using a pyrolyzer/thermal desorption accessory (Py/TD-GC-MS)*
- IEC 62321-3-3:2021, *Determination of certain substances in electrotechnical products - Part 3-3: Screening - Polybrominated biphenyls, polybrominated diphenyl ethers and phthalates in polymers by gas chromatography-mass spectrometry using a pyrolyzer/thermal desorption accessory (Py/TD-GC-MS)*
- 3) [Persistent, Bioaccumulative, and Toxic \(PBT\) Chemicals under TSCA](#)
- 4) [European Commission. Commission Delegated Regulation \(EU\) 2025/843 of 5 May 2025 amending Annex I to Regulation \(EU\) 2019/1021 of the European Parliament and of the Council as regards UV-328\(C/2025/2567\)](#)
- 5) [European Commission. Commission Delegated Regulation \(EU\) 2019/1021 of the European Parliament and of the Council as regards dechlorane plus \(C/2025/2887\).](#)
- 6) [European Commission. Commission Regulation \(EU\) 2015/2030 of 13 November 2015 amending Regulation \(EC\) No 850/2004 of the European Parliament and of the Council on persistent organic pollutants as regards Annex I \(Text with EEA relevance\)\(32015R2030\)](#)
- 7) [European Commission. Commission Delegated Regulation \(EU\) 2019/1021 of the European Parliament and of the Council as regards medium chain chlorinated paraffins\(Ares/2025/10131503\)](#)

<Related Application News Articles>

1. Analysis of PIP (3:1) in Plastic by Py/TD-GC-MS, [Application News No. 01-00191](#)
2. Analyzing UV-328 and Dechlorane Plus in Plastic by Pyrolyzer/ Thermal Desorption-GC-MS, [Application News No. 01-00488](#)

Shimadzu Corporation
www.shimadzu.com/an/
Shimadzu Scientific Korea
www.shimadzu.co.kr

For Research Use Only. Not for use in diagnostic procedures. Not available in the USA, Canada, and China.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Copyright © 2026 SHIMADZU group. All rights reserved.