

Differential Transformer Type Extensometers

DT Series

ISO Class 1 Compliant Differential Transformer Type Extensometers




UH-X/FX Series
Hydraulic Universal
Testing Machines
200 kN to 4,000 kN



AG-X Series Floor Type
Precision Universal Testers
250 kN/300 kN




AGS-X Series Floor Type
Precision Universal Testers
300 kN



The ISO 6892-1 standard for tensile testing of metallic materials prescribes Class 1 or higher extensometers, and JIS Z 2241 prescribes Class 2 or higher extensometers when measuring withstanding pressure. The DT series includes Class 1 compliant differential transformer type extensometers, which can be used to implement test methods compliant with metal tensile test standards.

Extensometer Class 1 Compliant

This instrument is compliant with Class 1 in ISO 9513 and JIS B 7741 ($\pm 3 \mu\text{m}$ or $\pm 1\%$ of the specified value, whichever is larger). It is also compliant with ASTM E83 Class C.




Compliant with Strain Rate Control Test Methods

ISO 6892

JIS Z 2241

Perform tests in compliance with "Test Rates by Strain Rate Control" from JIS B 7741 (Attachment JB (Reference)) and ISO 6892-1 (Method A).

Note: Requires an X type controller.



Compliance with Extensometer Class 1 and Strain Control

Please reference the table below when using DT series differential transformer type extensometers for Class 1 compliance and when implementing strain rate control.

	UH-X/FX/ AG-X/ AGS-X Series	UH-I Series / AG-I Series	Previous Testing Machines
Class 1	Use the amp built into the UH-X/FX and AG-X, or an ESA amp.	Use an ESA amp. Precision: Class 1 corresponds to the following measurement ranges. • DT50-10: 0 mm to 1 mm • DT50-5: 0 mm to 1.25 mm • DT50-2: 0 mm to 1 mm • DT50-50 is not compatible.	Retrofit the X type controller, and then use the amp built into the UH-X/FX and AG-X, or an ESA amp.
Strain control		Retrofit the X type controller, and then use the amp built into the UH-X/FX and AG-X, or an ESA amp.	

Specifications

Measurement range: Max. 25 mm (4 types)

Gauge length: 50 mm

Precision: Extensometer Class 1

(Within $\pm 3 \mu\text{m}$ or $\pm 1 \%$ of the specified value, whichever is larger)

*1) This Precision is the value when calibrated simultaneously with the amp. Existing amps will require on-site calibration.


*2) The latest version of amp built into the UH-X/FX and AG-X, or the latest ESA amp, is required for compliance with these DT extensometers. Existing amps may not be compliant.

*3) These DT extensometers do not accommodate fractures.


Size and Measurement Range for Applicable Samples

Model	Structure	Applicable Sample ^{Note)}		Measurement Range (mm)
		Rod Diameter (mm)	Flat Plate W x T (mm)	
DT50-50-45	A	$\varnothing 23$ to 45	W 40 (max.) x T 23 to 45	25
DT50-10-45	C			5
DT50-5-45	C			2.5
DT50-2-45	C			1
DT50-50-25	A	$\varnothing 6$ to 25	W 40 (max.) x T 6 to 25	25
DT50-10-25	C			5
DT50-5-25	C			2.5
DT50-2-25	C			1
DT50-10-10	B	$\varnothing 3$ to 10	W 40 (max.) x T 0.2 to 10	5
DT50-5-10	B			2.5
DT50-2-10	B			1


Note: Samples should have enough rigidity so they are not deformed when installing the extensometer.



DT50-50-45
(Structure A: 4-point holding system)



DT50-10-10
(Structure B: 4-point holding system)



DT50-2-25
(Structure C: 3-point holding system)



SHIMADZU Corporation

www.shimadzu.com/an/

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation and its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

For Research Use Only. Not for use in diagnostic procedures. The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.

© Shimadzu Corporation, 2013

Printed in Japan 3655-11313-10AIT